



# Technical / News Bulletin

UK

BULLETIN NO.:

TB 0444

DATE ISSUED

30/10/14

## BULLETIN DETAILS

**Product Group:** Power

**Model Affected:** Quickie Jive and Salsa HD family

**Effective Date:** 10/07/2014

**Serial No. Effective From:** N/A

**Bulletin Subject:** New single post armrest receiver

### Subject of Bulletin:

As part of Sunrise Medicals commitment to improving our products and listening to feedback from our customer base, a new improved single post armrest has been introduced on the Quickie Jive and Quickie Salsa HD seating system.

The new design incorporates a strengthened post, new stronger receiver and is now secured with an 8mm thumb screw to replace the old 6mm thumb screw. The shape of the receiver and the larger thumb screw ensure that there is no play at all in the mechanism. Additionally, a steel pin has been inserted into the aluminium post which engages with the thumb screw, ensuring that there is no wear over time and no need to tighten up the screw on a regular basis.

The new assembly is completely retro-fittable to existing chairs, and has new part numbers as detailed below. Parts manuals will be updated to include all new part numbers, which replace existing items.



New assembly

New receiver

New strengthened post with steel pin

| Old Part # | Description                   | New replacement part # | Description                      |
|------------|-------------------------------|------------------------|----------------------------------|
| 21004079   | Base plate fixing             | 21004106               | Lower single post fixing LH side |
|            |                               | 21004107               | Lower single post fixing RH side |
| 21004086   | Single post outer assembly LH | 21004113               | Single post outer assembly LH    |
| 21004087   | Single post outer assembly RH | 21004114               | Single post outer assembly RH    |

**Literature Affected:**

**Sunparts:** Updated with new part numbers

**Order Forms:** No change

**Website:** No change

**Owners Manual:** No change

If you have any questions regarding this bulletin, please do not hesitate to contact:

| Customer Services | Customer Support                         | Technical Service Centre                 | E-mail   |
|-------------------|--|--|--|
| <b>Telephone:</b> | <b>+44 (0)845 605 66 88<br/>Option 2</b> | <b>+44 (0)845 605 66 88<br/>Option 3</b> | <a href="mailto:enquiries@sunmed.co.uk">enquiries@sunmed.co.uk</a>           |
| <b>Fax:</b>       | <b>+44 (0)845 605 66 89</b>              | <b>+44 (0)845 605 66 89</b>              | <a href="mailto:help.technical@sunmed.co.uk">help.technical@sunmed.co.uk</a> |