Argon H

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Foreword (1.1.1)

Dear Customer,

We are very happy that you have decided in favour of a high-quality product from SUNRISE MEDICAL.

This user's manual will provide numerous tips and ideas so that your new wheelchair can become a trustworthy and reliable partner in your life.

Maintaining close links with our customers is of great importance to us at Sunrise Medical. We would therefore like to keep you up-to-date with our new and current developments. Keeping **close** to our customers also means fast service when you need replacement parts or accessories, or just have a question about your wheelchair — and with as little red tape as possible.

We want you to be satisfied with our products and service. Sunrise Medical therefore constantly works at continuous development of its products. For this reason, changes can occur in our range of products with regard to shape, technology, and fittings. Consequently, no claims can be construed from the data or pictures contained in this user's manual.

SUNRISE MEDICAL has been awarded the ISO 9001 Certificate, which affirms the quality of our products at every stage, from R & D to production.

As the manufacturer, SUNRISE MEDICAL, declares that the lightweight wheelchairs conform to the 93/42/EEC / 2007/47/EEC guideline.

Please contact your local, authorised SUNRISE MEDICAL dealer if you have any questions concerning the use, maintenance, or safety of your wheelchair.

In the case that there is no authorised dealer in your area or you have any questions, you can contact Sunrise Medical either in writing or by telephone (contacts are mentioned on the last page).

**Sunrise Medical Ltd.
Sunrise Business Park
High Street, Wollaston
West Midlands DY8 4PS
England**

Telephone: +44/1384-446622

Fax: +44/1384-446644

www.sunrisemedical.com

IMPORTANT:

DO NOT USE YOUR WHEELCHAIR UNTIL THIS MANUAL HAS BEEN READ AND UNDERSTOOD.

Table of contents

Foreword for Wheelchairs

General safety notes and driving restrictions 4-5

Wheelchair Components 5

Handling for transportation
Quick-release axles on rear wheel 6

Options
Step Tubes 6
Wheel locks 6
Suspension system 6
Footboard 7
Footrest 7
Castor 7
Seat 7
Seat height 8
Wheel Alignment 8
Back 9
Sideguard 9
Push handles 10
Lap belt instructions 10
Anti-Tip Tubes 11
Seat sling 11
Seat Depth 11
Crutch Holder 11
Travel Wheels 11
Tyres and Mounting 11

Trouble Shooting 12

Maintenance 12

Technical Data 12-13

Nameplates/Guarantee 13

Torque 13

Use

Wheelchairs are exclusively for a user who is unable to walk or has limited mobility, for their own personal use in- and outdoor.

The maximum weight limit (includes the user and any weight of accessories fitted to the wheelchair) is marked on the serial number label, which is affixed to the crossbar or stabiliser bar below the seat.

Warranty can only be taken on if the product is used under the specified conditions and for the intended purposes.

The expected life of the wheelchair is 5 years.

Please **DO NOT** use or fit any 3rd party components to the wheelchair unless they are officially approved by Sunrise Medical.

Area of application

The variety of fitting variants as well as the modular design mean that it can be used by those who cannot walk or have limited mobility because of:

- Paralysis
- Loss of extremity (leg amputation)
- Extremity defect deformity
- Joint contractures/joint injuries
- Illnesses such as heart and circulation deficiencies, disturbance of equilibrium or cachexia as well as for elderly people who still have strength in the upper body.

When considering provision, please also note the body size, weight, physical and psychological constitution, the age of the person, living conditions and environment.

General safety notes and driving restrictions (3.2)

The engineering and construction of this wheelchair has been designed to provide maximum safety. International safety standards currently in force have either been fulfilled or exceeded in parts. Nevertheless, users may put themselves at risk by improperly using their wheelchairs. For your own safety, the following rules must be strictly observed.

Unprofessional or erroneous changes or adjustments increase the risk of accident. As a wheelchair user, you are also part of the daily traffic on streets and pavements, just like anyone else. We would like to remind you that you are therefore also subject to any and all traffic laws.

Be careful during your first ride in this wheelchair. Get to know your wheelchair.

Before each use, the following should be checked:

- Quick-release axles on the rear wheels
- Velcro on seats and backrests
- Tyres, tyre pressure and wheel locks.

Before changing any of the adjustments of this wheelchair, it is important to read the corresponding section of the user's manual.

It is possible that potholes or uneven ground could cause this wheelchair to tip over, especially when riding uphill or downhill.

When riding over a step or up an incline frontally, the body should be leaning forward.

Danger!

NEVER exceed the maximum load of 125 kg for driver plus any items carried on the wheelchair. If you exceed the maximum load, this can lead to damage to the chair, or you may fall or tip over, lose control and may lead to serious injury of the user and other people.

Danger!

When it is dark, please wear light clothing or clothing with reflectors, so that you can be seen more easily. Make sure that the reflectors on the side and back of the wheelchair are clearly visible. We would also recommend that you fit an active light.

Danger!

To avoid falls and dangerous situations, you should first practice using your new wheelchair on level ground with good visibility.

Danger!

When getting on or off the wheelchair, do not use the footboards. These should be flipped up beforehand and swung to the outside as far as possible. Always position yourself as close as possible to the place where you wish to sit.

Danger!

Only use your wheelchair properly. For example, avoid travelling against an obstacle without braking (step, kerb edge) or dropping down gaps.

Danger!

The wheel locks are not intended to brake your wheelchair. They are only there to ensure that your wheelchair does not begin rolling unintentionally. When you stop on uneven ground, you should always use the wheel locks to prevent such rolling. Always apply both wheel locks; otherwise, your wheelchair could tip over.

Danger!

Explore the effects of changing the centre of gravity on the behaviour of the wheelchair, for example on inclines, slopes, all gradients or when overcoming obstacles. Do this with the secure aid of a helper.

Danger!

With extreme settings (e.g. rear wheels in the most forward position) and less than perfect posture, the wheelchair may tip over even on a level surface.

Danger!

Lean your upper body further forward when going up slopes and steps.

Danger!


Lean your upper body further back when going down slopes and steps. Never try to climb and descend a slope diagonally.

Danger!

Avoid using an escalator which may lead to serious injury in the event of a fall.

Danger!

Do not use the wheelchair on slopes $> 10^\circ$. The Dynamic safe slope is dependant on the chair configuration, the users abilities and the style of riding. As the users abilities and style of riding cannot be pre-determined then the max safe slope cannot be determined. Therefore this must be determined by the user with the assistance of an attendant to prevent tipping. It is strongly recommended that inexperienced users have Ant-tips fitted.

Danger!

It is possible that potholes or uneven ground could cause this wheelchair to tip over, especially when riding uphill or downhill.

Danger!

Do not use your wheelchair on muddy or icy ground. Do not use your wheelchair where pedestrians are not allowed.

Danger!

To avoid hand injuries do not grab in between the spokes or between the rear wheel and wheel lock when driving the wheelchair.

Danger!

In particular when using lightweight metal handrims, fingers will easily become hot when braking from a high speed or on long inclines.

Danger!

Only attempt stairs with the help of an attendant. There is equipment available to help you, e.g. climbing ramps or lifts, please use them. If there is no such equipment available, then the wheelchair must be tipped and pushed, never carried, over the steps (2 helpers). We recommend that users over 100 kg in weight do not use this stairway manoeuvre!

In general, any anti-tip tubes fitted must be set beforehand, so that they cannot touch the steps, as otherwise this could lead to a serious tumble. Afterwards the anti-tip tubes must be set back to their correct position.

Danger!

Make sure that the attendant only holds the wheelchair using securely mounted parts (e.g. not on the footrests or the sideguards).

Danger!

When using the lifting ramp make sure that the anti-tip tubes fitted are positioned outside the danger area.

Danger!

Secure your wheelchair on uneven ground or when transferring (e.g. into a car) by using the brakes.

Danger!

If and whenever possible, during a journey in a specially fitted vehicle for disabled people, vehicle occupants should use the seats in the vehicle and the appropriate restraint system. This is the only way to ensure that occupants will have the maximum protection if there is an accident. When using safety elements offered by SUNRISE MEDICAL and using a specially designed safety system, lightweight wheelchairs can be used as a seat when being transported in a specially fitted vehicle. (See the Chapter on "Transportation").

Danger!

Depending on the diameter and setting of the castors, as well as

the centre of gravity setting of the wheelchair, the castors may begin to wobble at high speeds. This can lead to the castors being blocked and the wheelchair may tip over. Therefore, please make sure that the castors are adjusted correctly (see the Chapter "Castors"). In particular, do not travel on an incline without brakes, travel at a reduced speed. We recommend that novice users use anti-tip tubes.

Danger!

Anti-tip tubes should prevent the chair tipping over backwards unintentionally. Under no circumstances should they take the place of transit wheels, and be used to transport a person in a wheelchair with the rear wheels removed.

Danger!

When reaching for objects (which are in front of, to the side or behind the wheelchair) make sure that you do not lean too far out of the wheelchair, as if you change the centre of gravity there is a risk of tipping or rolling over. The hanging of additional load (back pack or similar items) onto your chair backposts can affect the rearward stability of your chair, especially when used in combination with recliner backrests. This can cause the chair to tip backwards causing injury.

Danger!

For thigh amputees you must use anti-tip tubes.

Danger!

Before setting off, check that your tyre pressure is correct. For rear wheels it should be at least 3.5 bar (350 kPa). The max. pressure is indicated on the tyre. The knee-lever brakes will only work if there is sufficient tyre pressure and if the correct setting has been made (see the Chapter on "Brakes").

Danger!

If the seat and back sling are damaged, you must replace them immediately.

Danger!

Be careful with fire, in particular with burning cigarettes. Seat and back slings can be set alight.

Danger!

If the wheelchair is subject to direct sunlight for a long period of time, then parts of the wheelchair (e.g. frame, legrests, brakes and sideguard) may become hot (>41°C).

Danger!

Always make sure that the quick-release axles on the rear wheels are set properly and lock in. If the button on the quick-release axle is not pressed in, the rear wheel cannot be removed.

Note !

The effect of the knee-lever brake as well as the general driving characteristics are dependant on tyre pressure. The wheelchair is significantly lighter and easier to manoeuvre when the rear wheels are pumped up correctly and both wheels have the same pressure.

Note !

Make sure that your tyres have sufficient tread! Please note that you are subject to any and all traffic laws when driving in public traffic.

Note !

Always be careful with your fingers when working or adjusting the wheelchair!

The products shown and described in this manual may not be exactly the same in every detail as your own model. However, all instructions are still entirely relevant, irrespective of detail differences. The manufacturer reserves the right to alter without notice any weights, measurements or other technical data shown in this manual. All figures, measurements and capacities shown in this manual are approximate and do not constitute specifications.

Transportation

Your wheelchair is not suitable for use as a seat in vehicles.

Do not sit in the wheelchair during transportation. Use the vehicle seat or a properly secured transport approved seat.

Wheelchair Components (4.1)

We at SUNRISE MEDICAL have been awarded the ISO 9001 Certificate, which affirms the quality of our products at every stage, from R & D to production. This product complies with the standards set forth in EU directives. Optional equipment and accessories are available at extra charge.

Wheelchairs:


1. Push handles
2. Backrest upholstery
3. Sideguards
4. Seat sling
5. Footrest
6. Castors
7. Castor plates
8. Fork (6" castor)
9. Quick-release axles
10. Wheel locks
11. Handrim
12. Rear wheel

Handling

Quick-Release Axles for Rear Wheels (6.2)

The rear wheels are equipped with quick-release axles. The wheels can thus be installed or removed without using tools. To remove a wheel, simply depress the quick-release button on the axle (1) and pull it out.

CAUTION!

Hold the quick-release button on the axle depressed when inserting the axle into the frame to mount the rear wheels. Release the button to lock the wheel in place. The quick-release button should snap back to its original position.

Transporting the chair

Transporting the chair (6.5)

Removing the rear wheels will keep the chair as compact as possible. The backrest can be folded down by pulling the cord located on the backrest.

503

514

Options – Step Tubes

Step Tubes (7.1)

Step tubes are used by attendants to tip a wheelchair over an obstacle. Simply step on the tube to push a wheelchair, for example, over a curb or step.

519

Options – Wheel Locks

Wheel locks (7.10)

Your wheelchair is equipped with two wheel locks. They are applied directly against the tyres. To engage, press both wheel lock levers forward against the stops. To release the wheel, pull the levers back to their original positions.


Braking power will decrease with:

- Worn tyre tread
- Tyre pressure that is too low
- Wet tyres
- Improperly adjusted wheel locks.

520

The wheel locks have not been designed to be used as brakes for a moving wheelchair. The wheel locks should therefore never be used to brake a moving wheelchair. Always use the handrims for braking. Make sure that the interval between the tyres and wheel locks complies with the given specifications. To readjust, loosen screw (1) and set the appropriate interval. Tighten screw (see the page on torque).

520

CAUTION!

After each adjustment of the rear wheels, check the interval to the wheel locks and readjust if necessary.

Extension for the wheel lock Lever (7.11)

The longer lever helps to minimize the effort needed to set the wheel locks. The brake lever extension is screwed to the brakes. By raising this, it can be flipped forward.

CAUTION!

Incorrect mounting of the wheel lock will result in a higher effort to operate. This might cause the wheel lock extension lever to break.

521


Leaning onto the wheel lock extension lever while transferring will cause the lever to break. There is a risk it could break. Splashing water from tyres might cause the wheel lock to malfunction.

Options – Suspension system

Suspension system (7.17)

The function of the suspension system is determined by the buffer elements (1). Sunrise Medical offers a varied range of elements to suit the individual weight of the user.

To replace the elements, remove the 2 screws (2) (1 on each side). Turn the rocker arms (3) downwards; the buffer elements (1) can then be easily removed from the upper and lower openings.


522

To fit new sleeves, use the reverse procedure. Ensure that the buffer elements are securely located in the upper or lower openings.

Options – Footplate Adjustment

Adjusting the Footrest (7.25)

Removing screw (1) will allow you to adjust the footrest to correspond to the length of your lower leg. The angle of the footrest may be individually adjusted (not rigid footrest) by loosening screws (2). The side protection on the footrest prevents the feet slipping accidentally. When finished, make sure that all screws have been properly tightened (see the page on torque).


523

Options - Seat

Adjusting the seat height (7.53.2)

To adjust rear seat height, slacken and remove the 4 Torx screws (1) (2 each side) and the thread washer (2), which secures the clamp (3) for the camber tubes to the axle plates (4). Adjust the two camber tube clamps (3) to the required height and replace the 4 Torx screws. Before tightening the screws, please follow the instructions for setting the toe-in/toe-out to zero (section 7.61.1). Tighten the screws to 7 Nm.

NOTE – An adjustment to the castor angle may be necessary when adjusting the rear seat height.

Options - Castor

Castors, Castor Plates, Forks (7.40)

From time to time the wheelchair may veer slightly to the right or left, or the castors may flutter. This may be caused by the following:

- Forward and/or reverse wheel motion has not been set properly.
- The castor angle has not been adjusted properly.
- Castor and/or rear wheel air pressure is incorrect; wheels do not turn smoothly.

The optimum adjustment of the castors is required for the wheelchair to run in a straight line. Castors should always be adjusted by an authorized dealer. The castor plates must be readjusted, and the wheel locks must be checked any time the rear wheel position has been altered.

Options - Castor

Adjusting the castor (7.42.3)

To ensure that both forks are set parallel, simply count the teeth visible on both sides.

After setting the castor fork, the teeth will guarantee a secure position, allowing an adjustment of 16° in 2° increments.

Use the flat side to check for a right-angled position to the ground.

The patented design allows the castor fork to be turned, so that it can be reset at right-angles to the ground when the seat angle is changed.


524

525

526


195

Setting the directional stability

Remove the cover cap (1). Then release the Allen screw (2) inside. You can then open the screw (3). You can now turn the black socket (4) left or right.

Left – chair pulls to the left
Right – chair pulls to the right

Then tighten up the screws (3) again. Please set a 90° angle from the fork to the floor. Now tighten up the screws (2) again and close the tube with the cover cap (1).

Options – Wheel Alignment

Adjusting Wheel Alignment (7.61)

Important: Your chair will only roll properly if the rear wheel positions have been optimized, which means correctly adjusting the wheel alignment. To do this, measure the distance between both wheels front and rear to ensure that they are parallel to one another. The difference between both measurements should not exceed 5 mm. To adjust the wheels to make them parallel, loosen the screws and turn the axle sleeve accordingly. When finished, make sure that all screws have been properly tightened (see the page on torque).

Argon TI tracking adjustment

3. Setting the toe-in/toe-out to zero


NOTE: On a wheelchair with 0° camber cylinders it is not possible to set toe-in or toe-out. This setting is necessary only with 3°, 6° and 9° camber cylinders.

The term "toe-in or toe-out" defines how well the rear wheels of the chair are aligned in relation to the ground. This determines how well the chair will run. Normal resistance or rolling resistance is provided, when toe-in is set to zero.


To set toe-in/toe-out to zero:

Slacken the 2 screws (1) (1 each side) which secure the angle tube clamp. Check the ball in the horizontal (2) plane and turn the angle tube (3) until the ball is in the centre. Toe-in is now zero.

Before tightening screws (1), check that the camber tube is centred left-to-right. The gap should be the same on both sides, or there should be no gap at all. Tighten the screws to 7 Nm.

527

528


Width adjustment of the wheelbase (7.62)

K. REAR WHEELBASE

The rear wheelbase is defined as the distance between the upper side of the rear wheels and the backrest tubes, and is represented by measurement X. The factory setting is (1.25 cm). A bigger gap is usually required when a sufficient gap between the tyres and optionally height-adjustable armrests has to be created.

NOTE: When adjusting the rear wheelbase, adjust first one wheel then the other. If both sides are slackened at the same time, the adjustment of toe-in/toe-out will be altered.

To adjust the rear wheelbase, the parts of the camber (4) move telescopically into or out of the camber tube (5), and lock into place when they reach the end. Slacken screw (6) (located closest to the camber tube) on the left side of the chair. Move the quick-release axle inwards or outwards to achieve the desired wheelbase. Tighten the screws to 7 Nm. Repeat this procedure on the right side of the chair and adjust the gap so that it is the same amount as on the left side.

Options – Back

Angle adjustment of the folding backrest:

1. Open the upper screws and loosen the screw connection (Figure 1).
2. The hole in the back tube must be lined up with the hole in the connecting part, such that you can set the desired backrest angle (Figure 2).
3. Fit the screw connection with the nut and washer until they are hand tight, so that there is no lateral play between the components. But you must still be able to fold the back tube down easily.
4. Set the folding mechanism: To do this, loosen the nuts on the cam slightly (Figure 3). Set the cam (component 1) such that the folding mechanism locks into place with no play (if necessary use a 10 mm open-ended spanner to do this). Then hold the cam in position and tighten the nuts (5 Nm).
5. Repeat steps 1 to 4 for the second side.


Figure 1


Figure 2



Figure 3

Adjustable back sling (7.77.2)

The adjustable back-sling can be adjusted for tension by using several straps. The back-sling upholstery can be accessed from the inside via an opening and can be upholstered to suit individual tastes.

529

Height-adjustable backrest (7.77)


The back height may be set to various back heights, in 2.5 cm steps. Adjustment ranges (25-30 cm, 30-40 cm, 40-45 cm). Release the screw (1) and set the backrest to the desired height. Tighten up the screws again (see the page on torque).

506

Padded swingaway armrests (7.81.3)

The swing-away armrests can either be removed or swung away, so that it is easier to get into and out of the chair from the side. They can be height-adjusted (50 mm), by pushing the bolts (G) up or down in the pre-drilled holes on the armrest brackets. Fitting the swing-away and removable armrests is carried out by pushing the receiver tube (H) onto the rear frame tube.

530

Options – Sideguard

Quickie – height-adjustable sideguard (7.81.4) not part of the range in Germany

1. Fitting

- a. Push the outside armrest support into the clamp which is mounted on the chair frame.
- b. The armrest will automatically lock into place.

2. Height adjustment

- a. Turn the release lever to the second position.
- b. Push the arm support up or down to the required height.
- c. Turn the release lever back to the armrest locking position.
- d. Push the arm support in until the upper armrest locks into position.

3. Removing the armrest

- a. Turn the release lever to the first position and remove the armrest.

4. Changing the armrest

- a. Push the armrest back into the clamp.
- b. Turn the release lever back to the armrest locking position.

5. Adjusting the fit of the armrest receiver

To tighten or slacken the fit of the outside armrest in the clamp:

- a. Release the four bolts on the sides of the clamp.
- b. Hold the armrest in the clamp and squeeze the clamp together at the required fit.
- c. Tighten up the four bolts.

6. Adjusting the fit of the inner armrest

- a. The outside armrest is fitted with two stud bolts.
- b. Turn stud bolts in or out as required until you have found the desired fit.

516

Attention:

Sideguards and armrests are not designed to be used to lift or carry the wheelchair.

Options – Push handles


Height-adjustable Push Handles (7.90)

These handles are secured with pins to prevent them from sliding out unintentionally. Opening the quick-release lever (1) makes it possible to adjust the height of the push handles (3) to meet your individual needs. As you move the lever, you will hear a locking mechanism; you may now easily position the push handle as desired. The nut (2) on the tension lever determines how tightly the push handles are clamped into place. If the nut is loose after adjusting the tension lever, the push handle will also be too loose. Turn the push handle from side to side before use to make sure that it is clamped securely enough into place. After adjusting handle height, always clamp the tension lever (1) securely into place. If the lever is not secure, injuries could result when ascending stairs.

Fold-Down Push Handles (7.91)

If the push handles are not in use, they can be folded down depressing the button (2). When they are needed again, simply flip them back up until they click into place.

50

Options – Lap belt


Before using your wheelchair ensure the lap belt is worn.

The lap belt must be checked on a daily basis to ensure it is free from any obstruction or adverse wear.

Always make sure that the lap belt is correctly secured and adjusted prior to use. Too loose a strap could cause the user to slip down and risk suffocation or cause serious injury.

The lap belt is fitted to the wheelchair as shown in the illustrations. The seat belt comprises 2 halves. They are fitted using the existing seat stay retaining bolt fitted through the eyelet on the belt. The belt is routed under the rear of the side panel. (Fig. 24)

Adjust the belt position so buckles are in the centre of the seat. (Fig. 25)


Adjust lap belt to suit the user's needs as follows:

To reduce the belt length	To increase the belt length
Feed free belt back through male buckle and slide adjusters. Ensure belt is not looped at male buckle.	Feed free belt through slide adjusters and male buckle to provide more belt length.

When fastened check space between the lap belt and user. When correctly adjusted it should not be possible to insert more than the fat of the hand between the lap belt and the user. (Fig. 26)

The lap belt should be fixed so that the belt sits at an angle of 45 degrees across the users pelvis. The user should be upright and be as far back as possible in the seat when correctly adjusted. The lap belt should not allow the user to slip down in the seat. (Fig. 27)

To fasten buckle: Firmly push male buckle into female buckle.	To release belt: Press exposed sides of male buckle and push towards centre whilst gently pulling apart.

If in doubt about the use and operation of the lap belt then ask your healthcare professional, wheelchair dealer, carer or attendant for assistance.

Advice to client
The lap belt must only be fitted by an approved Sunrise Medical dealer / agent. The lap belt should only be adjusted by a professional, or a Sunrise Medical approved dealer / agent. The lap belt must be checked on a daily basis to ensure it is adjusted correctly and free from any obstruction or adverse wear. Sunrise Medical does not encourage the transportation of any person in a vehicle using this lap belt as a method of restraint. Please see Sunrise Medical transit booklet for further advice on transportation.

Maintenance:
Check lap belt, and securing components, at regular intervals for any sign of frays or damage. Replace if necessary.

NOTE:
The lap belt should be adjusted to suit the end user as detailed above. Sunrise Medical recommend that the length and fit of the belt be checked on a regular basis to reduce the risk of the end user inadvertently re-adjusting the belt to an excessive length.

Options – Anti-Tip Tubes

Anti-Tip Tubes

Sunrise Medical recommends anti-tip tubes for all chairs. When fitting anti-tip tubes, use a torque of 16 Nm.

1. Slotting the anti-tip tubes into the clamp
- a. Press the rear button on the anti-tip tube on the anti-tip tube adapter, so that both release pins are drawn inwards.
- b. Slot the anti-tip tubes (1) into the anti-tip tube adapter (2).
- c. Turn the anti-tip tubes downwards until the release pin locks into the clamp.
- d. Fit the second anti-tip tube in the same way.

53


2. Adjusting the anti-tip tubes
To achieve the correct ground clearance of approx. 1" to 2" (3.5 cm to 5.0 cm), the anti-tip tubes must be raised or lowered.

Press the anti-tip tube release button (3), so that both release pins are drawn inwards. Move the inner tube up or down to slot into the height holes provided (K). Release the button. Fit the second anti-tip wheel in the same way. Both wheels should be at the same height.

Options - Seat

Seat sling (7.101.3)

Remove the screws on the left side of the sling. Adjust the VELCRO® material to increase the tension of the seat sling. Tighten up the screws again. If the screws are difficult to fit back in place, try and locate the holes with a sharp object. Ensure also that the plastic base is in the correct position before the screws are tightened again.


Options – Seat Depth

Seat depth growth (Optional) (7.102.3)

Caution: This option is not available in Switzerland, due to legal patent reasons!

Using the optional cylinder for an offset seat, the backrest tubes can be moved either 1" (2.5 cm) or 2" (5 cm) further back than with the standard seat cylinder.

1. Before fitting, establish which back system (rigid back) is to be used and which seat offset is required.
2. If a seat depth increase of approx. 1" (2.5 cm) is required, screw (1) is screwed into hole 2.
3. If a seat depth increase of approx. 2" (5 cm) is required, screw (1) is screwed into hole 3.

Options – Crutch holder

Crutch holder (7.110)

This device permits crutches to be transported directly on the wheelchair. It has a Velcro loop to fasten crutches or other aids.

CAUTION!

Never try to use or even remove the crutches or other aids while moving.

306

Options – Travel Wheels

Travel Wheels (7.113)

Travel wheels should be used wherever your wheelchair would be too wide if the rear wheels were used. After the rear wheels have been removed with the aid of the quick-release axles, the transit wheels can immediately be used to continue riding. The transit wheels are mounted so that they are approx. 3 cm above the ground when not in use. They are thus out of the way when riding, transporting, or when tipping to pass over obstacles.

532

CAUTION!

Without rear wheels, your wheelchair has no wheel locks.

NOTE: When the chair is to be set up with the transit wheels and the anti-tip tubes, the transit mounting (G) must be installed between the camber tube clamp (I) and the anti-tip tube clamp mounting (not shown).

Tyres and Mounting

Tyres and mounting (8.1)

Always make sure that you maintain the correct tyre pressure, as this can have an effect on wheelchair performance. If the tyre pressure is too low, rolling resistance will increase, requiring more effort to move the chair forward. Low tyre pressure also has a negative impact on manoeuvrability. If the tyre pressure is too high, the tyre could burst.

The correct tyre pressure for a given tyre is printed on the surface of the tyre.

Tyres can be mounted the same way as an ordinary bicycle tyre.

Before installing a new inner tube, you should always make sure that the base of the rim and interior of the tyre are free of foreign objects. Check the pressure after mounting or repairing a tyre. It is critical to your safety and to the wheelchair's performance that regulation air pressure be maintained and that tyres are in good condition.

Trouble Shooting

Wheelchair pulls to one side (9.1.1)

- Check tyre pressure
- Check to make sure wheel turns easily (bearings, axle)
- Check to make sure both castors are making proper contact with the ground

Castors begin to wobble (9.2)

- Check angle of castors
- Check to make sure all bolts are secure; tighten if necessary (see the page on torque)
- Check to make sure both castors are making proper contact with the ground

Wheelchair squeaks and rattles (9.6)

- Check to make sure all bolts are secure; tighten if necessary (see the page on torque)
- Apply small amount of lubrication to spots where movable parts come in contact with one another

Wheelchair begins to wobble (9.7)

- Check angle of castors
- Check tyre pressure
- Check to see if rear wheels are adjusted differently

Maintenance and care

Maintenance (9.9)

- Check the tyre pressure every 4 weeks. Check all of the tyres for wear or damage.
- Check the brakes approximately every 4 weeks to make sure that they are working properly and easy to use.
- Tyres can be mounted the same way as an ordinary bicycle tyre.

Technical Data

Seat heights:

The choice of frames, forks and castors, as well as the rear wheel size (24", 26") determines the seat heights which can be achieved.

Possible seat heights

Important: Measurements without cushions!

Argon Ti – fixed castor receiver

Castors	Fork	Front seat height in cm	Rear seat height in cm	
			24"	26"
3" solid	72 mm	43	43-35	43-42
	72 mm	44	44-35	44-42
4" solid	118 mm	46	46-35	46-42
	118 mm	47	47-35	47-42
	118 mm	48	48-36	48-42
	118 mm	49	49-37	49-42
	118 mm	50	50-39	49-42
	138 mm	48	48-36	48-42
	138 mm	49	48-37	49-42
	138 mm	50	48-39	49-42
	138 mm	51	48-39	49-42
	138 mm	52	48-40	49-42
5" solid	118 mm	49	48-37	49-42
	118 mm	50	48-39	49-42
	118 mm	51	48-39	49-42
	138 mm	50	48-39	49-42
	138 mm	51	48-39	49-42
	138 mm	52	48-40	49-42
6" soft	138 mm	53	48-41	49-42
	118 mm	51	48-39	49-42
	118 mm	52	48-40	49-42
	138 mm	51	48-39	49-42
	138 mm	52	48-40	49-42
	138 mm	53	48-41	49-42

- All of the joints that are critical to using your wheelchair safely are fitted with self-locking nuts. Check every 3 months to make sure that all bolts are secure (see the page on torque).
 - Self-locking nuts should only be used once and should be replaced after single use.
 - Use only mild household cleansers when your wheelchair is dirty. Use only soap and water when cleaning the seat upholstery.
 - If your wheelchair should ever get wet, please dry it after use.
 - A small amount of sewing-machine oil should be applied to quick-release axles approximately every 8 weeks.
- Depending on the frequency and type of use, we recommend taking your wheelchair to your authorised dealer every 6 months, to have it inspected by trained personnel.

CAUTION!

Sand and sea water (or salt in the winter) can damage the bearings of the front and rear wheels. Clean the wheelchair thoroughly after exposure to these conditions.

Hygiene when being reused:

When the chair is to be reused, it should be prepared carefully, and be wiped and treated with spray disinfectant on all surfaces which could come into contact with the user. In doing this, you must use a disinfectant taken from the DGHM list, e.g. Antifect Liquid (Schülke & Mayr) for alcoholic quick disinfection for medical products and medical inventory when you require a quick effect. Please pay attention to the manufacturer's instructions of the disinfectant you are using. In general, a safe disinfectant cannot be guaranteed on seams. We therefore recommend, that you properly dispose of seat and back slings in the case of microbacterial contamination with active agents according to §6 infection protection law.

Technical Data

Overall width: 24" 26"

With 0° camber:	SW+20cm	SW+20cm
With 3° camber:	SW+22cm	SW+26cm
With 6° camber:	SW+28cm	SW+32cm
With 9° camber:	SW+34cm	SW+38cm

each with narrow-mounted handrims configuration

Overall length: 107 cm with seat depth 50

Overall height: 97 cm with back height 45

Weight in kg: from 7.6 kg

Maximum load:

Approved to load of 120 kg

Nameplates / Guarantee

Nameplates (11.1.2)

The nameplate is located on either the cross-tube assembly or the transverse frame tube. The name plate indicates the exact model designation as well as other technical specifications. Please provide the following information whenever you have to order replacement parts or to file a claim








- Serial number
- Order number
- Month/Year

Guarantee

We are very happy that you have decided in favour of a high-quality product from SUNRISE MEDICAL. As a sign of our gratitude, we are providing you with a 5 year guarantee on frame components. We are not responsible for any damage resulting from inappropriate or unprofessional installation and/or repairs, through neglect and wear, or from changes in any wheelchair components caused either by the user or by third parties. In such cases, this guarantee shall be considered null and void.

Custom wheelchairs cannot be exchanged.

Torque (12.1)

Torque for M6 screws is 7 Nm, unless otherwise specified.

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Sunrise Medical AB
Sunrise Medical B.V.
G2 A/S