HANDBIKE

# Shark RS/RT

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### Wheelchair Components

EN 


We at SUNRISE MEDICAL have been awarded the ISO-13485 certificate, which affirms the quality of our products at every stage, from R & D to production. This product complies with the standards set forth in EU and UK regulations. Options or accessories shown are available at extra cost.

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**HAND BIKE:**

1. REAR WHEEL
2. HEADREST
3. BACKREST UPHOLSTERY
4. GEAR-CHANGE
5. FRONT WHEEL
6. PARKING BRAKE
7. LEG SUPPORT
8. IMPACT PROTECTION

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
## Foreword

Dear Customer,

We are very happy that you have decided in favour of a high-quality product from SUNRISE MEDICAL.

This owner's manual will provide numerous tips and ideas so that your new hand bike can become a trustworthy and reliable partner in your life.

Maintaining close ties to our customers means a lot to us at Sunrise Medical. We like to keep you up-to-date about new and current developments at our company. Keeping close to our customers means: fast service, as little red tape as possible, working alongside our customers. When you need replacement parts or accessories, or if you just have a question about your hand bike – we are there for you.

We want you to be satisfied with our products and service. At Sunrise Medical we are constantly working to develop our products further. For this reason, changes can occur in our palette of products with regard to form, technology and equipment. Consequently, no claims can be construed from the data or pictures contained in this user's manual.

**The management system of SUNRISE MEDICAL is certified to ISO 13485 and ISO 14001.**

**UK  
CA** As the manufacturer, SUNRISE MEDICAL, declares that this product conforms to the UK Medical Devices Regulation 2002 No. 618.

**CE** As the manufacturer, SUNRISE MEDICAL, declares that this product conforms to the Medical Device Regulation (2017/745).

### NOTE:

#### General user advice.

Not following these instructions may result in physical injury, damage to the product or damage to the environment!

Notice to the user and/or patient: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

#### B4Me special adaptations

Sunrise Medical strongly recommends that in order to ensure that your B4Me product operates, and performs as intended by the manufacturer; all the user information supplied with your B4Me product is read and understood, before the product is first used.

Sunrise Medical also recommends that the user information is not discarded after reading it, but it is kept safely stored for future reference.

## Medical Device Combinations

It may be possible to combine this Medical device with one or more other Medical Device or other product. Information on which combinations are possible can be found at [www.Sunrisemedical.co.uk](http://www.Sunrisemedical.co.uk). All combinations listed have been validated to meet the General Safety and Performance Requirements, Annex I Nr. 14.1 of the Medical Device Regulation 2017/745.

Guidance on the combination, such as mounting, can be found at [www.SunriseMedical.co.uk](http://www.SunriseMedical.co.uk)

Please contact your local, authorised SUNRISE MEDICAL dealer if you have any questions concerning the use, maintenance, or safety of your hand bike.

If there is no authorised dealer in your area, or you have any questions, contact Sunrise Medical either in writing or by telephone.

### Sunrise Medical

Thorns Road

Brierley Hill

West Midlands

DY5 2LD

England

Phone: 0845 605 66 88

Fax: 0845 605 66 89

[www.SunriseMedical.co.uk](http://www.SunriseMedical.co.uk)

### IMPORTANT:

DO NOT USE YOUR HAND BIKE UNTIL THIS MANUAL HAS BEEN READ AND UNDERSTOOD.

## Area of Application

### Indications

The varieties of fitting variants, as well as the modular design, mean that it can be used by those who cannot walk or have limited mobility because of:

- Paralysis
- Loss of extremity (leg amputation)
- Extremity defect deformity
- Joint contractures/joint injuries
- Illnesses such as heart and circulation deficiencies, disturbance of equilibrium or cachexia as well as for elderly people who still have strength in the upper body.
- The Hand-bike device can be operated by people who are physically and mentally capable of safely operating the pushing device in all operating situations (uphill, downhill, etc.) and can follow the legal requirements for use on public roads.

### Contraindications

The wheelchair/ handcycle shall not be used in case of:

- Perception disorder
- Imbalance
- Loss of both arms
- Joint contracture or joint damage on both arms
- Seating disability

### LIFETIME

The intended lifetime of this product is 5 years. provided that: It is used in strict accordance with the intended use. All service and maintenance requirements are fulfilled.

## Table of Contents

Foreword	16
Area of Application	16
Table of Contents	17
Definitions	17
1.0 Safety and hazard notes	18
2.0 Warranty	19
3.0 Handling	20
4.0 Options	20
Brake	20
Footrest	21
RT gearshift	22
RS gearshift (Fig. 4.12)	22
Backrest adjustment	23
Height and depth adjustable fork	23
Rear axle	24
5.0 Maintenance and care (checks to be carried out before use)	24
6.0 Tyre maintenance	25
7.0.0 Disposal / Recycling of Materials	25
8.0 Nameplate	26
9.0 Torque	27

## Definitions

### Definitions of words used in this manual

Word	Definition
 <b>DANGER!</b>	Advice to the user of potential risk of serious injury or death if the advice is not followed
 <b>WARNING!</b>	Advice to the user of a potential risk of injury if the advice is not followed
 <b>CAUTION!</b>	Advice to user that potential damage to equipment may occur if the advice is not followed
<b>NOTE:</b>	General advice or best practice
	Reference to additional documentation

### NOTE:

Please keep a note of your local service agent's address and telephone number in the space provided. In the event of a breakdown, contact them and try to give all relevant details so they can help you quickly.

### NOTE:

Please note that driving a wheelchair requires sufficient cognitive, physical and visual skills. The user must be able to assess the effects of actions during the operation of the wheelchair and, if necessary, to correct them. These capabilities and the safe use of the additionally attached components cannot be assessed by Sunrise Medical as a manufacturer. We cannot accept any liability for any damage resulting from this.

Please refer to the operating instructions of the wheelchair and the additionally mounted components. Instruct the user in the safe use of the wheelchair and the additionally mounted components. Inform users of specific warnings that need to be read, understood, and respected.

Dealer signature and stamp

## 1.0 Safety and hazard notes

The technology and construction of this hand bike have been designed to provide maximum safety. International safety standards currently in force have either been fulfilled or exceeded in parts.

Nevertheless, users may put themselves at risk by improperly using their hand bikes. For your own safety, the following rules must be strictly observed.

Unprofessional or erroneous changes or adjustments increase the risk of accident. As a hand bike user, you are also part of the daily traffic on streets and pavements, just like anyone else. We would like to remind you that you are therefore also subject to any and all traffic laws. Be careful during your first ride in this hand bike. Get to know your bike.

Before each use, the following should be checked:

- Fixed axles or quick-release axles on the rear wheels.
- Velcro on seats and backrests
- Tyres, tyre pressure and wheel locks and main brake.
- Before each journey carry out a visual inspection of the frame, checking for cracks and other damage.

Before changing any of the adjustments of this hand bike, it is important to read the corresponding section of the Owner's manual.

NEVER exceed the maximum load of 125 kg for driver plus any items carried on the hand bike. If you exceed the maximum load, this can lead to damage to the hand bike, loss of control, or to serious injury to the user and other people.

- Exercise caution whenever you use the hand bike. For example, avoid travelling against an obstacle without braking (step, kerb edge) or dropping down gaps.
- Your hand bike has not been designed for transportation in a vehicle. DO NOT sit on the hand bike during any transportation in a vehicle; always use a properly secured vehicle seat!!
- This hand bike should be used exclusively to convey one person at a time. Any other use does not comply with the intended purpose.
- The bike may only be used by one person with a maximum weight of up to 125 kg.
- Use of the bike is only recommended for those who have the right physical and psychological constitution to be able to cope with public traffic.
- As a beginner, you should ride particularly carefully at the start.
- Familiarise yourself with the general traffic laws! In order to achieve a stable ride when travelling straight ahead, wherever possible, you should avoid hasty steering movements when travelling at speed or going round corners.
- When you are riding and braking, always hold the handles firmly with both hands, as otherwise, there is a high risk of an accident through steering error.

- The bike is a mechanical traction device. The speed must be adjusted to suit the technical capabilities of the bike and the terrain.
- At night and dawn/dusk, always use the lights. When you are travelling long distances, we recommend that you take a spare battery with you.
- Only ride on asphalt paths, or firm, smooth surfaces. When riding over curbs, steps as well as potholes, there is an increased risk of tipping over.
- Adjust your method of driving to your capabilities. Ride slowly and carefully, particularly when riding round corners.
- Do not put your hands in the chain area.
- Your bike should be serviced by your dealer at regular intervals.
- Be careful with fire, in particular with burning cigarettes. Seat and back slings can be set alight.
- With all moving parts, there is always a risk that fingers could become trapped. Please always handle carefully.
- If the hand bike is subject to direct sunlight / low temperatures for long periods, parts of your hand bike could become very hot (>41°C) or very cold (<0°).
- Always make sure that the quick tensioner on the rear wheels is set properly.
- The hand bike should not be used in heavy rain, in snow or on slippery or damaged surfaces. In addition, you should avoid any dangerous areas.
- Make sure that your tyres have sufficient tread!
- Please note that you are subject to any and all traffic laws when riding in public traffic.
- Never use the leg support or the impact protection to lift the hand bike
- Familiarise yourself with the braking distances from various speeds.
- You should only use product combinations which have been approved by Sunrise Medical.

The hand bikes shown and described in this user guide may not correspond in every detail exactly to your own model. However, all instructions are completely relevant, regardless of possible detail differences. The manufacturer reserves the right to alter without notice any weights, measurements or other technical data shown in this manual. All figures, measurements and capacities shown in this manual are approximate and do not constitute specifications.

## Before each use, you should carry out the following safety checks

- Before each journey, check the tyre pressure. If the air pressure is too low, there is an increased risk of tipping over when going round corners, it can lead to increased wear and tear and the rolling resistance will be increased.
- Avoid damage to the Bowden cables for the gearshift and braking units. There is a risk of kinking the cable when you transfer to/from the bike and if you operate the crank incorrectly.
- Check that the footrests are secure.
- Check that all components, in particular all screws, are secure.
- To avoid serious injury, you must secure your feet in the footrests, using the Velcro straps, before starting your journey.
- Always wear a safety helmet.
- You must always ride with a safety flag, as, due to the low height of the bike, you may be less visible to other traffic.

6) Consumable parts are normally excluded from the guarantee, except in the case that premature wear of the part is the direct result of a manufacturing fault. These items include amongst others upholstery, tyres, inner tubes and similar parts.

7) The above warranty conditions apply to all product parts for models purchased at full retail price.

8) Normally we do not accept responsibility if a repair or replacement of the hand bike is required for one of the following reasons:

a) The product or part has not been maintained or serviced in accordance with the manufacturer's recommendations as shown in the User Instructions and/or the Service Instructions. Accessories have been used which are not specified as original accessories.

b) The hand bike or part having been damaged by neglect, accident or improper use.

c) Alterations to the hand bike or parts, which are not in accordance with the manufacturer's specifications or the carrying out of repairs before informing the customer service agent.

## 2.0 Warranty

### Warranty

THIS DOES NOT AFFECT YOUR LEGAL RIGHTS IN ANY WAY.

### Warranty conditions

1) The repair or replacement is carried out by an authorised Sunrise Medical dealer/service agent.

-

2) To fulfil the guarantee conditions, should servicing need to be carried out on your hand bike under this agreement, contact the designated Sunrise Medical customer service agent immediately, with precise details on the type of difficulty. Should you be using the hand bike outside the area covered by the designated Sunrise Medical customer service agent, the work will be carried out under "guarantee conditions" by another agency as designated by the manufacturer.

3) Should any part of the hand bike require repair or replacement, as a result of a specific manufacturing or material defect, within 24 months (5 years on frame) from the date on which the possession of the product was transferred to the original purchaser, and subject to it remaining within that ownership, the part or parts will be repaired or replaced completely free of charge if returned to the authorised service agent.

Note: This guarantee cannot be transferred.

4) The guarantee also covers all repaired or replaced parts for the remaining period of the guarantee for the hand bike.

5) For spare parts which are fitted after the start of the original guarantee, we give a further 24-months guarantee.

### 3.0 Handling

#### Delivery

Your new Shark RS / RT is delivered fully assembled in a cardboard box. In order to avoid any damage during transportation, loose or parts to be fitted, such as the rear wheels are packed separately in the cardboard box.

#### Unpacking:

- Check the shipping cardboard box for external damage which may have been caused during transportation.
- Remove all packaging material.
- Take all bike parts (including wheels) carefully out of the box and check that they are complete.
- Check all parts for surface damage, scratches, cracks, dents, distortions and other defects.
- Should you find any damage, or if any parts are missing please let us know immediately.

#### Transporting the Shark RS / RT:

To keep the Shark RS / RT as compact as possible, the rear wheels are removed before delivery and packed separately in the box.

The wheels are mounted to the rear axle by simply pushing them on or by screwing them on, depending on your choice of axles.

#### CAUTION!

**Before each use, you must check that the screw and plug connections are secure!**

### 4.0 Options

#### Brake

#### Brakes

The bike has two brakes which work independently from one another.

#### Parking brake (Fig. 4.1 + 4.2)

The parking brakes, mounted on the fork serve as wheel locks. They can also be applied with the aid of a locking lever. The locking lever must be locked into place for these to work properly. The bike should not roll when the parking brakes are applied. This is the only way to guarantee safe transfer to and from the bike.

#### Main brakes (Fig. 4.3)

The brakes allow you to safely and comfortably brake your hand bike.

#### CAUTION

**If the main brakes should fail, then the parking brakes may be used once for emergency braking. After that, the main brakes must be repaired.**

Fig. 4.1


Fig. 4.2



Fig. 4.3

**⚠ CAUTION!**

**Braking power will decrease with:**  
**Worn tyre tread, air pressure which is too low, soiled and wet tyres.**  
**But the ground surface can also play a significant role in delaying the braking effect. The braking effect will be significantly reduced if the streets and paths are wet and dirty, as well as on sandy or stony ground. Dirty and wet rims and brake pads also reduce the braking effect when using rim brakes.**  
**For this reason, you should always adapt your method of riding to the environmental conditions in order to avoid serious accidents.**  
**Always check the efficiency of the brakes before starting your journey and ensure that there is sufficient brake function. The brakes should be checked regularly by your dealer and adjusted if required.**

**Back pedal brake (Fig. 4.4 - 4.6)**

The back pedal brake option allows you to brake by pedalling backwards, when used in combination with chain gearshift (Fig. 4.5). This option is required particularly for users with restricted or no finger function.  
The braking system can only be uncoupled for manoeuvring (Fig. 4.6).  
When riding or moving forward, use the rocker arm to put the system into an operational position.

**⚠ CAUTION!**

Riding with the back pedalling brake deactivated is not permitted and can lead to an increased risk of accidents.

Fig. 4.4



Fig. 4.5

Fig. 4.6

Fig. 4.7

**Footrest**

**Adjusting the footrest (Fig. 4.7 - 4.9)**

The footrests can be adjusted to suit the leg length. The adjustment can be made after releasing the screws (1). When you adjusting the length, you must make sure that the legs are slightly angled and that they touch the bar at the front. With the bike at a standstill, you should check that when turning at the maximum angle, the increased knee bend does not lead to contact with the drive crank.

There must be a gap of at least 5 cm between the heel and the floor. You will need this leeway when travelling round corners. When finished, make sure that all screws have been properly tightened (see the page on torque).




The footrest tubes for longer legs are always delivered at the maximum length. The customer can then shorten them to the required length.

After shortening the tube, you must always ensure that the footrest tube is pushed into the fork far enough so that both fixing screws touch the footrest tube.

### ⚠ CAUTION!

**Never use the leg support to lift the hand bike. Before use, you must secure your feet into the footrests.**

**If your feet accidentally come out when you are riding, this can lead to serious injuries.**

### RT gearshift

**Chain gearshift/hub gearshift (Fig. 4.10 - 4.11)**

Gear shift is carried out using the Grip Shift switch, which can be mounted either on the frame, the fork or the handle.

Shifting the upper chain rings on the crank is carried out using a lever which is mounted on the crank housing.

The chain gearshift enables you to easily keep moving, regardless of other road conditions.

**A description of the gear change unit is included as part of the delivery. If the gearshift suddenly does not operate smoothly, please contact your dealer!**

### RS gearshift (Fig. 4.12)

This is carried out using the trigger switch on the right handle or the Dual Control switch on the right-hand side of the fork.

## Backrest adjustment

### Backrest adjustment (Fig. 4.13 - 4.14)

During the seating test ensure that the cranks turn freely above the thighs, even when going round a corner. (Test backwards in freewheel, adjust seat cushion thickness if necessary ). Set the gap from the crank so that, when turning the crank, your arms stretch to a comfortable distance and can bend easily. The very best drive characteristics are achieved if you use the adjustment options for the backrest and the crank height, in order to establish the optimum position for you.

Adjusting the backrest is carried out lengthways by releasing and moving the fixing clamp. Adjusting the backrest angle is carried out by adjusting the height of the central tube behind the backrest.

### ⚠ CAUTION!

**After making the adjustment, re-tighten all screws. (see the page on torque)**


Fig. 4.14



Fig. 4.15

## Height and depth adjustable fork

### Diagonal adjustment (Fig. 4.15)

Undo the socket head screws (1) on the fork tube and the hex nuts (2) on the clamp. Push the crank bracket forwards or backwards. When the crank is in the required position, the chain tension will not be correct, so to correct the tension, chain links MUST be added or removed. This procedure should be carried out by a qualified technical engineer. Finally, re-tighten the socket head screws which fix the crank bracket to the fork.

### Vertical adjustment

Undo the socket head screws (3) which fix the crank bracket to the fork stem, but do not remove them. Move the hand crank unit up or down in the fork collar, until you have achieved the optimum position for you. Tighten the socket head screws (3) to fix the hand crank in the required position. Re-tension the chain if necessary.

## Rear axle

### Rear axle (Fig. 4.16 - 4.17)

There are 3 different axles in the range, which can each be screwed to the frame using the axle clamp.

Swap / change the axle:

Undo the 3 fixing screws on the left and right axle clamp. Remove the axle and exchange it.

To

complete the operation, place the half shell of the axle clamp back onto the axle and fix it with the three fixing screws.

Tighten the 3 fixing screws alternately until all 3 screws are fixed to a torque of 5 Nm.

**Note: When changing from aluminium to carbon tubes, the eccentric sleeve must also be changed.**

Fig. 4.16



Fig. 4.17

## 5.0 Maintenance and care (checks to be carried out before use)

### When

### What

### Comment

Every 4 weeks  
(depending on  
frequency of use)

Grease the cables where  
they leave the switches and  
the brake levers.

Check all screw  
connections are secure!  
Oil the drive chain!

Carry out yourself or with an attendant.  
Before oiling, please remove all old oil  
remnants from the components.  
Make sure that the surplus oil does  
not get the surrounding area  
(e.g. your clothes) dirty.  
Carry out yourself or with an  
attendant.  
Spray lightly with Teflon spray  
Be careful not to get it on your clothes!

Every 2-3 months  
(depending on fre-  
quency of use)

Check:  
• Cleanliness  
• General condition

Carry out yourself or with an attendant.  
If there are any cracks or general wear  
and tear, then change the rubbers or  
visit an authorised service centre  
for a refurbishment.

### Service:

If you have any questions or require any help, please contact your dealer, who will be pleased to advise you with regards to servicing and repairs.

**Depending on how much you have used your Shark RS / RT, we recommend that you take it to your dealer at least every 6 months, to have your Shark RS / RT inspected by trained personnel.**

## 6.0 Tyre maintenance

### Tyres

The Shark RS / RT is equipped with special 26" tyres on the rear wheels and a 26" drive wheel. Make sure that they are always at the recommended tyre pressure, as this will significantly affect the drive characteristics of your Shark RS / RT. If the tyre pressure is too low, movement and manoeuvrability is made more difficult due to the increased rolling resistance and this also increases the amount of effort required.

The optimum tyre pressure is printed on each tyre.

### Fitting and repairing tyres

Before fitting a new inner tube:

- Inspect the base of the rim and the tyre inner wall for foreign bodies and if necessary clean thoroughly.
- Check the condition of the rim tape, particularly at the valve hole. Only well fitting, strong rim tape offers long-term protection against inner tube damage through burrs and sharp spoke heads.

### Inflation:

Check all round on both sides, to make sure that the inner tube is not trapped between the tyre bead and the rim. Push the valve back gently and pull it back out, so that the inner tube achieves an optimum position around the valve.

### Tube tyres

As an option you can have tube tyres with full carbon rims.

With these, please note:

Tube tyres must be glued to the rim by a specialist. Should you have a puncture, you must not continue to ride (the immediate consequences are an increased risk of accident and the rim will be destroyed).

**CAUTION!**

**After fitting, check the tyre pressure.**

**It is essential that you maintain regulation tyre pressure and keep the tyres in perfect condition, in order to ensure your personal safety and keep your Shark RS / RT running properly.**

## 7.0.0 Disposal / Recycling of Materials

**NOTE:** If the hand bike has been made available to you as part of a charity or medical loans scheme, then it may not belong to you. If it is no longer required, then follow the instructions to return it as given by the organisation that made the hand bike available to you.

In the following section, there is a description of the materials used on the hand bike, in view of the disposal or recycling of the hand bike and its packaging. There may also be special local regulations in force with regard to disposal or recycling, these must be taken into account when disposing of your hand bike. (This can include the cleaning or decontamination of the hand bike prior to disposal).

**Aluminium:** Castor Forks, wheels, footrest,

**Steel:** Fixing points,

**Carbon:** Rear axle

**Plastic:** Handles, tube stoppers, armpads and tyres

**Packaging:** Made of soft polyethylene, cardboard

**Upholstery:** Woven polyester with PVC coatings and expanded combustion modified foam.

Disposal or recycling should be done through a licensed agent or authorised place of disposal. Alternatively your hand bike may be returned to your dealer for disposal.

### Hygiene measures when being re-used:

Prior to the wheelchair being re-used, it must be carefully prepared. All surfaces which come into contact with the user must be treated with a disinfection spray.

To do this, you must use a disinfectant as authorised/recommended in your country, for rapid alcohol-based disinfection for medical products and medical devices, which must be disinfected quickly. Please be aware of the manufacturer's instructions for the disinfectant you are using.

In general, a complete disinfection cannot be guaranteed on seams. We therefore recommend that you dispose of seat and back slings to avoid micro-bacterial contamination with active agents according to your local infection protection law.

## 8.0 Nameplate

### Nameplate


The nameplate is located on either the rear part of the frame or on the frame underneath the seat surface, as well as on a label in the owner's manual. The nameplate indicates the exact model designation and other technical specifications. Please provide the following pieces of information whenever you have to order replacement parts or to file a claim:


As we are committed to continually improving the construction of our wheelchairs and hand bikes, product specifications may vary slightly from the examples illustrated. All weight/dimensions and performance data is approximate and is provided solely for guidance.

All hand bikes must be used in accordance with the manufacturer's guidelines.

 Sunrise Medical GmbH  
Kahlbachring 2-4  
69254 Malsch/Heidelberg  
Deutschland  
Tel.: +49 (0) 7253/980-0  
Fax: +49 (0) 7253/980-222  
kundenservice@sunrisemedical.de  
www.SunriseMedical.de


<b>TYPE:</b>	Product Name/SKU Number.
	Seat width.
	Load Maximum.
	UKCA Mark.
	CE Mark.
	Consult instructions for use.
	Date of manufacture.
	Serial number.
	This symbol means Medical Device.
	Manufacturer's address.
	Importer's address
	UK Responsible Person
	Swiss Representative's address

## 9.0 Torque

Steering head bearings:	Handtight
Castor fork clamps:	7 Nm
Screw axle on drive wheel:	10 Nm
Footrest clamp screws:	7 Nm
Crank clamp screws:	7 Nm
Front backrest clamp screws:	7 Nm
Rear backrest clamp screws:	7 Nm

**Torque for M6 screws is 7 Nm, unless otherwise specified.**

Sunrise Medical S.r.l.  
Via Riva, 20 – Montale  
29122 Piacenza  
Italia  
Tel.: +39 0523 573111  
Fax: +39 0523 570060  
www.SunriseMedical.it

Sunrise Medical AG  
Erlenauweg 17  
CH-3110 Münsingen  
Schweiz/Suisse/Svizzera  
Fon +41 (0)31 958 3838  
Fax +41 (0)31 958 3848  
www.SunriseMedical.ch

Sunrise Medical AS  
Delitoppen 3  
1540 Vestby  
Norge  
Telefon: +47 66 96 38 00  
post@sunrisemedical.no  
www.SunriseMedical.no

Sunrise Medical AB  
Neogatan 5  
431 53 Mölndal  
Sweden  
Tel.: +46 (0)31 748 37 00  
post@sunrisemedical.se  
www.SunriseMedical.se

MEDICCO s.r.o.  
H – Park, Heršpická 1013/11d,  
639 00 Brno  
Czech Republic  
Tel.: (+420) 547 250 955  
Fax: (+420) 547 250 956  
www.medicco.cz  
info@medicco.cz  
Bezplatná linka 800 900 809

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Sunrise Medical Australia  
11 Daniel Street  
Wetherill Park NSW 2164  
Australia  
Ph: +61 2 9678 6600  
Email: enquiries@sunrisemedical.com.au  
www.SunriseMedical.com.au

Sunrise Medical  
North American Headquarters  
2842 Business Park Avenue  
Fresno, CA, 93727, USA  
(800) 333-4000  
(800) 300-7502  
www.SunriseMedical.com

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Sunrise Medical B.V.
Sunrise Medical HCM B.V.
Sunrise Medical SAS