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

© 2016 Sunrise Medical

All rights reserved.

The information provided herein may not be reproduced and/or published in any form, by print, photoprint, microfilm or any other means whatsoever (electronically or mechanically) without the prior written authorization of Sunrise Medical.

The information provided is based on general data concerning the construction known at the time of the publication of this manual. Sunrise Medical executes a policy of continuous improvement and reserves the right to changes and modifications.

The information provided is valid for the product in its standard version. Sunrise Medical can therefore not be held liable for any damage resulting from specifications of the product deviating from the standard configuration. Illustrations contained in the manual may deviate from the configuration of your product.

The information made available has been prepared with all possible diligence, but Sunrise Medical cannot be held liable for any errors contained in the information or the consequences thereof.

Sunrise Medical accepts no liability for loss resulting from work executed by third parties.

Names, trade names, trademarks etc. used by Sunrise Medical may not, as per the legislation concerning the protection of trade names, be considered as being available.

Table of contents

1	Introduction.....	20
1.1	This user manual	20
1.2	For further information	20
1.3	Symbols used in this manual	20
2.	The controller	21
2.1	DX2-REM420 / DX2-REM421 controller	21
3	Driving the wheelchair with the controller.....	24
3.1	Switching the controller on or off	24
3.2	Driving the wheelchair	24
3.3	Speed	24
3.4	Driving Profile	24
4	Operating the electrical adjustment options	25
4.1	Selecting the desired adjustment option	25
5	Lights.....	26
6.	Troubleshooting.....	27
7.	Locking the controller	28
8	Technical specifications.....	29
9.	Technical diagrams.....	30
9.1	Technical diagram	30
9.2	Technical diagram battery charger	31

1 Introduction

1.1 This user manual

This user manual will help you to use and maintain the controller of your power wheelchair safely. This user manual is a supplement to Sunrise Medical's general wheelchair user manual. When necessary this user manual refers to other manuals as shown below:

-  Wheelchair: Refers to the general wheelchair user manual.
-  Battery charger: Refers to the user manual for the battery charger.

Read this user manual and the other user manuals referred to carefully before using the product. If one of the user manuals was not included with your wheelchair, please contact your dealer immediately. In addition to this user manual, there is also a service manual for qualified specialists.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu or alternatively is available on request in large text.

1.2 For further information

Please contact your local, authorised Sunrise Medical dealer if you have any questions regarding the use, maintenance or safety of your wheelchair. In case there is no authorised dealer in your area or you have any questions, contact Sunrise Medical either in writing or by telephone.

For information about product safety notices and product recalls, go to www.SunriseMedical.eu

1.3 Symbols used in this manual

Note!

Pointing out possible problems to the user.

⚠ Caution!

Advice for the user to prevent damage to the product.

⚠ Warning!

Warnings for the user to prevent personal injury.

Not following these instructions may result in physical injury, damage to the product or damage to the environment!

2. The controller

2.1 DX2-REM420 / DX2-REM421 controller



A controller will usually have three basic functions:


- Driving and steering a wheelchair
- Operating electrical seat adjustments
- Charging the wheelchairs batteries

There are many different control systems for wheelchairs on the market. If the controller on your wheelchair does not resemble the one in figure 1 or 2 contact your dealer.

1. DX2-REM420
2. DX2-REM421

DX2-REM420 / DX2-REM421 by Dynamic Controls is a collective name for the entire control system of the wheelchair. The wheelchair is operated by means of a controller, which includes the following components:

From the left to the right ; DX2-REM420 and DX2-REM421

Part	Function
A. Joystick	In the 'drive' mode: driving and steering In the 'adjustment options' mode: • Left/right to select the adjustment options • Front/back to control the adjustment mode
B. On/off switch	Switching the controller on or off
C. Horn	Warning signal with sound
D.1 Battery indicator	Displays the power level of the battery
D.2 Maximum speed	Displays the maximum speed limit as set by the user
D.3 Light and/or indicator	Displays the status of lights and/or indicator
D.4 Fault indicator	Displays fault code in case of system error
D.5 Attendant profile LED	Displays if DK-ACU attendant controle module is set to 'attendant'
E. Adjustment option	To select the electrical seating adjustment 1, 2, 3 or 4. The actuator LEDs indicates the currently selected seating function
F. 1 Speed regulator	Reduce driving speed (slower)
F. 2 Speed regulator	Increase driving speed (faster)
G. Charge connector	Input for the battery charger
H. Charge plug for the battery charger	Connector for the battery charger
I. Lights button	Switching the lights on or off

Part	Function
J. Hazard lights	Warning signal with lights
K.1 Direction indicator left	Switches the left direction indicator on or off
K.2 Direction indicator right	Switches the right direction indicator on or off
L. Light sensor	Dims the display when it is dark
M Lock LED	Displays if the DX system is locked

3 Driving the wheelchair with the controller

3.1 Switching the controller on or off

To be able to drive or operate the electric adjustment options of the wheelchair, the controller must be switched on. Press the on/off button (B in figure 3).

3.2 Driving the wheelchair

Driving an electric wheelchair is done by operating a joystick. Move the joystick forwards and the wheelchair will also move forwards. Steer left and right and the wheelchair will turn.

3.3 Speed

The maximum speed can be controlled by the speed regulator on the controller (F1 and F2 in figure 3). De maximum speed is displayed by the LED lights (D2 figure 3). Speed can be controlled with the joystick while driving. If the joystick is moved a little, the wheelchair will move more slowly.

3.4 Driving Profile

This controller can also be set to make the wheelchair suitable for different driving profiles or environments. For example, selecting the profile for indoors will ensure that the wheelchair will react more ‘calmly’. Once outside, the profile can be adjusted to a more ‘robust’ setting. The currently selected drive profile is shown on the drive profile LED’s (D2 in figure 3). DX2-REM420 and 421 allows you to adjust the maximum speed within different profiles or environments

Selecting a profile

To switch to a different profile, just press the ‘speed regulator’ button (F1 and F2 in figure 3). Press the ‘speed regulator’ button until the most appropriate profile has been selected.

‘+’ to increase the drive profile.

‘-’ to decrease the drive profile

LEDs	Meaning	4
	Drive Profile 0*. The powerchair does not drive.	
	Drive Profile 1.	
	Drive Profile 2.	
	Drive Profile 3.	
	Drive Profile 4.	
	Drive Profile 5.	
	Attendant Drive Profile (Drive Profile 6).	

4 Operating the electrical adjustment options

Not every wheelchair has been equipped with electric adjustment options. We make a distinction between four different adjustment options to the seating system:

- Tilt adjustments
- High/low adjustments
- Backrest adjustments
- Legrest adjustments

4.1 Selecting the desired adjustment option

1. Press the ‘adjustment’ button (E in figure 3) to select the ‘adjustment options’ mode. The ‘adjustment option’ mode starts with the seating function that was last selected. The LED light indicates what electrical adjustment option is enabled (1, 2, 3 or 4).
2. Use the ‘adjustment option’ button or joystick to switch between the electrical adjustment options (‘+’ / right for the next seating function and ‘-’ / left for the previous seating function).
3. Moving the joystick forwards and/or backwards will activate the selected electric adjustment option (see table). Move the joystick forwards or backwards until the desired adjustment option has been attained.

	Move joystick backwards	Move joystick forwards
Tilt adjustment	The entire chair will tilt backwards	The entire chair will tilt forwards
High/low adjustments	The entire chair will be raised	The entire chair will be lowered
Backrest adjustments	The backrest will tilt backwards	The backrest will tilt forwards
Legrest adjustments	The legrest angle will increase, the footplate will be raised	The legrest angle will decrease, the footplate will be lowered

4. To return to the ‘drive’ mode: Press the ‘speed regulator’ button (F in figure 3).

Note!
If you adjust the seat height by using the high/low option, the speed will be reduced due to safety reasons.

5 Lights

Not every wheelchair has been equipped with lights. You can activate the following functions only by using the DX2-REM420:

- Lights (I in figure 3)
- Hazard Lights (J in figure 3)
- Direction indicators (K1 & K2 in figure 3)

6. Troubleshooting

If the wheelchair will not function while the batteries are fully charged, check the following points before consulting your dealer:

- Switch the controller off and then switch it on again. Check to see if the malfunction has been solved.
- Check if the free wheel switch was switched to Drive.
- Check if the joystick was in the 0 position when the controller was switched on. In other words, the joystick must not be moved when the controller is being switched on or off.

Malfunctions list

An extended malfunction list can be found in the appendix of the service manual (for qualified specialists only). The service manual can be found on www.SunriseMedical.eu



7. Locking the controller

To lock the controller (figure 5):

- Press the on/off button for 4 seconds when the DX system is turned on.
- The DX system turns off immediately when the on/off button is pushed.
- After 4 seconds the 'maximum speed' (D2 in figure 3) LEDs and the lock LED turn on for 1 second to indicate that the DX system is now locked.

To unlock the controller (figure 6):

- Press the on/off button.
- The lock LED will flash for 10 seconds to indicate that the DX System is locked.
- Press the horn button twice within 10 seconds. The DX System will turn on normally.
- If the horn button is not pressed within 10 seconds, the controller will turn off again.

8 Technical specifications

Operating voltage range:	18V - 32V DC (nom. 24V)
Charger rating:	max 12A RMS Continuous limited by DXBUS Rating
Quiescent Current:	< 0.25mA Off, typically 200mA On
Tube mount diameter	Min 15mm, Nom 19mm, Max 22mm
Force required to operate joystick DX-REM42x	Less than 2.2N
Force required to operate buttons DX-REM42x	Less than 5N
Protection rating:	IPx4
Operating Temperature Range:	-25°C to +50°C
Storage Temperature Range:	-40°C to +70°C
Operating Humidity Range:	0 to 95% RH



EMC tested on sample wheelchair:


	ANSI/RESNA WC/Vol:2 - 1998 Sec 21 EN12184: 1999 Sec 9.8.1
--	--

9. Technical diagrams

9.1 Technical diagram


The technical diagrams can also be found on the cover of the specific electronic component.

9.2 Technical diagram battery charger

The controller's standard configuration includes a '3-pin connection'. Ensure that the battery charger is properly connected so that the 'negative pole' and the 'inhibit' are connected, enabling the system to prevent the wheelchair from moving when the battery is being charged.

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Sunrise Medical AB