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
8 Technische specificaties

Supply Voltage:	24 Vdc
Operating Voltage:	16Vdc to 35Vdc
Peak Voltage:	35Vdc
Reverse Battery Current:	40Vdc
PWM Frequency:	20kHz $\pm$ 0.5%
Brake Voltage:	12/24Vdc
Brake Current:	200 μ A min. 1A max.
Charger Connector:	Use only Neutrik NC3MX
Battery Charging Current:	12Arms max.
Maximum Drive Current:	R-net 60 30A R-net 80 80A R-net 120 120A
Indicator Outputs:	45W per side
Lighting Outputs:	21W per side
Brake Light Output:	42W total
Actuator Current:	15A max at reduced speed. 12A max at full speed.
Moisture Resistance:	Electronics to IPX4
Operating Temperature:	Non LCD Modules -25°C to +50°C Modules with LCD Screens -10°C to +50°C
Storage Temperature:	Non LCD Modules -40°C to +65°C Modules with LCD Screens -20°C to +65°C

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



© 2016 Sunrise Medical

All rights reserved.

The information provided herein may not be reproduced and/or published in any form, by print, photoprint, microfilm or any other means whatsoever (electronically or mechanically) without the prior written authorization of Sunrise Medical.

The information provided is based on general data concerning the construction known at the time of the publication of this manual. Sunrise Medical executes a policy of continuous improvement and reserves the right to changes and modifications.

The information provided is valid for the product in its standard version. Sunrise Medical can therefore not be held liable for any damage resulting from specifications of the product deviating from the standard configuration. Illustrations contained in the manual may deviate from the configuration of your product.

The information made available has been prepared with all possible diligence, but Sunrise Medical cannot be held liable for any errors contained in the information or the consequences thereof.

Sunrise Medical accepts no liability for loss resulting from work executed by third parties.

Names, trade names, trademarks etc. used by Sunrise Medical may not, as per the legislation concerning the protection of trade names, be considered as being available.

Table of contents

1 Introduction.....22

1.1 This user manual 22

1.2 For further information 22

1.3 Symbols used in this manual 22

2 The controller23

2.1 R-net controller 23

2.2 Direct access 23

2.3 R-net attendant steering 23

2.4 Main components 24

2.5 Display screen 25

3 Driving the wheelchair26

3.1 Switching the controller on or off 26

3.2 Driving the wheelchair 26

3.3 Adjusting the push bar height 26

3.4 Speed 27

3.5 Driving Profile 27

3.6 Selecting a profile 27

4 The electrical adjustment options.....28

4.1 The controller and the attendant steering 28

4.2 The Direct access 29

4.3 Mounting of the Direct access 30

5 Lights.....31

6 Troubleshooting.....32

7 Locking the controller33

7.1 Keypad Locking 33

7.2 Key Locking (optional) 33

8 Technical specifications.....34

9 Technical diagrams.....35

9.1 Technical diagram 35

9.2 Technical diagram battery charger 36

1 Introduction

1.1 This user manual

This user manual will help you to use and maintain the controller of your power wheelchair safely. This user manual is a supplement to Sunrise Medical's general wheelchair user manual. When necessary this user manual refers to other manuals as shown below:

-  Wheelchair: Refers to the general wheelchair user manual.
-  Battery charger: Refers to the user manual for the battery charger.

Read this user manual and the other user manuals referred to carefully before using the product. If one of the user manuals was not included with your wheelchair, please contact your dealer immediately. In addition to this user manual, there is also a service manual for qualified specialists.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu or alternatively is available on request in large text.

1.2 For further information

Please contact your local, authorised Sunrise Medical dealer if you have any questions regarding the use, maintenance or safety of your wheelchair. In case there is no authorised dealer in your area or you have any questions, contact Sunrise Medical either in writing or by telephone.

For information about product safety notices and product recalls, go to www.SunriseMedical.eu

1.3 Symbols used in this manual

Note!

Pointing out possible problems to the user.

⚠ Caution!

Advice for the user to prevent damage to the product.

⚠ Warning!

Warnings for the user to prevent personal injury.

Not following these instructions may result in physical injury, damage to the product or damage to the environment!

2 The controller

2.1 R-net controller

A controller will usually have three basic functions:

- Driving and steering a wheelchair
- Operating electrical seat adjustments
- Charging the wheelchairs batteries

There are many different control systems for wheelchairs on the market. If the controller on your wheelchair does not resemble the one in figure 1, contact your dealer.


R-net by PG Drives Technologies is a collective name for the entire control system of the wheelchair. The wheelchair is operated by means of a controller. Main components of the controller can be found in chapter 2.4.


2.2 Direct access

Direct access is an extra keypad that can be connected to the controller (figure 2). This gives the user the possibility to operate the adjustments without scrolling through a menu. There are different ways of mounting the direct access to suit the user, see chapter 3.

2.3 R-net attendant steering

It is possible to connect an extra controller on the back (figure 3). This so the attendant can control the wheelchair and operate electrical options.

2.4 Main components

Part	Function
A. Joystick	In the 'drive' mode: driving and steering In the 'adjustment options' mode:: • Left/right to select the adjustment options • Front/back to select the adjustment options
B. On/off switch	Switching the controller on or off
C. Horn	Warning signal with sound
D. LCD colour display screen	Display and feedback
E. "Mode" button	Changing between the 'driving' and the 'adjustment options' mode
F1: Speed regulator	Reduce driving speed (slower)
F2: Speed regulator	Increase driving speed (faster)
G. Charge connector	Input for the battery charger
H. Charge plug for the battery charger	Connector for the battery charger
I. Lights button	Switching the lights on or off
J. Hazard lights	Warning signal with lights
K1: Direction indicator left	Switches the left direction indicator on or off
K2: Direction indicator right	Switches the right direction indicator on or off
L. "Profile" button	Select driving profile
M. Direct access	Buttons to directly control electrical adjustments (see 4.2)
N. Maximum speed	Displays the maximum speed limit
O. User/attendant switch	Switching between user and attendant controller
P. User/attendant indicator	Displays the active controller

2.5 Display screen

When the R-net is powered-up the screen as in figure 5 will appear.

1. Battery indicator
Displays the power level of the battery
2. Speed indicator
Displays a graph of the actual speed
3. Kilometre (Mileage) counter
Displays the total distance driven
4. Kilometre (Mileage) an hour counter
Displays the actual speed
5. Maximum speed
Displays the maximum speed limit as set by the user
6. Current profile
Displays the selected profile
7. Clock
Displays the time
8. Adjustment option (figure 6)
Displays the selected adjustment options

Display settings

The settings menu (figure 7) allows you to adjust the display to suit your preferences: clock adjustment, clock display format, odometer, brightness of the backlight, background colour.


The settings menu is accessed by depressing simultaneously the 'Speed Down' and 'Speed Up' buttons (F1 en F2 in figure 4) for a few seconds.

In the menu move the joystick:

- front/back to go through the setting options.
- left/right to select a setting option.

Once you selected a setting option, use front/back and left/right to set the desired values/options.

Leave the menu by selecting 'Exit'.

3 Driving the wheelchair

3.1 Switching the controller on or off

To be able to drive or operate the electronic adjustment options of the wheelchair, the controller must be switched on. Press the on/off button (B in figure 4).

Note!

The joystick should be released (neutral position) when switching on the controller otherwise there will be an error code

If the attendant steering needs to be activated, the user/attendant switch needs to be pressed (O in figure 4).

3.2 Driving the wheelchair

Driving an electric wheelchair can be done by operating the controller or the attendant steering.

Driving with the controller or the attendant steering on push bar


Move the joystick forwards and the wheelchair will also move forwards. Steer left and right and the wheelchair will turn.

Driving with the attendant steering mounted in the Z-steering mechanism

To operate the controller in the Z-steering mechanism the push handles needs to be squeezed to the metal bar below. The chair will drive in the forward direction. When moving the bars to the left, the chair will turn right and when moving them right, the chair will turn left. Lift the push handles if you want to drive backwards.

3.3 Adjusting the push bar height

The height of the push bar with the Z-steering mechanism is adjustable in 6 steps and can easily be set. Pull A (figure 8) and set the height by releasing the strap on the desired height. You should hear a click when setting the push bar height. This indicates that the bar has locked into place. It is also possible to flip down the total Z-steering.

3.4 Speed

The maximum speed can be controlled by the speed regulator on the controller (F1 en F2 in figure 4). The speed bar on the display screen will display the maximum speed. Speed can be controlled with the joystick while driving. If the joystick is moved a little, the wheelchair will move more slowly.

3.5 Driving Profile

This controller can also be set to make the wheelchair suitable for different driving profiles or environments (L in figure 4). For example, selecting the profile for indoors will ensure that the wheelchair will react more 'calmly'. Once outside, the profile can be adjusted to a more 'robust' setting. The name and number of the current profile is displayed below the speed bar. R-net allows you to adjust the maximum speed within different profiles or environments.

3.6 Selecting a profile

To switch to a different profile, just press the 'profile' button (L in figure 4). Press the 'profile' button until the most appropriate profile has been selected.

4 The electrical adjustment options

Not every wheelchair has been equipped with electronic adjustment options. We make a distinction between four different adjustment options to the seating system (figure 9):

- Tilt adjustments
- High/low adjustments
- Backrest adjustments
- Legrest adjustments

Note

If you adjust the seat height by using the high/low option or use the electrical tilt adjustment 0 - 45°, the speed will be reduced for your safety.




Operating electrical adjustment options can be done by the controller, the attendant steering or the Direct access:

4.1 The controller and the attendant steering

1. Switch on the controller (A in figure 10)
2. To activate the attendant steering the user/ attendant switch needs to be pressed (A in figure 11 + 12)
3. Press on the 'mode' button (B) to select the 'adjustment options' mode. The controller will now be in the 'adjustment options' mode. It is possible to follow the status on the display screen.

The joystick or the push handles (C) are used to select and operate the electric adjustment options.

4. Move the joystick / handles to the left or the right to select the desired adjustment option. The actuator LED's indicates the currently selected seating function.
5. Moving the joystick forwards / backwards or the handles up / down will activate the selected electric adjustment option (see table on next page). Move the joystick / handles until the desired adjustment option has been attained.

Adjustment	Move joystick backwards / Pushing the push handles down
Tilt	The entire chair will tilt backwards
High/low	The entire chair will be raised
Backrest	The backrest will tilt backwards
Legrest	The legrest angle will increase, the footplate will be raised



Adjustment	Move joystick forwards / Lifting the push handles up
Tilt	The entire chair will tilt forwards
High/low	The entire chair will be lowered
Backrest	The backrest will tilt forwards
Legrest	The legrest angle will decrease, the footplate will be lowered

6. Press on the ‘mode’ button to return to the ‘drive’ mode.

4.2 The Direct access

With the Direct access option the electrical adjustments can be operated directly, without scrolling through a menu (figure 13).

7. Switch on the controller (A in figure 10)
8. The upper buttons will operate the specific adjustment up / forwards. The lower buttons will operate the adjustment down / backwards (see table below). Push the button on the Direct access key pad to operate the electrical adjustment option until the desired adjustment position has been attained.


Adjustment	Upper buttons	Lower buttons
1 Tilt	The entire chair will tilt forwards	The entire chair will tilt backwards
2 High/low	The entire chair will be raised	The entire chair will be lowered
3 Backrest	The backrest will tilt forwards	The backrest will tilt backwards
4a Central legrest	The central legrest angle will increase, the footplate will be raised	The central legrest angle will decrease, the footplate will be lowered
4b Left legrest	The left legrest angle will increase, the footplate will be raised	The left legrest angle will decrease, the footplate will be lowered
5 Right legrest	The right legrest angle will increase, the footplate will be raised	The right legrest angle will decrease, the footplate will be lowered

4.3 Mounting of the Direct access

The Direct access should be mounted to the armrest in the way that suits best to the user:

- A. On joystick side, in front of the armrest
- B. On the opposite joystick side, below the armrest, left
- C. On the opposite joystick side, below the armrest, right
- D. On the opposite joystick side, below the armrest, front

5 Lights

With the R-net it is possible to control the following lights:

- Lights (I in figure 4)
- Hazard Lights (J in figure 4)
- Direction indicators (K1 & K2 in figure 4)

6 Troubleshooting

If the wheelchair will not function while the batteries are fully charged, check the following points before consulting your dealer:

- Switch the controller off and then switch it on again. Check to see if the malfunction has been solved.
- Check if the free wheel switch was switched to Drive.
- Check if the joystick was in the 0 position when the controller was switched on. In other words, the joystick must not be moved when the controller is being switched on or off.

Malfunctions list

An extended malfunction list can be found in the appendix of the service manual (for qualified specialists only). The service manual can be found on www.SunriseMedical.eu

7 Locking the controller

The control system can be locked in one of two ways. Either using a button sequence on the keypad or with a physical key. How the control system is locked depends on how the wheelchair manufacturer has programmed the system.

7.1 Keypad Locking

To lock the wheelchair using the keypad;

- While the control system is switched on, depress and hold the on/off button.
- After 1 second the control system will beep. Now release the on/off button
- Deflect the joystick forwards until the control system beeps.
- Deflect the joystick in reverse until the control system beeps.
- Release the joystick, there will be a long beep.
- The wheelchair is now locked.

To unlock the wheelchair using the keypad;

- If the control system has switched off, press the on/off button.
- Deflect the joystick forwards until the control system beeps.
- Deflect the joystick in reverse until the control system beeps.
- Release the joystick, there will be a long beep.
- The wheelchair is now unlocked.


7.2 Key Locking (optional)

To lock the controller with the special key (H1 in figure 15):

- Insert and remove the special key into the Charger Socket on the Controller (G in figure 15). The controller is now locked and the following screen will be displayed.

To unlock the controller with the special key:

- If the control system has switched off, press the On/off (B in figure 4).
- Insert and remove the special key into the Charger Socket on the Controller (G in figure 15). The controller is now unlocked.

8 Technical specifications

Supply Voltage	24 Vdc
Operating Voltage	16Vdc to 35Vdc
Peak Voltage	35Vdc
Reverse Battery Current	40Vdc
PWM Frequency	20kHz $\pm$ 0.5%
Brake Voltage	12/24Vdc
Brake Current	200 μ A min. 1A max.
Charger Connector	Use only Neutrik NC3MX
Battery Charging Current	12Arms max.
Maximum Drive Current	R-net 60 30A R-net 80 80A R-net 120 120A
Indicator Outputs	45W per side
Lighting Outputs	21W per side
Brake Light Output	42W total
Actuator Current	15A max at reduced speed. 12A max at full speed.
Moisture Resistance	Electronics to IPX4
Operating Temperature	Non LCD Modules -25°C to +50°C Modules with LCD Screens -10°C to +50°C
Storage Temperature	Non LCD Modules -40°C to +65°C Modules with LCD Screens -20°C to +65°C



EMC tested on sample wheelchair:


Applicable aspects	Applicable standards
Radiated emissions Electrostatic discharge immunity Radiated r.f. field immunity	EN12184: 2009 and underlying norms

9 Technical diagrams

9.1 Technical diagram


The technical diagrams can also be found on the cover of the specific electronic component.

9.2 Technical diagram battery charger

The controller's standard configuration includes a '3-pin connection'. Ensure that the battery charger is properly connected so that the 'negative pole' and the 'inhibit' are connected, enabling the system to prevent the wheelchair from moving when the battery is being charged.

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The logo features a stylized swoosh above the word 'QUICKIE' in a bold, sans-serif font, followed by a registered trademark symbol (®).
Sunrise Medical AB