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© 2016 Sunrise Medical

All rights reserved.

The information provided herein may not be reproduced and/or published in any form, by print, photoprint, microfilm or any other means whatsoever (electronically or mechanically) without the prior written authorization of Sunrise Medical.

The information provided is based on general data concerning the construction known at the time of the publication of this manual. Sunrise Medical executes a policy of continuous improvement and reserves the right to changes and modifications.

The information provided is valid for the product in its standard version. Sunrise Medical can therefore not be held liable for any damage resulting from specifications of the product deviating from the standard configuration. Illustrations contained in the manual may deviate from the configuration of your product.

The information made available has been prepared with all possible diligence, but Sunrise Medical cannot be held liable for any errors contained in the information or the consequences thereof.

Sunrise Medical accepts no liability for loss resulting from work executed by third parties.

Names, trade names, trademarks etc. used by Sunrise Medical may not, as per the legislation concerning the protection of trade names, be considered as being available.

Table of contents

1	Intro	oduction	19
	1.1	This user manual	19
	1.2	For further information	19
	1.3	Symbols used in this manual	19
2	The	controller interface	20
	2.1	Omni interface	20
	2.2	Main components	21
	2.3	Display screen	22
	2.4	User menu	22
3	Drivi	ing the wheelchair	23
	3.1	Switching the wheelchair on or off	23
	3.2	Driving the wheelchair	23
	3.3	Speed	23
	3.4	Driving profile	23
	3.5	Selecting a profile	23
4	Ope	rating electrical adjustments	24
	4.1	Selecting the desired adjustment option	24
5	Inpu	It devices	
	5.1	External joystick	25
	5.2	External joystick lite	25
	5.3	Mini flex switch joystick	25
	5.4	Micro pilot joystick	26
	5.5	Micro guide joystick	26
	5.6	Proximity head array	26
	5.7	Spec switch	27
6	Ligh	ts	28
7	Trou	Ibleshooting	

1 Introduction

1.1 This user manual

This user manual will help you to use and maintain the controller of your power wheelchair safely. This user manual is a supplement to Sunrise Medical HCM's general wheelchair user manual. When necessary this user manual refers to other manuals as shown below:

- Wheelchair: Refers to the general wheelchair user manual.
- Battery charger: Refers to the user manual for the battery charger.
- Controller: Refers to the Omni operating manual from PGDT.

Read this user manual and the other user manuals referred to carefully before using the product. If one of the user manuals was not included with your wheelchair, please contact your dealer immediately. In addition to this user manual, there is also a service manual for qualified specialists.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu or alternatively is available on request in large text.

1.2 For further information

Please contact your local, authorised Sunrise Medical dealer if you have any questions regarding the use, maintenance or safety of your wheelchair. In case there is no authorised dealer in your area or you have any questions, contact Sunrise Medical either in writing or by telephone.

For information about product safety notices and product recalls, go to www.SunriseMedical.eu

1.3 Symbols used in this manual

Note!

Pointing out possible problems to the user.

▲ Caution!

Advice for the user to prevent damage to the product.

▲ Warning!

Warnings for the user to prevent personal in jury.

Not following these instructions may result in physical injury, damage to the product or damage to the environment!

2 The controller interface

2.1 Omni interface

A controller will usually have three basic functions:

- Driving and steering the wheelchair
- Operating electrical adjustments
- Charging the wheelchairs batteries

There are many different control systems for wheelchairs on the market. If you don't have the controller interface on your wheelchair as in figure 1, contact your dealer.

To drive the electric wheelchair with Omni, a special type of input device for connection is needed. This can for example be an external joystick, a mini-joystick or a head array. Sunrise Medical offers several external input devices to connect to the Omni. More information about these input devices can be found in chapter 5.

Omni is a controller interface by PG Drives Technologies. More information about the Omni interface can be found in this manual or in the Omni operating manual from PGDT (Controller).

2.2 Main components

Part	Function
A. LCD colour display screen	Display and feedback
B. On/off switch	Switching the controller on or off. There is a facility for an optional, externally mounted switch which can be operated by the user ('External on/off input')
C. "Profile" button	Changing between the available drive profiles
D. "Mode" button	Changing between the available modes
E. "-/+" buttons	These buttons make adjustments to the highlighted On Board Programming (OBP) parameter. Holding both buttons down at the same time launches the Omni into the Settings menu screen.
F. Navigation	To navigate the OBP screens
G. Communication cable	Connector to link the Omni to the R-net system
H. IR transmitter	Allows the Omni to replicate commonly used IR devices, such as remote controls for TV's, DVD's, Cable/Satellite or environmental controls such as automatic door openers
I. Charge connector	Input for the battery charger
J. 1+2 User switch jacks	There are two jack sockets (3,5 mm/1/8") which provide connection to the user switches.
K. 1+2 D-Types SID connectors	Specialty Input Device (SID) connectors, any specialty type of input device for connection to the Omni, e.g. head control or mini joysticks
L. External on/off input	There is a jack socket (3,5 mm/1/8") which provide connection to a user operated on/off button (optional)
M. Sip and puff input	This input will accept a 3.5mm/1/8" pipe connected to a Sip and Puff mouthpiece

2.3 Display screen

Power-up can be done via the on/off button on the Omni interface (B in figure 2). When switching on the Omni the external controller will also become active. Some of the external input devices can also be used to power-up the system. Depending on the type of input device the wheelchair has. Adjustments and settings can be operated by the Omni and by the external controller.

When the Omni is powered-up the screen as in figure 3 will appear.

- 1. Battery indicator Displays the power level of the battery
- Speed indicator Displays a graph of the actual speed
- Kilometre (Mileage) counter Displays the total distance driven
- 4. Kilometre (Mileage) an hour counter Displays the actual speed
- Maximum speed Displays the maximum speed limit as set by the user
- Current profile Displays the selected profile
- 7. Clock Displays the time
- 8. Adjustment option (figure 4) Displays the selected adjustment options

2.4 User menu

By holding the - and + button down at the same time, the user menu will open (figure 5). In the user menu several functions can be provided. Forward and revers commands will change the selected function. At the right hand of each line, there will be either a > or a number (surrounded by < >). A > symbol indicates that a right command will enter that function.<3> indicates that left and right commands will change the selection for that function.

3 Driving the wheelchair

3.1 Switching the wheelchair on or off

To be able to drive or operate the electronic adjustment options of the wheelchair, the controller must be switched on. Press the on/off button (B figure 2) or the external on/off button.

3.2 Driving the wheelchair

Driving and operating the wheelchair only can be done by the external input device (chapter 5)

3.3 Speed

The speed bar on the display screen will display the maximum speed. Speed can be controlled with the external input device while driving. Responsiveness of the wheelchair will depend on the type of input device connected to the Omni.

3.4 Driving profile

The Omni can also be set to make the wheelchair suitable for different driving profiles or environments (C figure 2). For example, selecting the profile for indoors will ensure that the wheelchair will react more 'calmly'. Once outside, the profile can be adjusted to a more 'robust' setting. The name and number of the current profile is displayed below the speed bar of the display screen.

3.5 Selecting a profile

To switch to a different profile, just press the 'profile' button (C figure 2). Press the 'profile' button until the most appropriate profile has been selected.

4 Operating electrical adjustments

Not every wheelchair has been equipped with electronic adjustment options. We make a distinction between four different adjustment options to the seating system (figure 6):

- Tilt adjustments
- High/low adjustments
- Backrest adjustments
- Legrest adjustments

Note!

If you adjust the seat height by using the high/low option or use the electrical tilt adjustment 0 - 45°, the speed will be reduced for your safety.

4.1 Selecting the desired adjustment option

- 1. Switch on the controller
- Press the 'mode' button of the Omni interface (D figure 2) to select the 'adjustement options' mode. It is possible to follow the status on the display screen.

To select and operate the electronic adjustment options the external input device should be used (figure 7). Functionallity of the input device depends on the programming of the Omni. In general:

- Move the joystick to the left or the right to select the desired adjustment option (1 + 2). The selected adjustment option will be visible on the display screen.
- Moving the joystick forwards and/or backwards to activate and control the selected electronic adjustment option (3 + 4). Move the joystick forwards/backwards until the desired adjustment option has been attained.

	Move joystick backwards	Move joystick forwards
Tilt adjustment	The entire chair will tilt backwards	The entire chair will tilt forwards
High/low adjustments	The entire chair will be raised	The entire chair will be lowered
Backrest adjustments	The backrest will tilt backwards	The backrest will tilt forwards
Legrest adjustments	The legrest angle will increase,	The legrest angle will decrease,
	the footplate will be raised	the footplate will be lowered

5. Press the 'mode' button to return to the 'drive' mode.

5 Input devices

Sunrise Medical offers several input devices to connect to the Omni:

- 1. External joystick
- 2. External joystick lite
- 3. Mini flex switch joystick
- 4. Micro pilot joystick
- 5. Micro guide joystick
- 6. Proximity head array

5.1 External joystick

The external joystick (figure 8) is a proportional joystick. This means that the further you push the joystick from the standard position, the faster the wheelchair will move. Releasing the joystick will stop the wheelchair.

- Deflection and 100 grams pressure
- Sunrise Medical offers the external joystick mounted in the worktop (fixed, mechanical or electrical flip-down), chin arm (mechanical or electrical) and on the armrest.

5.2 External joystick lite

The external joystick lite (figure 9) is comparable with the external joystick. The reaction is also proportional, but the lite version requires 50 percent less force than the external joystick.

- Deflection and 50 grams pressure
- Sunrise Medical offers the external joystick lite mounted in the worktop (fixed, mechanical or electrical flip-down), chin arm (mechanical or electrical) and on the armrest.

5.3 Mini flex switch joystick

The mini flex switch joystick (figure 10) is a small, black joystick with 5 different switches: to the front, rear, right and left and the fifth is activated by pressing down on the joystick shaft. Pressing down the shaft will turn the controller on or off.

The shaft is available with a hard plastic tip or a flexible rubber tip. This rubber tip can be cut to the desired size.

 Small pressure required to activate the switch

 Sunrise Medical offers the mini flex switch joystick mounted on the chin arm (mechanical or electrical).

5.4 Micro pilot joystick

Micro pilot is a miniature joystick based on a different technology than other mini "throw" based joysticks (figure 11). As little as 10 grams of force with virtually no joystick deflection, will activate the chair. This force is adjustable from 10 - 50 grams. All internal metal construction means outstanding durability.

- Approximate joystick dimensions: diameter 5/8", height without knob 2 ¹/₂"
- Non deflection and 10 to 50 grams pressure
- Sunrise Medical offers the micro pilot mounted in the worktop (fixed) and on the chin arm (mechanical or electrical)

5.5 Micro guide joystick

Micro guide (figure 12) has the same size and durability as the micro pilot but with joystick deflection. For clients that prefer the feedback of movement. Deflection with a standard knob will be 1/4" neutral to forward joystick. For activation a force of 40 grams is needed.

- Approximate joystick dimensions: diameter 5/8", height without knob 2.8"
- Deflection (1,5 cm movement) and 40 grams
 pressure
- Sunrise Medical offers the micro guide mounted in the worktop (fixed) and on the chin arm (mechanical or electrical)

5.6 Proximity head array

The proximity head array (figure 13) is a head control with 3 sensors. All functions are operated by movements of the head. The sensors are mounted inside the right and left wings. Control movements will also be in those directions. The sensor mounted inside the back pad controls movement in a forward or a forward and reverse direction. Depending on how the chair is programmed.

5.7 Spec switch

In case there are spec switches connected to the Omni interface, they have the following functions (figure 14):

- Black button: On/off switch
- Red button: 'mode' button

• Green button: D-types SID connectors More information about these functions can be found in 2.2 Main components.

In case the external (lite) joystick is used as chin control, there are two black switches, one left of the joystick and one at the right side (figure 15). Functionality of these switches is depending on programming.

6 Lights

Depending on programming, with the Omni it is possible to control the following lights:

- Lights
- Hazard Lights
- Direction indicators

7 Troubleshooting

If the wheelchair will not function while the batteries are fully charged, check the following points before consulting your dealer:

- Switch the controller off and then switch it on again. Check to see if the malfunction has been solved.
- Check if the free wheel switch was switched to Drive.
- Check if the joystick was in the 0 position when the controller was switched on. In other words, the joystick must not be moved when the controller is being switched on or off.

Malfunctions list

An extended malfunction list can be found in the appendix of the service manual (for qualified specialists only). The service manual can be found on www.SunriseMedical.eu

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