ADD ON BIKE  
AANKOPPELFIETS  
KIT HANDBIKE


Attitude  
Attitude Junior  
Attitude Hybrid

DIRECTIONS FOR USE  
GEBRUIKERSHANDLEIDING  
NOTICE D'UTILISATION

  
QUICKIE®

000690707.02

**IF YOU ARE VISUALLY IMPAIRED, THIS DOCUMENT CAN BE VIEWED IN PDF FORMAT AT [WWW.SUNRISEMEDICAL.CO.UK](http://WWW.SUNRISEMEDICAL.CO.UK)**

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#### **Add-on Bike Components**

EN

We at SUNRISE MEDICAL have been awarded the ISO-9001 certificate, which affirms the quality of our products at every stage, from R&D to production. This products meet the requirements in accordance with EC guidelines. Options or accessories shown are available at extra cost.

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04

**Attitude:**

1. Front wheel
2. Fork
3. Handle
4. Hubmotor
5. Crank with handle
6. Gear change
7. Battery
8. Brake
9. Parking brake
10. Frame
11. Docking system
12. Release Lever
13. Docking Pins (wheelchair).

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

**Contents**

<b>Definitions</b>	<b>4</b>
<b>Foreword</b>	<b>5</b>
<b>Use</b>	<b>5</b>
<b>Area of application</b>	<b>5</b>
<b>1.0 General safety notes and driving limits</b>	<b>6</b>
Safety instructions – prior to each ride.	7
Safety instructions – while driving	7
Safety instructions - after your ride:	7
<b>2.0 Warranty</b>	<b>8</b>
<b>3.0 Handling</b>	<b>9</b>
Delivery:	9
Unpacking:	9
Manual version:	9
Hybrid version:	9
Start Up:	9
Storage:	9
Transport:	9
Security:	9
<b>4.0 Adjusting the Pedal Position</b>	<b>10</b>
Adjustment:	10
Handles:	10
<b>5.0 Adjusting the Connecting Frame</b>	<b>12</b>
Adjustment:	12
<b>6.0 Docking and Dedocking</b>	<b>13</b>
Docking: Fig.6.1 - Fig.6.5	13
De-Docking: Fig.6.6 - Fig.6.10	14
<b>7.0 Functional Components</b>	<b>15</b>
Connecting frame (Fig.7.1):	15
Brakes:	15
Parking brake (Fig.7.2):	15
Main brakes (Fig.7.3):	15
Performance brake lever (Fig.7.4):	15
Coaster brake (Fig.7.5):	17
Coaster brake integrated in the gear hub (Standard on the Attitude Junior)	17
Derailleur gears/gear hub (Fig.7.6-Fig.7.7):	17
Trigger switch (Fig.7.6-Fig.7.7):	17
Grip shift switch (Fig.7.8):	17
Angle- and depth adjustment of the crank position	18
Adjusting the angle of the frame (Fig.7.9):	18
Adjusting the depth of the frame (Fig.7.10):	18
Adjusting the chain-tension (Fig.7.11)	18
Length-adjustment of the cranks (Attitude Junior only)	19
<b>8.0 Maintenance</b>	<b>20</b>
Checks to be carried out before use:	20
<b>9.0 Service Period</b>	<b>21</b>
Servicing:	21
<b>10.0 Tyre Maintenance</b>	<b>21</b>
Fitting and repairing tyres:	21
Inflation:	21
<b>11.0 Maintenance Brakes</b>	<b>21</b>
Maintaining the Brakes:	21
<b>12.0 Cleaning/Hygiene</b>	<b>21</b>
Hygiene when being reused:	21
<b>13.0 Disposal / Recycling of Materials</b>	<b>22</b>
Materials Used:	22
<b>14.0 Nameplate</b>	<b>22</b>
<b>15.0 Torque</b>	<b>23</b>
Torque Settings:	23
<b>16.0 Technical Data</b>	<b>23</b>

**Definitions**

**Definitions of words used in this manual**

Word	Definition
 <b>DANGER!</b>	Advice to the user of potential risk of serious injury or death if the advice is not followed
 <b>WARNING!</b>	Advice to the user of a risk of injury if the advice is not followed
 <b>CAUTION!</b>	Advice to user that potential damage to equipment may occur if the advice is not followed
<b>NOTE:</b>	General advice or best practice
	Reference To Additional Documentation

**NOTE:**

- Please keep a note of your local service agent's address and telephone number in the space provided.
- In the event of a breakdown, contact them and try to give all relevant details so they can help you quickly.
- The Add-on Bikes shown and described in this user guide may not correspond in every detail exactly to your own model. However, all instructions are still entirely relevant, irrespective of detail differences.
- The manufacturer reserves the right to alter without notice any weights, measurements or other technical data shown in this manual. All figures, measurements and capacities shown in this manual are approximate and do not constitute specifications.

Dealer signature and stamp

## Foreword

### Dear Customer,

We are very happy that you have decided in favour of a high-quality product from SUNRISE MEDICAL.

This owner's manual will provide numerous tips and ideas so that your new Add-on Bike can become a trustworthy and reliable partner in your life.

For Sunrise Medical, it is very important that we have a good relationship with our customers. We like to keep you up-to-date about new and current developments at our company. Keeping close to our customers means: fast service, as little red tape as possible, working closely with customers. When you need replacement parts or accessories, or if you just have a question about your Add-on Bike – we are there for you.

We want you to be satisfied with our products and service. At Sunrise Medical we are constantly working to develop our products further. For this reason, changes can occur in our palette of products with regard to form, technology and equipment. Consequently, no claims can be construed from the data or pictures contained in this user's manual.

**The management system of SUNRISE MEDICAL is certified to EN ISO 9001, ISO 13485 and ISO 14001.**

As the manufacturer, **SUNRISEMEDICAL**, declares that the add-on bikes conform to the 93/42/EEC directive.

Please contact your local, authorised SUNRISE MEDICAL dealer if you have any questions concerning the use, maintenance, or safety of your Add-on Bike.

In case there is no authorised dealer in your area, or you have any questions, contact Sunrise Medical either in writing or by telephone.

**Sunrise Medical**  
**Thorns Road**  
**Brierley Hill**  
**West Midlands**  
**DY5 2LD**  
**England**  
**Phone: 0845 605 66 88**  
**Fax: 0845 605 66 89**  
**www.SunriseMedical.co.uk**

Do not use your Add-on Bike until this Owner's Manual has been read and understood. Please follow the instructions contained within the Owner's Manual.

## Use

The ATTITUDE is a mechanical propelled / powered drive unit, which is attached to a wheelchair, to support the mobility and integration of the wheelchair user in daily life. It is designed for personal use indoors or outdoors.

**The maximum weight limit (includes both the user and any weight of accessories fitted to the Add-on Bike) is marked on the serial number label. The serial number label is located on the frame as well as in the owner's manual.**

Warranty can only be taken on if the product is used under the specified conditions and for the intended purposes.

### Lifetime

The expected life of the add-on bike is 5 years, provided that:

- It is used in strict accordance with the intended use.
- All service and maintenance requirements are fulfilled.

### DANGER!

- Do not fit any unauthorized electronic equipment.

## Area of application

The add-on bike provides the user with the potential to safely cover long distances in an ergonomically and ecologically efficient manner, using their own wheelchair. It expands the activity radius significantly.

In general we recommend the use of a rear wheel extension or additional weights to improve the traction of the driving wheel especially when driving uphill. The attached wheelchair must be equipped in accordance with road traffic regulations when used on public roads, pavements and spaces.

The max user weight is 100 kg. The max load for transported objects (e.g. in bike panniers) is 5kg each.

Please also check the Owner's Manual of the bike panniers.

The add-on bike can be operated by people who:

- are physically and mentally capable of safely operating the device in all operating situations (steering, braking, pedalling..) and can meet the legal requirements for use on public roads.
- Have been advised and trained in its use by Sunrise Medical or an authorized Sunrise Medical dealer

### Operating Conditions

The ATTITUDE should be able to cope with tarmac surfaces and firm paths with different topography. Avoid driving on unpaved or loose surfaces (e.g. on loose gravel, in sand, mud, snow, ice or through deep puddles of water), as this may expose the user to unforeseen risks.

## 1.0 General safety notes and driving limits

The technology and construction of this add-on bike have been designed to provide maximum safety. International safety standards currently in force have either been fulfilled or exceeded in parts.

Nevertheless, users may put themselves at risk by improperly using their add-on bikes. For your own safety, the following rules must be strictly observed.

Unprofessional or erroneous changes or adjustments increase the risk of accident. As an add-on bike user, you are also part of the daily traffic on streets and pavements, just like anyone else. We would like to remind you that you are therefore also subject to any and all traffic laws. Be careful during your first ride in this add-on bike. Get to know your add-on bike.

### **WARNING!**

- NEVER exceed the maximum load of 100 kg for driver excluding any items carried on the add-on bike.
- If the max user weight of the attached wheelchair is less than 100 kg, this lower limit will be valid for the complete system.
- Exceeding the maximum load can lead to damage to the add-on bike, loss of control, or serious injury to the user and other people.
- Always supervise children using the Attitude junior. Familiarize them with all functions and options on the bike, in particular with the function of the driving- and parking-brake on the bike. It's the responsibility of the supervising person to make sure the docking is safely locked before the ride and maintenance is being conducted as specified in the maintenance plan.
- Exercise caution when using the add-on bike. For example, avoid negotiating obstacles such as, steps, curbs, path edges or dropping down gaps without slowing down or braking first.
- Your add-on bike has not been designed for transportation in a vehicle. DO NOT sit on the add-on bike during any transportation in a vehicle; always use a properly secured vehicle seat!!
- This add-on bike should be used exclusively to convey one person at a time. Any other use does not comply with the intended purpose.
- Use of the add-on bike is only recommended for those who have the right physical and psychological constitution to be able to cope with public traffic.
- Please note that you are subject to any and all traffic laws when riding in public traffic.
- As a beginner, you should ride particularly carefully at first.
- Familiarise yourself with the braking distances from various speeds.
- Familiarise yourself with the general traffic laws!
- In order to achieve a stable ride when travelling straight ahead, wherever possible, you should avoid hasty steering movements when travelling at speed or going round corners.
- A three-wheeled device is always less stable than a four-wheeled device.
- When you are riding, braking or manoeuvring, always hold the handles firmly with both hands. Failure to do so could increase the risk of an accident through steering error.
- The add-on bike is a mechanical propelled / powered drive unit. The speed must be adjusted to suit the technical capabilities of the bike, terrain and operator. The max speed is 25 km/h. Never exceed this max speed.
- At night and dawn/dusk, always use the lights.
- When you are travelling long distances, we recommend that you take a spare battery with you.
- Only ride on asphalt paths, or firm, smooth surfaces.
- Use caution when riding over curbs, steps, potholes, etc. as there is an increased risk of tipping over.
- Adjust your method of driving to your capabilities.

- Ride slowly and carefully, particularly when riding round corners.
- Do not put your hands in the chain area.
- Do not touch the motor, because it can get very hot.
- Your bike should be serviced by your dealer at regular intervals. (min. every year)
- As with all moving parts, there is a risk that fingers could become trapped. Please always handle carefully.
- If the add-on bike is subject to direct sunlight / low temperatures for long periods, parts of your add-on bike could become very hot (>41°C) or very cold (<0°).
- Always make sure that the quick tensioner on the drive wheel is set properly.
- The add-on bike should not be used in storm conditions, heavy rain/snow or on slippery or damaged surfaces.
- You should only use product combinations which have been approved by Sunrise Medical.

### **DANGER!**

**CHOKING HAZARD** – This mobility aid uses small parts which under certain circumstances may present a choking hazard to small children.

The wheelchairs shown and described in this user guide may not correspond in every detail exactly to your own model. However, all instructions are completely relevant, regardless of possible detail differences.

**NOTE:** The add-on bikes shown and described in this user guide may not correspond in every detail exactly to your own model. However, all instructions are completely relevant, regardless of possible detail differences. The manufacturer reserves the right to alter without notice any weights, measurements or other technical data shown in this manual. All figures, measurements and capacities shown in this manual are approximate and do not constitute specifications.

## Safety instructions – prior to each ride.

### CAUTION!

- Check the conditions of your wheels, (including the wheels of your wheelchair) visually, (spokes and rims, damages of tyres) as well as the pressure of the tyres. The correct pressure is indicated on the tyre and for the add-on bike is normally 3-4 bar.
- Low tyre pressure increases the risk of tipping over especially in bends, it also effects the braking action, the wear and tear and the driving performance due to increased rolling resistance.
- Always check that the tread depth of your tyres is Okay.
- Check the function of the brakes; the add-on bike may not be used if the brakes are not 100% functioning!
- Avoid damage to the Bowden cables for the gear-shift and braking units. There is a risk of kinking the cable when you transfer to/from the bike and if you operate the crank incorrectly.
- Check the fork by visual inspection for damage (e.g. cracks and breakages)
- Check that all components, especially the bolts, are tight.
- Check the connection between the add-on bike and your wheelchair.
- Check the frame and the components of your wheelchair for any damage, (e.g. cracks)
- Check that the lights are functioning and the reflectors are sound; they must be visible all the time during use.
- Ensure that your feet cannot slip off the foot plate of your wheelchair, if necessary use a special fixation device.
- Check all cables and electrical connections.
- Remember, the max load for the saddle bags is 5 kg each.
- Check that the battery is locked correctly and fully charged.
- Switch the battery on and check the light.
- Always wear a helmet.

## Safety instructions – while driving

### WARNING!

- It is recommended Make yourself familiar with the product in flat terrain before driving up or down slopes.
- Always hold the handles with both hands firmly as otherwise, there is a high risk of an accident through steering error.
- Always adapt your speed to your driving capability and to the traffic - / terrain conditions
- Use particular caution when approaching stairs, edges, drops or other hazard areas.
- Reduce your speed to walking pace and lean your body inwards, when driving round a bend.
- When waiting e.g. at a pedestrian crossing, a traffic light, on hills or slopes or at ramps of any type, always activate the brake.
- Due to the increased turning circle, turning can be sometimes impossible e.g. in corridors and other confined spaces.
- Only use the defined brackets to fix objects, (additional weights, bike panniers)
- When driving on pedestrian areas, observe the maximum permitted speed (walking speed)
- Always follow the road traffic regulations when driving on public roads and footpaths,
- Avoid driving on unpaved or loose surfaces (e.g. on loose gravel, in sand, mud, snow, ice or through deep puddles of water).
- Avoid abrupt steering actions
- Avoid driving in parallel to slopes
- Avoid turning on slopes
- On slopes the traction of the driving wheel is lowered and the braking effect is considerably reduced. The driving style and speed should be always adjusted so that the product can be safely stopped at all
- Times without any risk.
- The traction can be improved by using an additional weight or the rear wheel extension.
- Towing or using a trailer is not permitted.

- Driving on stairs is not permitted.
- Obstacles e.g. curbs must be run over always frontal / rectangular with very low speed to avoid tipping over and potential injuries. The max obstacle height is 50.0 mm, but this depends on the distance of the stand wheels and castors to the ground and therefore of the adjustments of the add-on bike – wheelchair system.
- It must be avoided that any wheelchair components or parts of your body get trapped when driving over or passing obstacles. This could result in falls and serious injuries as well as damage to the
- Add-on bike and the wheelchair.
- When driving on wet surfaces, there is an increase slip hazard because the traction of the tyres is reduced. Adapt your driving style accordingly.
- Always use the main brake to slow down your bike. In case the main brake will fail the parking brake can be used as emergency brake once.
- Sliding of the driving wheel during braking can be avoided by moving your body weight forward.
- Avoid strong braking actions while driving through bends.

## Safety instructions - after your ride:

- Switch off the product immediately when it is not in use to avoid accidental operation / movements caused by accidentally touching the crank handle.
- Switching off helps to preserve the battery charge.

## 2.0 Warranty

### **THIS GUARANTEE DOES NOT AFFECT YOUR LEGAL RIGHTS IN ANY WAY.**

Sunrise Medical\* provides a guarantee, as set out in the warranty conditions, for products to its customers covering the following.

#### **Warranty conditions:**

1. Should a part or parts of the product require repair or replacement as a result of a manufacturing and/or material fault within 24 months, then the affected part or parts will be repaired or replaced free of charge. The warranty will only cover manufacturing defects
2. To enforce the warranty, please contact Sunrise Medical Customer Service with the exact details of the nature of the difficulty. Should you be using the product outside the area covered by the Sunrise Medical customer service agent, repairs or replacement will be carried out by another agency as designated by the manufacturer. The product must be repaired by a Sunrise Medical designated Customer Service agent, (dealer).
3. For parts, which have been repaired or exchanged within the scope of this warranty, we provide a warranty in accordance with these warranty conditions for the remaining warranty period for the product in accordance with point 1).
4. For original spare parts which have been fitted at the customer's expense, these will have a 12 months guarantee, (following the fitting), in accordance with these warranty conditions.
5. Claims from this warranty shall not arise, if a repair or replacement of a product or a part is required for the following reasons:
  - a. Normal wear and tear, which include but is not limited to the following parts where fitted; batteries, armrest pads, upholstery, tyres, brakes shoes, ferrules, etc.
  - b. Any overloading of the product, please check the EC label for maximum user weight.
  - c. The product or part has not been maintained or serviced in accordance with the manufacturer's recommendations as shown in the user instructions and/or the service instructions.
  - d. Accessories have been used which are not specified as original accessories.
  - e. The product or part having been damaged by neglect, accident or improper use.
  - f. Changes/modifications have been made to the product or parts, which deviate from the manufacturer's specifications.
  - g. Repairs have been carried out, before our Customer Service has been informed of the circumstances.
6. This guarantee is subject to the law of the country in which the product was purchased from Sunrise Medical“

\* Means the Sunrise Medical facility from which the product was purchased.

#### **Additional Notes For Australia Only:**

- i. For goods provided by Sunrise medical Pty Ltd in Australia, our goods come with a guarantee by Sunrise Medical that cannot be excluded under Australian Consumer Law.
- ii. You are entitled to a replacement or refund for a major failure and for compensation for any foreseeable loss or damage.
- iii. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.
- iv. The benefits to you given by this warranty are in addition to your other rights and remedies under a law in relation to the goods to which the warranty relates.

## 3.0 Handling

### Delivery:

Your new add-on bike is delivered fully assembled in a cardboard box. In order to avoid any damage during transportation, loose parts or parts to be fitted, are packed separately in the cardboard box.

### Unpacking:

- Check the shipping cardboard box for external damage which may have been caused during transportation.
- Remove all packaging material.
- Take the add-on bike carefully out of the box and check that it is complete.
- Check all parts for surface damage, scratches, cracks, dents, distortions and other defects.
- The scope of delivery consist usually of the ATTITUDE add-on bike:

### Manual version:

Connecting frame with complete docking system.

### Hybrid version:

Connecting frame with complete docking system, display, Lithium Ion battery and charger.

**NOTE:** Please make a note of the Battery Locking Key identification number. This number can be used to order on-line spares at:  
<https://keyservice.axa-stenman.com/>

### Start Up:

If the Attitude Hybrid has not been use for 48 hours or more, please turn the battery on.

**NOTE:** Should you find any damage, or if any parts are missing please let us know immediately.

### CAUTION!

- The add-on bike may only be used with the delivered docking system and the suitable docking pins for your wheelchair. The add-on bike must be handed over ready to use from Sunrise Medical or an authorized Sunrise Medical dealer.
- Before using for the first time, the battery should be fully charged.

### Storage:

The add-on bike should always be stored in an environment with low humidity, on an easy to clean surface at ambient room temperature (+15°C to + 25°C).

### CAUTION!

To avoid creating a hazard the Lithium Ion battery, must be stored according the description given in the section Part 2 Neodrives Battery.

### Transport:

Your add-on bike has not been designed for transportation in a vehicle.

### WARNING!

- DO NOT sit on the add-on bike during transportation in a vehicle.
- Always use a properly secured vehicle seat.
- The product and all its components must be secured during transport so that they are not damaged (e.g. by falling over) and do not become any risk to other people.

### CAUTION!

- When loading make sure that the cables are not caught up, become kinked or otherwise damaged.
- The add-on bike may not be used with damaged cables.

### Security:

- To protect your Attitude against theft and also against unauthorized use, please use a suitable lock, such as those designed for bicycles.
- For extra security on the Attitude Hybrid, please remove the SMMI and the battery.

Key Number:

## 4.0 Adjusting the Pedal Position

### Adjustment:

The default position of the crank is based on data which depicts the average best fit. If this is not suitable for you, please contact your authorized Sunrise Medical dealer.

The pedal position can be adjusted to the user, (arm length and body position), by an authorized dealer in accordance with available adjustment options,

### ⚠️ WARNING!

- The cranks must not touch the knee or thigh while pedalling, especially through bends.
- Do not allow the elbows to completely extend when the crank handles are furthest away from the body.
- If you have poor seating stability, or weak trunk muscles, a suitable restraint must be used. Please contact your therapist, doctor or authorized dealer.

The pedal position can be adjusted at least 0- 50.0 mm forward and 0 - 10 mm upward, (Fig.4.2).

### Handles:

### ⚠️ WARNING!

- The handles must be held firmly with both hands whilst driving and held so that the cables always are oriented upwards. Holding the cranks in a horizontal position or upside-down will result in cable-damage.
- Cables should never be kinked or trapped, as this could cause damage and loss of function.

Fig.4.1


Fig.4.2


Fig.4.3


Fig.4.4




Fig.4.5

**Fig.4.6**

Docking top	Docking bottom	(X) mm	(Y) mm	(S) mm	(a) mm	(L) mm
position 1	position 1	640	370	25	75 - 95	35
position 1	position 1	650	370	25	75 - 95	38
position 1	position 1	660	370	25	75 - 95	41
position 1	position 1	670	370	35	75 - 95	44
position 1	position 1	680	370	35	75 - 95	47
position 1	position 1	690	370	35	75 - 95	50
position 2	position 2	660	365	15	75 - 95	36
position 2	position 2	670	365	25	75 - 95	39
position 2	position 2	680	365	25	75 - 95	42
position 2	position 2	690	365	25	75 - 95	45
position 2	position 2	700	365	25	75 - 95	48
position 2	position 2	710	365	35	75 - 95	51

Frame angle:	95°
Seat depth:	440 mm
Seat height rear:	450 mm
Seat height front:	500 mm
COG:	80 mm

## 5.0 Adjusting the Connecting Frame

### Adjustment:

The connecting frame can be extended by approximately 20mm to fit a wheelchair with a wider front frame.

Therefore the 12 allen screws on the top and bottom cover of the adjuster clamp need to be loosened. The connecting frame can then be pulled out on both sides, (Fig.5.1). Tighten all allen screws with a torque of 10Nm.

### ⚠ CAUTION!

- Ensure that the connecting frame is extended equally on both sides
- Never extend the connecting frame too far.
- Ensure that the ends of the frame tube are within the "MAX" marking printed on top of the adjuster clamp, (see Fig.5.2).

Fig.5.1



Fig.5.2


## 6.0 Docking and Dedocking

### ⚠️ WARNING!

- Always make sure that the device is switched off before docking and dedocking.
- Ensure the wheelchair is fully connected with the add-on bike before use, by using the handgrip A.
- As with all moving parts, there is a risk that fingers and / or clothing could become trapped. Please always handle carefully.

Refere to the docking check list below:




<p><b>Docking Pins Fully locked In</b></p>		
<p><b>Docking Levers Fully Locked In</b></p>		

### ⚠️ CAUTION!

- Always keep the docking system clean.
- Clean it after every trip paying particular attention to sand and mud.

### Docking: Fig.6.1 - Fig.6.5

- Drive your wheelchair as close as possible to the add-on bike.
- The 2 docking pins on your wheelchair should stay directly in front of the lower docking connector.
- The add-on bike should stand on an even surface.
- Lock the brakes of your wheelchair.
- Pull it towards you by using the handgrip, until the 2 docking pins lock into A.
- Now push the add-on bike upwards/forwards until the stand wheels and the castors leave the ground B.
- Swing the upper docking connector downwards until it stops, C.
- Lower the add-on bike until the upper docking connector is in contact with the front frame of your wheelchair.
- Pull the red locking lever backwards until it locks firmly and completely, D.
- Ensure the wheelchair is fully connected with the add-on bike.





**De-Docking: Fig.6.6 - Fig.6.10**

**⚠ WARNING!**

The de-docking procedure should take place on an even surface.

**NOTE:** De-docking follows the opposite sequence as described for Docking.

- Lock the wheelchair brakes.
- Open the red locking levers, A.
- Push the add-on bike forward/upwards until the upper docking connector has no contact with the front frame of your wheelchair, B.
- Swing the upper docking connector upwards, C.
- Lower the add-on bike until the stand wheels and castors have contact to the ground.
- Pull the red release lever by the handgrip, D.
- Pull the add-on bike downwards towards yourself until the docking pins are totally disengaged.
- Hold this position and push the add-on bike forward so that it is clear of your wheelchair, E.

## 7.0 Functional Components

### Connecting frame (Fig.7.1):

The connecting frame acts as a parking stand, and facilitates easy docking/manoeuvring when the add-on bike is disconnected from the wheelchair. The height of the stand wheels can be adjusted using the bolts, (A).

#### **⚠ WARNING!**

- The stand wheels must always stay min 30.0 mm above the ground, when the add-on bike is connected to the wheelchair.
- The stand wheels must be always adjusted so that the stand of the add-on bike remains stable when the add-on bike is disconnected from the wheelchair.

### Brakes:

The add-on bike has two brakes that work independently from one another.

### Parking brake (Fig.7.2):

The parking brake is assembled to the fork and locks in automatically when activated. Therefore the red lever must lock in. The add-on bike should not move with activated parking brake to make sure parking is safe. To release the parking brake just push the red lever.

### Main brakes (Fig.7.3):



The brakes allow you to safely and comfortably brake your add-on bike.

### Performance brake lever (Fig.7.4):

The Performance brake lever (optional) gives you the possibility to adjust the lever-feel and the way that the main brake is applied.

Turning the red adjusting screw (A) on the brake lever clockwise results in a harder lever feel and a better, linear, controllability of the brake-force.

Turning the red adjusting screw (A) on the brake lever counter-clockwise will result in a softer lever-feel and more progressive brake-force.

 **WARNING!**

- If the main brakes should fail, then the parking brakes may be used once for emergency braking.
- The main brakes must be repaired immediately before next use.
- Hard and sudden application of the brakes may cause you to fall forward with your upper body, resulting in possible injury. Therefore always stabilise yourself by holding onto the handles firmly with both hands.
- Braking power can be severely decrease by any or all of the following:
  1. Worn tyre tread.
  2. Air pressure which is too low.
  3. Soiled and wet tyres.
  4. Worn brake pads.
- Ensure that the brake pads and brake discs do not become contaminated with oil and/or grease. If this happens the brake pads must be replaced. The brake disc must be professionally degreased with a brake cleaner.
- The disc brake is not self-adjusting and needs to be re-adjusted when the brake pads wear.
- Be aware that the condition of the ground surface can significantly affect braking efficiency.

 **WARNING!**

- The braking effect will be greatly reduced if the surface is wet, dirty, sandy, stoney, uneven or otherwise compromised by environmental conditions. It is best to avoid such conditions, but if this is not possible adapt your driving style to avoid the risk of potentially serious accidents.
- Dirty and wet rims and brake pads also reduce the braking effect when using rim brakes. For this reason, you should always adapt your method of riding to the environmental conditions in order to avoid serious accidents.
- Always check the efficiency of the brakes before starting your journey to ensure that they are working correctly. The brakes should be checked regularly by your dealer and adjusted if required.

### Coaster brake (Fig.7.5):

The coaster brake option allows you to brake by pedalling backwards. This option is required particularly for users with restricted or no finger function.

The braking system can only be uncoupled for manoeuvring (Fig.7.5) by moving the lever A.

When riding or moving forward, the system will be put automatically back into an operational position.

#### **⚠ WARNING!**

- Riding with the back-pedalling brake deactivated is not permitted and can lead to an increased risk of accidents.

### Coaster brake integrated in the gear hub (Standard on the Attitude Junior)

The gear hub Shimano Nexus has an integrated coaster brake function, which can be activated by pedalling backwards.

### Derailleur gears/gear hub (Fig.7.6-Fig.7.7):

Gear shifting is carried out using a trigger or grip shift switch, which are attached to the handgrips.

The Derailleur gears/gear hub enables you to easily keep moving, regardless of other road conditions.

The derailleur gears may only be switched while pedalling. Reducing the pedaling force will allow a fast change of gears. The gear hub can also be operated when standing still.

### Trigger switch (Fig.7.6-Fig.7.7):

The trigger switch on the left handle, activates the 3 gear hub and the trigger switch on the right handle, activates the 10 derailleur gears.

When operating the big lever the next lower gear will be activated.

When operating the smaller lever the next higher gear will be activated.


The chosen gear is indicated at the operating unit.

### Grip shift switch (Fig.7.8):

The grip shift switch at the right handle operates the 10 derailleur gears. Turning the gripshift switch to the left activates the next lower gear. Turning the grip shift switch to the right activates the next higher gear.

#### **⚠ CAUTION!**

Please contact your authorized dealer in case the derailleur gears/ gear hub does not functioning correctly

### Angle- and depth adjustment of the crank position

If your Attitude is equipped with the optional angle- and depth adjustment, you can adjust the angle and the length of the frame using a hinge and a telescopic tube on the upper part of the frame.

#### Adjusting the angle of the frame (Fig.7.9):

- Loosen the 4 allen screws (B, 2 on each side)
- Whilst holding the upper part of the frame (A)
- Adjust the upper half of the frame to the desired angle and hold
- Tighten all 4 allen screws (B) evenly with a torque of 12Nm

#### Adjusting the depth of the frame (Fig.7.10):

- Loosen the allen screw (B)
- Extend or insert the telescopic tube (A) to the Desired length
- Tighten the allen screw with a torque of 20Nm




#### Adjusting the chain-tension (Fig.7.11)

Adjusting the angle and/or the depth has an impact on the tension of the chain. Please check chain-tension and extend or shorten the chain if necessary (additional chain-links included in delivery).

On the Attitude manual equipped with 8spd hub-gearing and on the Attitude Junior, it's possible to adjust the chain tension within a certain range using the chain adjuster mounted to the right hand side dropout.

- To tension the chain, loosen the allen screw (A) and slide the pulley (B) upwards or downwards until the correct tension of the chain is achieved. Hold the pulley and tighten the allen screw (A) with a torque of 7Nm.

The chain adjuster can compensate a certain range of adjustment, but not the full range of adjustability. Therefore, an additional piece of chain with a quick-lock is included in delivery. You can extend your chain without tools by simply opening the quicklock (A) which is already in your chain.

- Locate the quick-lock (A) on your chain, clean the area around the quick-lock, hold the chain-links left and right of the quick-lock firmly and push them against each other (Fig.7.12).
- Insert the chain extension, re-assemble the quick-locks (A) on both sides and close them by pulling them apart (Fig.7.13)

The tension of the chain is correct when it can be tensioned manually in the middle of the chain by approx. 10mm. (Fig.7.14)

On the Attitude manual or Hybrid with 10/30 derailleur gears, the derailleur functions as chain tensioner and can compensate a certain range of adjustment. However, chain tension must never be too tight.

To check chain tension, shift the derailleur in the highest position (chain on the largest sprocket) and check tension: The tension of the chain is correct when it can be tensioned manually in the middle by approx. 10mm (Fig.7.14).

#### ⚠ CAUTION:

- Always tighten the allen screws to the correct torque. Over- or undertorquing the allen screws might result in failure of the handbike frame.
- Adjusting the angle and/or the depth has an impact on the tension of the chain. Please check chain-tension and extend or shorten the chain if necessary (additional chain-links included in delivery).



#### Length-adjustment of the cranks (Attitude Junior only)

The cranks of the Attitude Junior can be adjusted in length to suit the users arm-length.

To change the position of the handles on the cranks, the handles need to be removed from the cranks using a spanner.

Mind that there is a right-handed thread on the right-hand side crank, left-handed thread on the left-hand side crank. Unscrew the handles counter-clockwise on the right-handed crank but clockwise on the left-hand crank.

Position the handles in the desired position and tighten them with a torque of 35Nm.

## 8.0 Maintenance

Before any maintenance is started, turn the battery off and remove it.

**NOTE:** All Service and Maintenance procedures carried out on the Add-on Bike require the Battery Locking Key. Without this key no maintenance is possible.

### Checks to be carried out before use:

When	What	Comments
<ul style="list-style-type: none"> <li>Every 4 weeks, (depending on use).</li> </ul>	<ul style="list-style-type: none"> <li>Grease the cables where they leave the switches and the brake levers.</li> <li>Check all screw connections are secure.</li> <li>Oil the drive chain.</li> </ul>	<p>The following can be done by user or attendant:</p> <ul style="list-style-type: none"> <li>Before oiling, please remove old oil remnants from the components.</li> <li>Spray lightly with Teflon spray.</li> <li>Ensure that surplus oil/spray does not contaminate the surrounding area, (Clothes, brakes, etc).</li> </ul>
<ul style="list-style-type: none"> <li>Every 2 - 3 months, (depending on use).</li> </ul>	<ul style="list-style-type: none"> <li>Clean regularly.</li> <li>Check the unit is clean and in good condition.</li> </ul>	<p>The following can be done by user or attendant:</p> <ul style="list-style-type: none"> <li>Check for cracks in the rubber components.</li> <li>Replace any cracked rubber components or take it to your nearest authorised service centre for refurbishment.</li> </ul>

Maintenance plan	1st inspection	2nd inspection	Annual inspections
	"250 km or after 6 month"	"1000 km or after 1 year"	"every 1000 km / or every year"
Check all screws / fastening elements are firmly seated	X	X	X
Functional / safety check of all lights (if fitted)	X	X	X
Functional / safety check of the docking system (add-on bike and wheelchair)	X	X	X
Clean and oil / grease all pivot points and bearings		X	X
Check function of gearing	X	X	X
Check brake-pads for wear and if necessary, re-adjust the distance from brake-pad to brake-disc.	X	X	X
"Check function of main and parking brake (replace cables and brake pads if needed)"	X	X	X
"Check chain, gear, gear rim for wear and tear (and replace if needed)"		X	X
"Check redirection chain wheel and replace it if needed (and replace if needed)"		X	X
"Visual inspection of the fork and connecting frame for any damage (e.g. crack, formation, deformation, etc.)"	X	X	X
"Visual inspection of the wheelchair frame and back tubes for any damage (e.g. crack, formation, deformation, etc.)"	X	X	X
Check tyres (tread, cracks, damages) and replace if needed		X	X
Check, clean and oil the gear components (including chain)	X	X	X
Check the spoke tension and the rim, hub for damages on the add-on bike and the wheelchair	X	X	X
Check all cable and electrical connectors	X	X	X
Check tightness of steering shaft	X	X	X

## 9.0 Service Period

### Servicing:

If you have any questions or require any help, please contact your authorized Sunrise Medical dealer, who will be pleased to advise you with regards to servicing and repairs. Depending on how much you have used your add-on bike, we recommend that you take it to your dealer at least every 6 months, to have your add-on bike inspected by trained personnel.

## 10.0 Tyre Maintenance

The add-on bike is equipped with a 20" drive wheel. Make sure that they are always at the recommended tyre pressure, as this will significantly affect the driving characteristics. If the tyre pressure is too low, movement and manoeuvrability is made more difficult due to the increased rolling resistance and this also increases the amount of effort required.

The optimum tyre pressure (usually 3-4 bar) is indicated on each tyre.

### Fitting and repairing tyres:

Before fitting a new inner tube:

- Inspect the base of the rim and the tyre inner wall for foreign bodies and if necessary clean thoroughly
- Check the condition of the rim tape, particularly at the valve hole. Only well fitting, strong rim tape offers long-term protection against inner tube damage through burrs and sharp spoke heads.

### Inflation:

Before inflation check all round on both sides, to make sure that the inner tube is not trapped between the tyre bead and the rim. Push the valve back gently and pull it back out, so that the inner tube achieves an optimum position around the valve.

### WARNING!

After fitting, check the tyre pressure. It is essential that you maintain regulation tyre pressure and keep the tyres in perfect condition, in order to ensure your personal safety and keep your add-on bike properly.

## 11.0 Maintenance Brakes

### Maintaining the Brakes:

Check regularly the condition of the brake pads and the brake disc. After 1,500 km the brake cables should be replaced completely. The brake pads need to be replaced if the thickness measured at the thinnest point is 2,5 mm. The distance between brake-pads and disc needs to be re-adjusted as the brake-pads wear.

## 12.0 Cleaning/Hygiene

- The add-on bike should be cleaned regularly and after each bigger journey with commercially available, household cleaning agents.
- If needed only use a little water.
- Handles, display and battery may only be cleaned by wiping with a damp (not wet) cloth.
- The drive wheel can be cleaned with a soft sponge or a soft brush.
- The docking pins at the wheelchair and the docking system can be cleaned with water, a soft sponge or a soft brush (recommended after each ride).
- After cleaning please use a penetrating oil (e.g. Ballisto WD 40, Caramba...) to lubricate the moving parts of the docking.

### CAUTION!

- Never use steam or high pressure to clean the product.
- Never use any abrasive cleaning agents or aggressive, acidic cleaners.
- Make sure that water will not get into contact with the electrical connectors.

### Hygiene when being reused:

When the add-on bike is to be reused, it should be prepared carefully, and be wiped and treated with spray disinfectant on all surfaces which could come into contact with the user. If you need to do this quickly, you must use a liquid, alcohol based disinfectant suitable for medical products and devices. Please pay attention to the manufacturer's instructions of the disinfectant you are using.

### 13.0 Disposal / Recycling of Materials

**NOTE:** If the add-on bike has been made available to you as part of a charity or medical loans scheme, then it may not belong to you. If it is no longer required, then follow the instructions to return it as given by the organisation that made the add-on bike available to you.

#### Materials Used:

In the following section, there is a description of the materials used on the add-on bike, in view of the disposal or recycling of the add-on bike and its packaging. There may also be special local regulations in force with regard to disposal or recycling, these must be taken into account when disposing of your add-on bike. (This can include the cleaning or decontamination of the add-on bike prior to disposal).

- Aluminum:** Fork, wheels, connecting frame, crank, main brake lever, gearing components.
- Steel:** Docking pins, pivots knee lever docking system.
- Composite:** Handles, plugs, parking brake lever and tyres.
- Packaging:** Made of soft polyethylene, cardboard
- Battery:** Lithium Ion battery (hazardous product) Please see battery section for more detailed information.


**NOTE:** Disposal or recycling should be done through a licensed agent or authorised place of disposal. Alternatively your add-on bike may be returned to your dealer for disposal.

### 14.0 Nameplate

The nameplate is located on the frame as well as on a label in the owner's manual. The nameplate indicates the exact model designation and other technical specifications. Please provide the following pieces of information whenever you have to order replacement parts or to file a claim:

SAMPLE

ATD XXXXX-XXX	Product Name/SKU Number.
maxXX	Maximum safe slope, depends on wheelchair/ add-on bike setting, posture and physical capabilities of the user.
Kg Max XXX kg	Load Maximum.
CE	CE Mark
i	User's Guide.
<del>u</del>	Not Crash tested
XX.XX.XXXX	Date of manufacture.
SN ATD :XXXX	Serial No

As we are committed to continually improving the construction of our wheelchairs and add-on bikes, product specifications may vary slightly from the examples illustrated. All weight/ dimensions and performance data is approximate and is provided solely for guidance. All products from Sunrise Medical comply with the 93/42/EEC directive for medical products.

**CE** All add-on bikes must be used in accordance with the manufacturer's guidelines.

Sunrise Medical GmbH  
 Kahlbachring 2-4  
 69254 Malsch/Heidelberg  
 Deutschland  
 Tel.: +49 (0) 7253/980-0  
 Fax: +49 (0) 7253/980-222  
 kundenservice@sunrisemedical.de  
 www.SunriseMedical.de

## 15.0 Torque

### Torque Settings:

Torque in general is 7 Nm, unless otherwise specified.

The screws used are fitted with threadlock at the factory and can be adjusted up to 5 times after fitting. After that they should be replaced or alternatively re-secured with threadlock to medium tightness..

### ⚠ WARNING!

Secured nuts may only be used once.

Fig.15.1



Fig.15.2

Fig.15.3

## 16.0 Technical Data

General	
max user weight	100 kg
Seat width range wheelchair	280 - 460 mm
Weight (cpl)	21,5kg
Wheelsize	20" (508 mm)
Min Turning radius	2700 mm*
Min slope	5° (9%)*
Max. inclination (static)	10°*
Mix. obstacle height	35 mm*
Max. obstacle height	50 mm*
Gearing	SRAM X5 10 speed • or Dual Drive (3x10 speed)
Gear switch	SRAM trigger /grip shift switch
Motor (see also motor section 9)	Hybrid only
Motor (peak)	250 W (650W)
Torque (full load)	12 Nm (40 Nm)
Range*	120 km
Voltage	36 Volt
Speed	25 km/h
Efficiency	80%
Battery (see also battery section)	Hybrid only
Battery type	Lithium ion
Battery capacity	14,5 Ah
Battery protection rating	IP 54
Rated voltage	36,2 V
End-of-charge voltage	42 V
Total energy	515 Wh
Max discharge current	30 A
Ambient operating temperature	-20 °C to 60 °C
Ambient storage temperature	+10°C to +30°C
Ambient chargingg temperature	0 °C to 40 °C
Battery weight	3,5 kg
Display (MMI)	
Display	Monochrome
Display Size (diagonal)	2,4 " (61 mm)
Display resolution	240 x 320 pixel
Protection Rating Display	IP 65
Operating temperature	-20° to + 50°

\* Depending on the wheelchair – add-on bike combination values can vary e.g. rear wheel extension, surface conditions, user weight, driving capability, add-on weight.

# Neodrives System

## Part 1 Neodrives sMMI and Motor

## Part 2 Neodrives Battery

<b>Part 1 Neodrives sMMI and Motor</b>	<b>26</b>	<b>Part 2 Neodrives Battery</b>	<b>43</b>
<b>1.0 Safety Advice/Instructions</b>	<b>26</b>	<b>4.0 Introduction</b>	<b>43</b>
Intended use of the neodrives components:	26	Signs and symbols:	43
Permissible operating conditions / operating sites:	26	Permissible operating conditions / operating sites:	43
<b>2.0 Standard Scope of Delivery</b>	<b>27</b>	Standard scope of delivery (neodrives components):	43
Neodrive Components:	27	<b>5.0 Technical Data</b>	<b>44</b>
Technical Data:	27	Key Elements Overview:	44
<b>3.0 smart Man-Machine Interface (sMMI)</b>	<b>28</b>	Safety instructions and precautions on the use of the battery:	45
Controller:	28	<b>6.0 Safety Instructions/Precautions</b>	<b>45</b>
Fitting and removing the sMMI	29	Safety instructions and precautions on how to store the battery:	45
Functions of the sMMI	29	Safety instructions and precautions on the charging process:	45
Start menu:	30	Safety information and precautions on how to transport and	
Selecting the assistance level:	30	ship the battery:	46
Activating recuperation:	31	Safety information and precautions on the charger:	46
Braking assistant:	32	Starting up:	47
Important information:	32	Information on the operating modes:	47
Pushing aid:	33	Inserting the battery:	47
Using the pushing aid:	33	Connecting the battery to the motor cable:	47
For add-on bikes:	33	Switching on the battery:	47
For add-on bikes:	33	<b>7.0 Operation</b>	<b>47</b>
Date and time:	33	Battery capacity indicator:	48
Travel modes:	34	Detaching the battery:	49
Displaying journey information:	34	Switching off the battery:	49
Information and warning symbols:	35	Disconnecting the cable:	49
Battery capacity display (standard display):	35	Removing the battery:	49
Warning:	35	Charging the battery:	50
Braking assistant active:	35	Connecting the battery charger:	50
Service reminder:	35	Charging process:	50
Warning – temperature problem:	35	LED displays during the charging process:	50
Warning – system fault:	35	Keys:	51
Warning – general fault:	35	Cleaning the battery:	51
Easy Display mode:	36	Disposal:	51
Activating Easy Display:	36	Liability:	51
Deactivating Easy Display:	36	<b>8.0 Disposal</b>	<b>51</b>
Other settings:	36	<b>9.0 Liability</b>	<b>51</b>
USB port: (Fig.3.14)	36		
Programming options by the specialist dealer:	37		
Firmware updates and relaying them to the motor and battery:	37		
Note on the sMMI plugs: (Fig.3.15)	38		
Thermal management:	38		
Motor:	39		
Removing the drive wheel:	39		
Attaching the drive wheel:	40		
Cleaning the motor and sMMI:	41		
Motor:	41		
sMMI Unit:	41		
Transportation:	41		
Safety precautions:	41		
Error indications and possible remedies:	42		
Introduction:	43		
Important instructions :	43		
Proper use of the neodrives battery:	43		

## Part 1 Neodrives sMMI and Motor

### 1.0 Safety Advice/Instructions

#### WARNING!

It is not a legal requirement at present to wear a helmet when using a add-on bike, nevertheless for your own safety, wearing a helmet is recommended!

#### Intended use of the neodrives components:

- Your add-on bike, equipped with neodrives components by the specialist dealer on delivery, is designed as a hybrid bicycle for ordinary personal transport in public road traffic use.
- Adjustments and repairs to the add-on bike and its individual components only apply as intended use insofar as they are explained and permitted in this operating manual, in the operating manual of the add-on bike manufacturer, the instructions of the component manufacturers or other documents included when purchasing your add-on bike.
- The manufacturer accepts no liability for damage caused by negligence as a result of misuse, improper maintenance or repairs or improper use.
- It is the responsibility of the rider to check the add-on bike as required, to have any work carried out on it and to use it responsibly.
- This operating manual only describes the use of the neodrives components fitted to your add-on bike and corresponds to the state of the art at the time of print.
- The manufacturer reserves the right to make changes resulting from further development of the mechanics, software or legal requirements.

The manufacturer regards the following cases as examples of misuse of the neodrives components fitted to your add-on bike:

- Use of the drive system that contravenes the instructions and recommendations in this operating manual.
- Exceeding the technical limits laid down in this operating manual.
- Technical modifications to the neodrives components.
- Modifications to the software of the neodrives components.
- Unauthorised attachments or use of the neodrives components on bicycles or a different add-on bike to the one supplied to you.

The manufacturer shall not be liable for any damage caused by misuse of the components.

#### WARNING!

Before using the device, carefully read all safety and hazard information contained in the individual chapters of this operating manual and all other enclosed documents.

#### Permissible operating conditions / operating sites:

Neodrives components can be used at temperatures between -20 °C and +50 °C.

#### WARNING!

- Observe the information on permissible operating conditions in the ad-on manufacturer's operating instructions.
- Any limits regarding the permissible operating conditions (e.g. maximum climbing capability, maximum permissible height of obstacles, max. load, etc.) must be observed when using the add-on bike!
- Observe the safety and hazard information provided in the individual chapters of the operating manual.

## 2.0 Standard Scope of Delivery

### Neodrive Components:

- Drive motor.
- Smart MMI (display) incl. dock.
- This operating manual.

### Technical Data:

#### Drive

Range*:	120 km
Speed	25 km/h
Power rating (peak)	250 Watt (650 Watt)
Operating voltage	36 Volt
Rated torque	12 Nm
Peak torque	40 Nm
Level of efficiency	80 % (incl. electronics)
Control system performance electronics	Integrated into the wheel hub
Cassette receiver commercially available	Plug-in cassette, up to 10 times
Brake disc	From 160 mm diameter
Torque receiver	Variable torque supports can be adjusted to the dropout
Weight	4.36 kg (just the drive including plug and cable, without brake disc, free wheel, cassette)

#### Smart MMI

Display	control monochrome
Display diagonal, resolution	2.4 inch, 240 x 320 pixels
Dimensions sMMI without dock (W X L X H)	53 mm x 85 mm x 14 mm
Connectivity	Micro-B 1.1 USB, 5 volts power supply, 500 mA connection to PC with diagnosis and parametrisation software
Mechanical / electrical contacting	twist-to-lock mechanism, corrosion-protected contacts, spring loaded
Illumination	LED backlight, 70–350 cd/m <sup>2</sup>
Display screen	Scratchproof, hardened acrylic glass screen
Weight of the sMMI (detached)	55 g

#### Smart MMI dock

Controller	23 mm inside diameter, 3 buttons (up, down, menu), hard wired
Mounting plate	Stem and handlebar mounting, angle can be adjusted in 10° steps, height can be adjusted using spacers
Weight (incl. cable and remote control)	60 g

#### Entire system

Operating temperature	- 20 °C to + 50 °C (below 0 °C recuperation or the braking assistant is automatically deactivated)
Protection rating	IP65

**NOTE:** (\*) The range depends on the battery used and on the terrain and the prevailing travelling conditions. The specified range can be achieved at optimum travelling conditions (for instance level terrain, recently charged batteries, ambient temperature of 20 °C, smooth journey etc.), a driving power of 100 watts and a pedalling efficiency of 100 watts. We reserve the right to modify the design and technology of our products to incorporate the latest developments. This operating manual is available for download on our website [www.neodrives.de](http://www.neodrives.de). If you require a version that is written in a larger font, contact our Alber Service Center.

### 3.0 smart Man-Machine Interface (sMMI)

#### Controller:

The controller fitted to the handlebar of your Pedelec is used to access menus and activate functions in the sMMI. The following functions are stored:

- Button 1 = UP (one step upwards)
- Button 2 = confirm menu key or selection
- Button 3 = DOWN (one step downwards).

#### Interface:


The sMMI is fitted to the handlebar or stem of your Pedelec. Using the controller's buttons, (Fig.3.1), you can access various functions and activate or deactivate parameters.

For some functions you can also permanently store various parameters in the software of the sMMI (see section "Firmware Updates").

Please contact your specialist dealer in this regard, they will be happy to advise.


The following provides an overview of the menu structure of the sMMI:

Fig.3.1

**NOTE:**

- **Short-cut :** Activating the Pushing Aid: push up button (1) for 2 seconds.
- **Short-cut:** Tour reset: push down button (3) for 2 seconds.

### Fitting and removing the sMMI

#### Fitting:

- Place the sMMI [4] in the correct position, (the neodrive logo faces the rider) and skewed at an angle of about 30 degrees to the dock [5] (Fig.3.1).
- Turn the sMMI [4] on the dock [5] using slight pressure by 30 degrees clockwise so that both components are in alignment, (Fig.3.2).
- The electrical connections to the controller, motor and battery pack are thereby automatically established.

#### Removing:

- Turn the sMMI [4] on its dock [5] by about 30 degrees anticlockwise.
- The electrical connections are thereby disengaged and the sMMI [4] can be removed.
- Before removing, the sMMI needs to be switched off, (see "Functions").

#### ⚠ CAUTION!

- To protect your add-on bike from unwanted use by a third party or against theft, the sMMI should always be removed from the handlebar when not in use.
- Removing the sMMI, though, does not replace securing your add-on bike against theft by other suitable means (using a bicycle lock, a safety chain or similar device).

### Functions of the sMMI

#### Switching On:

To switch on the sMMI, tap the menu button [3] on the controller. After a few seconds a welcome screen appears, followed by the start menu, which is shown to the side. If functions are already activated or the battery is not fully charged, the display of your sMMI may differ from the illustration in some parts.

#### Switching off:

To switch off your add-on bike you need to press down button [2] on the controller for about 2 seconds when the start menu is displayed. This takes you to the sub menu depicted to the side, in which display you can jump up or down using buttons [1] and [3] of the controller.

The selected field in each case is displayed within a U-shaped border.

Select the field "Turn off" and tap button [2]. Your add-on bike is now switched off.

#### Automatic switch-off:

If your add-on bike is not used for 10 minutes, the system is automatically switched off. Pressing the menu button again switches the system on once more.

#### ⚠ CAUTION!

Do not switch your Pedelec off by removing the sMMI, the electronics could be damaged by this.

Fig.3.2

Fig.3.3

Fig.3.4



Fig.3.5

**Start menu:**

As described previously, when starting up, the start menu illustrated in the accompanying image appears, (Fig.3.6).

**Explanation:**

- 6. Approximate value of the distance in km, which can be covered at the preselected assistance level (remaining range).
- 7. Indication of the current speed travelled.
- 8. Assistance level selection.
- 9. Set assistance level.
- 10. Recuperation selection.
- 11. Diverse varying information and warning indicators, see section 3.2.11 (shown here as remaining battery capacity)
- 12. Display of various travelling information.

The letter “M” to set the menu functions in the left-hand lower third of the display, is only visible when at rest. For safety reasons, different functions cannot be selected when travelling.

**⚠ CAUTION!**

- All parameters that you change will be permanently stored and are available each time you switch on the sMMI. The display of your add-on bike may therefore differ from the illustration.

**Selecting the assistance level:**

(Fig.3.6-Fig.3.7).

An assistance level that you wish to use (there are 5 possible levels) can be permanently stored by your specialist dealer as part of a travel profile in the sMMI. This is instantly available to you after switching on the sMMI and is shown in the sMMI display as a bar element [8].

You can also manually change the level of assistance at any time using the controller buttons [1] UP and [3] DOWN, whereby the number of individual white bar elements in the field [8] increases or decreases according to the level of assistance you select. When switching off the add-on bike, any changes made using the controller buttons are not accounted for. In this respect, when restarting the add-on bike, only the assistance level stored in the travel profile is available.

According to the selected assistance level, field [6] gives an indication of the distance that can be travelled using motorised support. The higher the selected assistance level, the more energy is consumed by the motor.

The range to be achieved reduces accordingly.

**NOTE:** When the add-on bike is idle please note:



- If the pushing aid is activated, you can adjust the assistance levels from a travel speed of about 8 km/h.
- Motorised support is also immediately available when activating the pushing aid by moving the pedals.

**Exceptions:**

After switching on the Add-on bike, the wheel must complete 2–3 wheel revolutions before the motorised support will start. If the pushing aid is deactivated, the motorised support is immediately available when moving the pedals.

Similarly, the assistance level can also be increased or decreased when the add-on bike is idle.

The start menu field [12] has several features. By tapping the controller button [2], various travel information can be displayed.

**Activating recuperation:**

By activating the recuperation function you can recover energy when travelling and store it in the battery. This is possible or useful from a speed of 15 km/h. Activating recuperation and adjusting is carried out using buttons [1] and [3] of the controller.

- One white bar element [10] means: 50 % energy recovery in the battery (factory setting, configurable)
- Two white bar elements mean: 100 % energy recovery in the battery (factory setting, configurable)

Depending on the battery being used and the speed, at 100 % energy recovery a maximum of 6A–8A is recovered.

If recuperation needs to be deactivated, this is also performed using button [1] of the controller.

During energy recovery “0 A” is displayed in field [12] and the charging process [a] symbol in field [11]. If energy cannot be recovered due to a battery that is already charged above 90 %, recuperation levels cannot be selected (shown in the display by symbol [b]). As soon as the battery is partially discharged, recuperation can be switched on again (shown in the display by symbol [a]).


**Automatic recuperation (optional)**

Depending on what features the manufacture has equipped the Add-on bike with, automatic recuperation can also be activated by pressing the rear wheel brake. By doing this, every braking process automatically stores energy in the battery. To ensure controlled and safe braking, recuperation during braking is set at 40 % energy recovery.

**NOTE:**

- Recuperation can only be activated at temperatures above 0 °C. If the temperature drops below 0 °C, recuperation is automatically deactivated.
- At speeds of less than 15 km/h the motor is not at its optimum operating point, which is why recuperation cannot be activated.
- Activating recuperation is not possible with a fully charged battery because this could otherwise be damaged as a result of overcharging. Once at a battery capacity  $\leq 90\%$ , recuperation can be activated.

**Fig.3.8**

**Braking assistant:**

Before travelling, if you wish, the braking assistant can be activated. This supports you when braking and ensures energy recovery to the battery (if this is not charged above the limit of 90 % or the temperature is not below 0 °C).

To activate braking assistance, in the start menu (see section 3.2.2), press button [2] on the controller to move to the next sub menu. There, select “menu” and in the next selection screen that appears, select “braking assistant”.

When selecting “on”, a new menu item opens in which you can preselect using buttons [1] and [3] of the controller, any speed from which motor-supported braking should take place.

Settings between 10 and 25 km/h are possible. Once you have set the required speed, go back to the start menu by pressing button [1] of the controller several times.

If you wish to deactivate the braking assistant, repeat the process and select “off”

instead of “on”. Please ensure that when the braking assistant is deactivated, no automatic braking or energy recovery to the battery takes place – except when selecting manual recuperation (see section 3.2.5).

**Important information:**

When the battery is fully charged, the braking assistant can in fact be activated, however, not actively used. This is only possible at a battery capacity of < 90 %. For this reason, in field [11] of the display, the symbol of the fully charged battery [a] is shown instead of the symbol for the activated braking assistant [b].

The sMMI software now waits until the battery is partially discharged, then automatically switches on the braking assistant and changes to symbol [b].

Important information on how the braking assistant works

If you have pre-set, for instance, 20 km/h, the system maintains this maximum speed irrespective of higher or lower gradients, provided that the incline is also steep enough to achieve the set maximum speed. The drive achieves this until it reaches the maximum motor torque.

If this is exceeded, the braking action gradually subsides and you need to brake yourself so as not to get faster.

Whilst the motor regulates the vehicle speed under the conditions described above, power is fed into the battery, charging it.

As soon as the pedals are pushed, the braking assistant automatically deactivates. It reactivates if the pedals are not pushed any more and therefore no force is acting on the chain or on the force sensor in the wheel hub. However, the assistant is only effective if, once adjusting the pedalling, the speed is not more than 25 km/h. By manually braking the assistant can be taken back to the speed window where it automatically activates.

**Fig.3.9**

(a)

(B)

### Pushing aid:

Should you require motorised support when pushing, for example on steep hillsides, then you can activate the pushing aid. To do this, while you are in the start menu (see section 3.2.2), press button [2] on the controller to move to the next sub menu. There, select “menu” and in the next selection screen that appears, select “pushing aid”. You can activate or deactivate the pushing aid in the next window that opens. Once you have set the required function, go back to the start menu by pressing button [1] of the controller several times. An activated pushing aid is shown in the start menu by the symbol [13].

### Using the pushing aid:

#### For add-on bikes:

The pushing aid is used by pressing button [1] for a pushing aid forwards or button [3] for a pushing aid backwards. This turns on the motor and moves your add-on bike at a maximum speed of 6 km/h (forwards) or 4 km/h (backwards) for as long as you are holding down the respective button ([1] or [3]). This is shown in field [13] as a white bar element.


#### For add-on bikes:

If buttons [1] or [3] are released, the motor switches off. It can only cut in again if add-on bike is idle. If the pedals are moved during pushing, the sMMI automatically switches to the mode of the pre-set assistance level (see section 3.2.4) so that the selection activated for this level is available from a speed of about 8 km/h. If the pedals are not moved, the sMMI switches back to the pushing aid below a speed of 8 km/h. If the sMMI is switched off, the activated pushing aid is retained in the program and is immediately ready for use upon start-up. However, after start-up the drive wheel must complete 2 to 3 revolutions before the motor moves the add-on bike when pressing button [1] or [3]. The maximum speed of the pushing aid can be set according to your requirements by your specialist dealer.

### Date and time:

The date and time can also be set individually. The set time is shown in field [12] of the start menu whereas the date is only used for internal calculations in the sMMI. As already seen in the functions described in previous sections, you can also set the date and time in the known way from the start menu in the various sub menus (see also the overview in section 3). You can modify the required parameters here.

Fig.3.10

**Travel modes:**

There are three travel modes stored in the sMMI – BOOST, TOUR and ECO.

In Eco mode the torque and therefore the maximum available motor power is automatically reduced by about 40 %, thereby reducing power consumption. At the same time the agility is also modified so that the drive responds better. The Eco travel mode is particularly suited to tours in which the battery charge needs to last as long a distance as possible.

In Tour mode, 75 % of the maximum motor torque is available. Performance and range are both at a high level. Similarly in this mode the heat generation in the motor (see section 4 Thermal management) is moderate so that long, steep inclines can be accomplished very well in Tour mode.

In Boost mode the full drive performance is accessed. It is suited to quick city travel, including moving off powerfully at traffic lights. The full “Boost performance” is not available all the time under certain circumstances. In difficult conditions, for example steep ascents, the drive power may be reduced because of heat generation (see section 4 Thermal management). The range in Boost mode is also less than in the other modes.

To activate the required mode, in the start menu (see section 3.2.2), press button [2] on the controller to move to the next sub menu. There, select “menu” and in the next selection screen that appears, select “travel mode”. This opens another window in which you can activate the required mode. Once you have activated the mode, go back to the start menu by pressing button [2] several times.

The Eco and Sport travel modes can only be selected when idle and not when travelling.

**NOTE:** The parameters stored in the Boost, Tour and Eco modes can be adapted to your travelling requirements. Contact your specialist dealer in this regard.

**Displaying journey information:**

(bike computer functions)

Before, during and after a journey, various values and information are displayed in field [12] and regularly saved. You can adjust this in the normal way by tapping menu button [2] on the controller to display the subsequent function.

**The following denote:**

Current time (needs a pre-setting, see section 3.2.8).

Displays the distance you have covered (needs a “tour reset”). Navigate using the controller buttons (see section 2) to the function “tour reset” (start menu – menu – tour reset; see diagram in section 3). Each time after activating “tour reset” the display returns to “zero”.

The distance covered by you since first using your add-on bike.

Your average speed per trip. The average speed is calculated from every “tour reset” and is displayed after a 10 minutes of travelling time.


The time in which you travel one more more trips; any idle time of your add-on bike is not included. The journey time is calculated from every “tour reset”.

The current power consumption in amperes.

Your driving performance in watts

The travel mode you activated (Boost, Tour or Eco).

**Fig.3.11**

**Information and warning symbols:**

As standard, in field [11], the remaining battery capacity is displayed. But depending on the travel situation, information and warnings as shown in the following may also be displayed in this field.

**Battery capacity display (standard display):**

The remaining battery capacity is shown in stages by a decreasing white bar.

**Charging process of the battery:**

1. The battery charger is plugged in and the battery is being charged (only for battery models which are charged via a second charging socket which is why the power connection to the sMMI from the battery does not need disconnecting).
2. The battery is being charged by energy fed by the motor (recuperation, see also section 3.2.5).

**Warning:**

The battery is empty. No power can be taken from the battery, motorised support for the add-on bike is no longer available. Please charge the battery as soon as possible using the supplied battery charger.

**Braking assistant active:**

During a downhill journey your Pedelec brakes automatically within the pre-set limits (see section 3.2.6).

**Service reminder:**

The service interval of your Pedelec has lapsed. Please arrange a service appointment with your specialist dealer. The display can be reset by your dealer using diagnosis software.

**Warning – temperature problem:**

For extremely long and steep ascents (primarily in Boost mode) the drive generates heat, which results in an automatic reduction in performance from a temperature from +80 °C inside the motor. As a result of the intelligent multipoint thermal management monitoring system (see section 4), the motor will never overheat. The performance is reduced to such an extent that damage is excluded. In very rare cases (e.g. heat built-up from external heat sources) the system may shut down completely until temperatures return to within the permissible operating range. When switching off the display will show the thermometer symbol.

**Warning – system fault:**

Inside the warning symbol, instead of the exclamation mark, various letters are shown (for instance the letter “B” in the accompanying diagram), which mean the following:


- B = battery fault
- C = communication fault in the system
- M = motor fault
- R = controller fault

When a fault occurs, motorised support is not available. Please contact your specialist dealer.

**Warning – general fault:**

A fault has occurred in the system, motorised support is no longer available. Please contact your specialist dealer.

**Fig.3.12**

**Fig.3.13**

### Easy Display mode:

By activating Easy Display, the following graphical symbols presented in section 3.2.2 are enlarged on the display.

- indication of the current speed travelled
- indication of the daily kilometres travelled (requires prior reset, see section 3.2.10)
- indication of the remaining battery capacity including indication of the expected remaining range

The relevant display can be selected by tapping button [2] on the controller.

If buttons [1] or [3] are tapped, the display shows the assistance level for about 3 seconds that can be increased or decreased using the two buttons.

### Activating Easy Display:

From the start menu you go via Menu – Settings – Easy Display to display “on/off”.

If Easy Display is activated, the sMMI automatically goes into this mode every time it is switched on.

### Deactivating Easy Display:

To deactivate the mode you need to press button [2] on the controller for about 2 seconds. Via Menu – Settings – Easy Display you arrive at the display “on/off”.

**NOTE:** When Easy Display is activated the functions “pushing aid” and “recuperation” are not available.

### Other settings:


As standard, your sMMI has a German user interface. However, if required, an English / Dutch / French interface can be activated. The software status of the sMMI, the motor and the battery pack can also be retrieved from the menu item “information”. Language activation or status requests are performed similarly to the procedures described in the previous sections via various menus and sub menus.

### USB port: (Fig.3.14)

The sMMI has a USB port [14], which is primarily used by your specialist dealer to carry out maintenance and diagnosis work. Ensure that the rubber cover is always fully inserted and seals the port. If the sMMI is not fully sealed, moisture can get inside or fog the display from the inside.

In addition to English we can also set Dutch and French. However, if required, an English / Dutch / French interface can be activated.

Fig.3.14

### Programming options by the specialist dealer:

Your specialist dealer can tailor the travel characteristics of your add-on bike to your requirements using a diagnosis and programming software. Generally, the factory settings are ideally tailored to the respective model and do not need to be changed. The factory settings may vary depending on your bicycle model (mountain bike / trekking). Nevertheless, if adjustments need to be made, the following parameters can be programmed:

<b>Speed of the pushing aid forwards:</b>	Set to 4 km/h as default.
<b>Standard assistance level:</b>	The assistance level that is always available when the Pedelec is switched on. Can be set in the range 0–5, the value 3 is pre-programmed as default.
<b>sMMI standard settings:</b>	Language settings, time format (12/24 h)
<b>sMMI lock:</b>	Optionally, the sMMI can be permanently connected to the motor. This results in the sMMI only being operational with this one motor/system and cannot be used with another motor.
<b>Service interval:</b>	Depending on the specified suitable maintenance interval, this can be set or reset by date or kilometres travelled (depending on which comes first).
<b>Wheel circumference:</b>	To correctly display the speed and to comply with statutory requirements with regard to speed limitation, the bicycle dealer can change the wheel circumference. This is only necessary if the drive wheel is subsequently fitted with a tyre which increases or reduces the original wheel circumference, or the motor has been reassigned to a larger or smaller rim.

### WARNING!

- This parameter may only be changed if the statutory requirements are complied with (max.25 km/h for add-on bikes) and in the event of misuse this leads to the warranty and product liability being void.
- Similarly, unlawful tampering can have criminal consequences if investigated by the police.
- All changes to the travel parameters are recorded in the sMMI data storage device.

### Firmware updates and relaying them to the motor and battery:

As part of product maintenance and to expand the range of functions, firmware updates are available to download to specialist dealers from time to time. Your specialist dealer will be happy to advise you. If you have loaded an update to your sMMI from your specialist dealer, the new firmware will install when you next start up your add-on-bike.

To do this, proceed as follows:

- Place the sMMI on its dock [5] as described in section 3.1.
- The connection to the motor and battery is automatically established and instead of the start menu, a warning note appears “firmware update is running”. A bar chart also shows the progress of the data transfer.
- Once the data is fully transferred, the screen display automatically changes to the start menu (see section 3.2.2).
- Check the settings stored by you in the sMMI if applicable. This may have changed as a result of the update.

### CAUTION!

Do not interrupt the update process, for example, by removing the sMMI from its dock. This may lead to damage.

**Note on the sMMI plugs: (Fig.3.15)**

If you need to remove the sMMI plugs leading to the battery at any time, please note the following information when connecting them again:

- The plug of the sMMI and the socket on the additional cable to the battery must be correctly aligned when connecting.
- Both parts are a little tricky to connect but this prevents moisture ingress during operation.
- Do not kink the attached cables at any time when connecting the plug and the socket. There is a risk of the cable breaking.

**Thermal management:**

A combination of three temperature sensors, an intelligent software control system and patented air circulation ensure the motor is optimally cooled. In practice this means greater and longer output on climbs or under high (attachment) loads.

Advantage: protection against premature overheating on long climbs and high loads – longer support on hills, higher level of efficiency and thereby a lower battery consumption as the motor is optimally cooled.

About the theory


As with all drives, gearless wheel hub motors are also optimised to an operating point consisting of speed, load and output. Our wheel hub motors are designed for operation in the speed range between 15 km/h and 25 km/h and a normal drive output of 250 watts. In this speed and performance range they achieve the greatest efficiency and range, which means that the supplied energy is optimally converted into drive energy.

Whenever a motor is operated outside of the optimum operating point, its level of efficiency decreases. This leads to the fact that the energy is no longer optimally converted, rather part of the supplied energy is converted into heat. The range therefore decreases and the heat needs to be discharged. In the neodrives motors this heat discharge is achieved via a large contact area inside the motor (stator carrier) to the dropout or chain stays of the bicycle frame. In addition, cooling ribs inside and outside the drive housing ensure the greatest possible heat exchange with the surroundings. Heat which cannot be discharged causes the drive motor to heat up.

The neodrives wheel hub motors monitor both the supplied energy and the temperatures in the motor. This enables damage to be prevented as a result of overheating in an overload situation. However this also results in the motor performance that is accessible by the rider. Being reduced to prevent overheating.

**Fig.3.15**

If a temperature of 80 °C is exceeded in the motor electronics, the motor control system reduces the input power and thereby the support. This means that the higher the temperature increase in the motor, the less drive output can be accessed and the lower the support available. When the motor cools down, the energy input is increased again and the drive output rises. Important: The motor cannot be damaged by heat build-up. The temperature symbol (Fig.3.16) only appears when the power is completely reduced.

**Fig.3.16**

This regulation of the drive output as a function of the motor temperature is progressive so that there is always support but the motor is not damaged by overheating.

In practice

As a result of the points explained above, the day-to-day practice is dependent on the outside temperature, total weight, incline, the terrain, cadence, air pressure and speed. These factors may lead to a temperature being reached which causes the output or support to be reduced.

However, this does not mean a fault or failure of the drive, you can continue to cycle at lower support. In an extreme case it may cause a brief complete shutdown.

Extreme example: An incline of 10 – 12 % at an altitude of 500 metres, a total weight of 120 kg, loose terrain, a maximum assistance level, a speed of < 10 km/h and a cadence of 60 rpm mean operation in an unfavourable range at low efficiency and range at simultaneously high heat generation. This will result in a reduction in the drive output.

Tip: Ideally, by selecting a lower gear with higher cadence, the travel mode "Tour" or "Eco", a reduced assistance level and/or a short break (in which the drive can cool down again), you can continue to cycle.

**⚠ CAUTION!**

The motor must never be “forcibly cooled” externally with water. This can lead to damage and does not help especially with cooling as it is primarily the inside of the motor that is hot.

**Motor:**

The drive wheel of your Pedelec can be removed from the bicycle frame at any time, for example, for cleaning purposes or in case of a puncture. Proceed extremely carefully when removing and subsequently installing the drive wheel and note, in particular, the instructions and information of the manufacturers of the various components attached to the wheel, especially the brake disc.


**NOTE:** For clarity reasons the following diagrams only show the drive motor integrated in the wheel but not the entire drive wheel.

**Removing the drive wheel:**

Note or mark the cable routing and the fixing points of the cable ties before removing the drive wheel. First loosen and remove all cable ties which are securing the cables [16] coming from the motor and cables and leads of other components to the bicycle frame. Then disconnect the plug [15] on the cable of the motor [16] from the socket [17] on the battery cable. Loosen the two nuts [18] or the quick-release mechanism with which the wheel is fixed to the frame so that the entire wheel can be removed from the frame of your add-on bike.

**⚠ CAUTION!**

- Note or mark the position of the torque support [20]. This must be re-fitted in exactly the same position as it was before removal, when you later attach the wheel.
- Never hold or transport the wheel you have removed by the cable, [16] coming from the motor. There is a risk of the cable breaking.

**Fig.3.17**

### Attaching the drive wheel:

Make sure that all the components attached to the wheel have been installed in accordance with the instructions and specifications of the respective manufacturer. This relates in particular to the brakes and the gear box. Do not forget to reinstall the torque support [20] in the same position as it was taken from.

Then push the wheel into the slot of the frame and tighten it using the hub axle nuts [18] in the following sequence:

- first tighten on the side of the gear box (diagram A)
- then on the side of the brakes (diagram B)

In each case, the tightening torque of the two nuts is between 30 and 40 Nm.

Ensure that the tooth lock washer [19] is located beneath the axle hub nut [18] otherwise there is the risk of the axle hub nut [18] coming loose.

If your wheels are equipped with a quick-release mechanism, please observe the manufacturer's specifications on installation and the tightening torque.

If the wheel is correctly attached to the frame, the motor with the end of the cable that leads to the battery can be connected. In doing so, make sure the plug [15] is correctly aligned with the socket [17]. The curved surfaces ([A] and [B]) must be aligned with each other. Then fasten all cables and leads with cable ties to the bicycle frame and perform a final function test.


#### **⚠ WARNING!**

- Make sure the cables are laid correctly. Failure to do so could lead to the cable getting caught in the brake disc, the drive or the spokes and by locking the wheel cause a fall.
- Never install the motor without the torque support [20]. This could cause a total loss (twisting off the cable). This voids all warranty or guarantee claims.

#### **⚠ CAUTION!**

- It is essential to observe the instructions and specifications of the manufacturer of the various components attached to the wheel in all your installation work. This particularly applies to the brakes, the gear box and quick-release mechanism.
- Carry 5 cable ties with you along with your repair tool, to securely fasten any loose cables during a trip.
- Installing or removing the drive wheel is best carried out when the add-on bike is upside down (standing on its crank and connecting frame). Remove the sMMI attached to the handlebars first so that it does not get damaged.
- Always use the original gear cluster built by the bicycle manufacturer. Using other makes may lead to restricted function or to the gear cluster scraping the chainstays

**Fig.3.18**

## Cleaning the motor and sMMI:

When cleaning the device, never use cleaning benzine, thinner, acetone or similar agents. Never use abrasive detergents or aggressive cleaning agents. Instead, use only conventional household cleaning agents and disinfectants (isopropyl alcohol).

### Motor:

- Your add-on bike motor should be regularly cleaned of dirt, ideally with a dry brush or a damp (not wet) cloth. Never clean the motor under running water, for example a garden hose or even a high-pressure water jet. However, travelling in the rain and on wet lanes is entirely possible.
- Water ingress can destroy the motor. Therefore always ensure that liquid or moisture does not get into the motor.
- Do not clean the motor when it is warm, for instance, immediately after a trip. Wait until it has cooled down. Otherwise this could cause damage.
- If the motor is dismantled, for example for cleaning purposes, it must never be held or transported over the cables which would run the risk of the cable breaking.
- If the motor has been removed from the add-on bike (see section 4.1), the plug from the motor and the socket of the cable to the battery pack need to be checked for any contamination or cleaned before connecting.

### sMMI Unit:

- The contacts of the sMMI dock are spring loaded and from time to time cleaned with a contact spray to ensure they function perfectly and for a long time.
- The sMMI housing may only be cleaned using a damp (not wet) cloth.

uninterrupted vehicle operation or downtime as a result of continuous, direct sunlight)  
cause system downtime, leave the motor to cool down for about 10 minutes before continuing your journey.

- The maximum speed (non-motorised operation) of the system is 75 km/h. Exceeding this speed compromises the electronic components which could be damaged in the worst-case scenario. The maximum speed is logged by the system.

## Transportation:

The following information needs to be noted when transporting the Pedelec by car.

- Protect all components of your Pedelec against moisture and dirt by suitable means.
- Remove the battery and the sMMI from the bicycle before you attach the Pedelec to the bicycle rack of your car. This also reduces the weight that you need to lift, particularly when using a roof rack system.
- Always transport the battery and the sMMI inside your car.
- Even when transporting inside the car (e.g. in an estate car), the sMMI and battery should be removed to avoid damage occurring when loading and during the journey.
- If you have a roof carrier system with down tube clamps ensure that when tightening the clamping device, the battery mounting slide rail is not crushed/damaged.
- Ensure that the ends of the cables cannot cause damage to the Pedelec or your car during transportation.
- After reaching your destination check all the Pedelec's contacts for possible foreign bodies or moisture. To ensure reliable function all plug connections in particular must be free from dirt and foreign bodies and be completely dry.
- Never lay your Pedelec on the gear box side during transportation, for example in a car boot. This could be damaged.

## Safety precautions:

- When not in use, do not expose your add-on bike to strong sunlight for long periods of time. This could cause the motor to heat up and in an extreme case, not be able to operate at full power. Even plastic parts age quicker under intense sunlight.
- If increased temperatures (caused for instance by uninterrupted vehicle operation or downtime as a result of continuous, direct sunlight), cause system downtime, leave the motor to cool down for about 10 minutes before continuing your journey.
- The maximum speed (non-motorised operation) of the system is 75 km/h. Exceeding this speed compromises the electronic components which could be damaged in the worst-case scenario. The maximum speed is logged by the system.

**Error indications and possible remedies:**

The system cannot be turned on.	<ul style="list-style-type: none"> <li>• Is the battery correctly inserted in it's holder?</li> <li>• Are all the plugs connected correctly?</li> <li>• Are there ant deposits, (e.g.metal filings), on the magnetic plug on the battery?</li> <li>• Note: Please check this carefully.</li> <li>• Is the battery activated?</li> <li>• After 48 hours of non-use the battery goes into standby mode and needs to be reactivated by pressing the battery button.</li> <li>• Has the sMMI lock been activated during servicing?</li> <li>• Note: The sMMI only works with the designated motor, (see section 3.2.14).</li> <li>• Do the contacts of the sMMI spring back properly?</li> <li>• Check this by using your fingers and pressing each each pin individually into the dock. Each pin should spring back. Treat sticking pins with contact spray.</li> </ul>
The battery will not charge.	<ul style="list-style-type: none"> <li>• Are there ant deposits, (e.g.metal filings), on the magnetic plug on the battery charger?</li> <li>• Note: Please check this carefully.</li> <li>• Is the ambient temperature higher than 0°C?</li> <li>• Note: Below 0°C the battery charger will not charge the battery. Always charge the battery at room temperature.</li> <li>• Follow the detail of the charging process, in particular the fault codes in the charger operating manual.</li> </ul>
No Motorised Support. (sMMI in operation, but motor support unavailable)	<ul style="list-style-type: none"> <li>• Check the motor cable and plug are correctly aligned, (section 4.1).</li> <li>• Does the display show an error message?</li> <li>• If so, follow the relevant recommendations in section 3.2.11.</li> <li>• Was the start up routine observed for the light? (section 3.2.3).</li> <li>• Is the system permanently in recuperation mode?</li> <li>• If so, check the switch on the rear wheel brake lever is fitted correctly, (sMMI with brake cable only).</li> <li>• Is the sMMI correctly fitted on the dock, (section 3.1).</li> </ul>
Recuperation / hill climbing system is not working.	<ul style="list-style-type: none"> <li>• Is the battery capacity more than 90%?</li> <li>• Note: Recuperation only works when the battery is at 90% or more.</li> <li>• Is the current travelling speed less than 15 km/h?</li> <li>• Note: Recuperation only works above 15 km/h.</li> <li>• Is the current travelling speed more than 28 km/h?</li> <li>• Note: Recuperation is not available above 28 km/h.</li> </ul>
The assistance levels cannot be changed whilst standing idle.	<ul style="list-style-type: none"> <li>• You have activated the pushing aid in the menu.</li> <li>• As soon as you step onto the pedals you will be able to select the assistance levels.</li> <li>• Alternatively, you can deactivate the pushing aid via the menu, (section 3.2.7).</li> </ul>
The motor is not achieving maximum output.	<ul style="list-style-type: none"> <li>• It is possible that the scooter is in high temperature range.</li> <li>• When the motor becomes hotter than 80C, performance is gradually reduced.</li> <li>• Leave the motor to cool down for about 10 minutes in a cool shaded area, then try again.</li> <li>• Battery discharged.</li> <li>• As the battery voltage decreases the performance/maximum speed also drops.</li> <li>• The speed difference between a fully charged battery and a nearly discharged battery can be in the region of 2-3 km/h.</li> <li>• Charge the battery as soon as you can.</li> </ul>
The display shows the symbol for the service reminder, (section 3.2.11).	<ul style="list-style-type: none"> <li>• You may continue to use your bike normally, but please arrange for your authorised dealer to service it at your earliest convenience.</li> <li>• After servicing they will reset the display.</li> </ul>

**Part 2 Neodrives Battery**

**4.0 Introduction**

**Introduction:**

**Important instructions :**

 Along with this operating manual, your add-on bike comes with additional documents. Please observe the specifications and instructions in these documents.

**Proper use of the neodrives battery:**

The neodrives battery is used exclusively to power the neodrives drive system. No other components may be connected to it. Any other use requires the written authorisation of the manufacturer.

This operating manual only describes the use of the neodrives battery and corresponds to the state of the art at the time of print.

The manufacturer reserves the right to make changes resulting from further development of the mechanics, software or legal requirements. The manufacturer regards uses including the following as misuse of the battery:

- Use of the battery that contravenes the instructions and recommendations in this operating manual.
- Exceeding the technical limits laid down in this operating manual.
- Technical changes to the battery.
- Changes to the software of the battery.
- Unauthorised attachments or use of the battery.

The manufacturer shall not be liable for any damage caused by improper use of the battery.

 Before using the device, carefully read all safety and hazard information contained in the individual chapters of this operating manual and all other enclosed documents.

**Signs and symbols:**

Important information for your safety are identified in this operating manual as follows:

Word	Definition
 <b>DANGER!</b>	Advice to the user of potential risk of serious injury or death if the advice is not followed
 <b>WARNING!</b>	Advice to the user of a risk of injury if the advice is not followed
 <b>CAUTION!</b>	Advice to user that potential damage to equipment may occur if the advice is not followed
<b>NOTE:</b>	General advice or best practice
	Reference To Additional Documentation
	Indicates magnetic forces

Observe these indications at all times to avoid injury to people and damage to the product.

**Permissible operating conditions / operating sites:**

 **CAUTION!**  
The battery must only be used at ambient temperatures between -20 °C and +60 °C.

 Observe the information on permissible operating conditions in the add-on bike manufacturer's operating instructions. The manufacturer's restrictions on the permissible operating conditions must be observed. Observe the safety and hazard information provided in the individual chapters of the operating manual.

**Standard scope of delivery (neodrives components):**

- Battery including 1 key pair.
- Battery slide rail to fit the battery (already installed on the add-on bike).
- This operating manual

5.0 Technical Data

Category	ICR18650MG1
Battery type:	Lithium ion
Rated capacity:	14.5 Ah
Rated voltage:	36.2 V
End-of-charge voltage	42 V
Total energy:	515 Wh
Maximum discharge current:	30 A
Ambient charging temperature:	0 °C to 40 °C
Ambient operating temperature:	-20 °C to 60 °C
Ambient storage temperature:	+10°C to +30°C
Number of cells:	50
Protection rating:	IP54
Weight:	Approximately 3.5 kg

**NOTE:** We reserve the right to modify the design and technology of our products to incorporate the latest developments. Please retain this operating manual for future reference. This operating manual is also available to download from our website [www.neodrives.de](http://www.neodrives.de).

**Key Elements Overview:**

**Battery**

- Battery (housing) 1
- Key 2
- Charger socket / motor connection 3
- On/off button 4
- LED display 5

**Motor**

- Motor cable connector 7
- Motor 8


**Charging device**

- Charging cable connector 9

**On the add-on bike**

- Battery slide rail 6

**Fig.5.1**

## 6.0 Safety Instructions/Precautions

### WARNING!

- Read and observe the following safety instructions and precautions before activating the battery and before starting the charging process.
- Failure to comply with the safety precautions and instructions may damage the product or result in electric shock, fire and/or serious injuries.
- The lithium-ion battery contains chemical substances, which may cause hazardous reactions if the safety instructions specified here are disregarded.
- The manufacturer accepts no liability for damage resulting from non-compliance with these instructions.

### Safety instructions and precautions on the use of the battery:

- Before using for the first time, the battery should be fully charged.
- The battery must only be used at temperatures between -20 °C and 60 °C.
- The battery must not be exposed to heat (e.g. radiators) or fire. External heat exposure can lead to the battery exploding.
- In the (unlikely) event of the battery overheating or catching on fire, you must ensure that the battery does not come into contact with water or other liquids. The only suitable extinguishing agent that is recommended by the cell manufacturers is sand.
- Your add-on bike uses power in all operating modes. Therefore, if possible, charge the battery after every use.
- The battery must only be used to supply energy to the neodrives components. Any other use requires the written authorisation of the manufacturer.
- The battery must not be opened or taken apart. Improper opening or wilful destruction of the battery involves the danger of serious injury. In addition, opening the battery voids any warranty claim.
- Never connect the battery contacts in the socket [3] to metallic objects, or ensure that the contacts never come into contact with metallic objects (for example with metal filings).
- If the socket [3] is dirty, clean it with a clean and dry cloth.
- Never immerse the battery in water.
- The useful life of the battery depends, among other things, on its storage location. Consequently, never leave the battery (irrespective of whether the battery is installed in or has been removed from the add-on bike) in hot locations for prolonged periods. In particular, use the boot of a car parked in the sun only for transport rather than as a storage location.
- The battery must not be subjected to mechanical impact. If, for example, the add-on bike overturns and the battery thereby hits the ground, the battery must be checked by the manufacturer. Contact your specialist dealer in this regard. You must not continue to use a damaged battery.
- If damaged or defective, the battery must be singled out and checked. Please contact your specialist dealer and discuss with them the next steps with regard to return and repair. The defective/damaged battery must not be used again or opened.
- Always ensure that the battery is kept clean and dry.

### Safety instructions and precautions on how to store the battery:

- Protect the battery immediately upon separating from the battery charger or motor. Never allow any moisture or foreign particles (e.g. metal fragments, small nails, filings or other conductive metals) to get into the battery.
- Do not expose the battery to moisture of any kind during storage (water, rainwater, snow, etc.).
- Before storing it, charge the battery and check its charge status every 3 months.
- Store the battery in a cool and dry location where it is safe from damage and unauthorised access.
- To achieve the optimum battery service life, store the battery at a temperature from 18°C to 23°C and at a humidity of 0 to 80 per cent. The charge status under these conditions should be 70 per cent.
- Check the charge status of the battery every 3 months during storage and charge it to 70 per cent if necessary.

### Safety instructions and precautions on the charging process:

- Only charge the battery in a ventilated, dry and dust-free environment.
- Do not charge the battery in the presence or vicinity of flammable liquids or gases.
- Do not expose the battery to moisture of any kind during charging (water, rainwater, snow, etc.).
- Do not carry out the charging process in rooms where moisture may affect the battery.
- The battery must only be charged at temperatures between 0 °C and 40 °C. If you attempt to carry out a charging process outside of this temperature range, the battery mechanism automatically switches off the charging process. The battery reaches its maximum service life if it is charged at temperatures between 10 °C and 30 °C.
- Only use the designated charger to charge the battery. Your specialist dealer can provide the relevant information.
- Using an unsuitable battery charger can lead to malfunction and result in the battery having a limited service life. There is also a danger of fire and explosion.
- When the charging process is complete, disconnect the charger from the mains socket first before disconnecting it from the battery.
- Ensure adequate air circulation as soon as the battery is charged. In principle, only charge the battery under supervision.
- Damaged batteries must not be recharged or used any further.
- Damaged battery chargers (damage to the plug, housing, cable) must not be used.

### Safety information and precautions on how to transport and ship the battery:

Lithium-ion cells are used in the neodrives battery. Transport and shipping of the battery are, therefore, subject to all applicable statutory requirements, which must be strictly observed. For instance, a defective battery must never be transported by plane.

If your battery is defective, take it personally to your specialist dealer as shipping lithium-ion batteries by post or other carriers is subject to strict rules and regulations. Once again, we recommend that you contact your specialist dealer.

As transport regulations may change annually, we strongly recommend that you consult your travel operator, airline or shipping line before you set out on your trip in order to learn about the current applicable regulations. A defective battery must not be taken on a plane or put in your luggage.

If your battery is fitted to the add-on bike when being transported, relaxed transport regulations apply in accordance with UN3171.

#### NOTE:

- Make sure you keep the packaging container of the battery in case you need to transport it.
- Discuss transportation with your specialist dealer prior to shipping.

### Safety information and precautions on the charger:

Before starting the charging process, read and follow all instructions and precautions enclosed with the charger and the following precautions and safety information.

- Only use the designated charger to charge the battery. Your specialist dealer can provide the relevant information.
- Using an unsuitable battery charger can lead to malfunction and result in the battery having a limited service life. There is also a danger of fire and explosion.
- The charging process ends automatically as soon as the battery is fully charged. This avoids overcharging.
- When the charging process is complete, we recommend disconnecting the charger from the mains socket first before disconnecting it from the battery.
- Never use a charger that is not recommended by your specialist dealer.
- Do not expose the charger to moisture of any kind during charging (water, rain water, snow, etc.).
- Do not carry out the charging process in rooms where moisture may affect the charger.
- Be aware of condensation. If the charger is brought from a cold into a warm room, condensation may form. In this case, refrain from using the charger until all condensation has dissipated. Please note that this may take several hours.
- Never carry the charger by its power cable or the charger cable.
- Never tug at the power cable to disconnect the charger from the mains socket.
- Never subject cable and plug to any pressure. Overstretching or bending the cable, pinching a cable between a wall and a window frame or placing heavy objects on a cable or a plug may result in electric shock or fire.
- Lay the power cable and the attached charging cable so that nobody can step on it or trip over it and so both cables are protected against any other harmful effects or stress.

- Do not operate the charger if the power cable, the charging cable or the plugs attached to the cables are damaged. Damaged parts must be replaced immediately by the authorised specialist dealer.
- Do not use or disassemble the charger when it has received a hard blow or was dropped or damaged in another way. Take the damaged charger to a specialist dealer who has been authorised to perform repairs.
- The charger must not be used by children.
- Never attempt to disassemble or modify the charger.
- Do not cover the charger during the charging process or place any objects on top of the device.
- Never connect the terminals of the charging plug with any metal objects.
- Ensure that the plug is firmly inserted in the socket.
- Never touch plugs with wet hands.
- Do not use the plug of the charger and/or the mains plug if they are wet or dirty. Before inserting it, clean the plug using a dry cloth.

## 7.0 Operation

### Starting up:

#### Information on the operating modes:

Essentially, the battery has two operating modes. It is either in "Active mode" or in "Deep Sleep mode".

In Active mode the battery consumes at least 5 mA per hour (consumption of the battery's own electronics).

To keep the battery's own consumption as low as possible, the battery automatically switches to the so-called Deep Sleep mode after 48 hours.

#### Inserting the battery:

- Place the battery [1] on the battery slide rail [6] mounted on the add-on bike.
- Push the battery [1], as shown in the diagram, up to the front edge of the battery slide rail [6].
- Lock the battery [1] by carefully turning the key [2] clockwise until it stops. The battery can no longer be removed from the battery slide rail.
- Remove the key [2] from the battery [1].

#### Connecting the battery to the motor cable:

- Insert the plug [7] of the cable coming from the motor into the socket [3] on the battery [1].
- The two parts are correctly aligned and interlocked automatically by means of a magnetic latch.

Before inserting the plug [7] in the socket [3], ensure that both parts are clean and there are no metallic particles on them. If you detect such particles, use a dry, clean cloth to remove them.

#### Switching on the battery:

If the battery has been used within 48 hours it does not need to be switched on. The Add-on bike is ready for use and can be switched on via the sMMI and used.

If the battery is first being activated now, or if it has not been used for more than 48 hours (Deep Sleep mode), it needs to be switched on.

- Tap button [4].
- All the LEDs [5] flash three times to show it is switched on.
- Your add-on bike is now ready for use and can be switched on via the sMMI and used.

#### **CAUTION!**

If the motor cable is not yet connected to the battery, the battery will still be in Active mode when switched on.

If the battery cannot be switched on, the cell voltage may be too low. In this case connect the charger and then press the on/off button [4]. The battery will be charged for a minute.

Fig.7.1


Fig.7.2



Fig.7.3


**Battery capacity indicator:**

You can check the battery capacity at any time using the LED display.

If the battery has not been used for more than 48 hours:

- Tap button [4].
- The battery is switched on, all LEDs (a – e) flash three times.
- Tap button [4] again.
- Now the LEDs indicate the capacity of the battery as shown in the following table.
- If the battery has been used in the past 48 hours:
- Tap button [4].
- Now the LEDs indicate the capacity of the battery as shown in the following table.

**Fig.7.4**

LED lights Steadily	LED Flashes	Battery Charge
-	a	- 19 %
a	-	20 - 39 %
a, b	-	40 - 59 %
a, b, c	-	60 - 79 %
a, b, c, d	-	80 - 99 %
a, b, c, d, e	-	100 %

### Detaching the battery:

#### Switching off the battery:

The battery is switched off using the sMMI controller (see the sMMI and motor operating manual). This puts the battery initially in Active mode for 48 hours. This means that, within this period, the sMMI can be reactivated at any time without having to first switch on the battery. The power required for this is minimal.

#### Disconnecting the cable:

- Disconnect the add-on bike on the sMMI.
- Then remove the plug of the motor cable [7] from of the socket [3].

#### ⚠ CAUTION!




Ensure that the plug of the motor cable [7] does not come into contact with any metallic particles when being set down (risk of contamination).

#### Removing the battery:

- Insert the key [2] in the lock on the battery [1].
- Carefully turn the key [2] anticlockwise until it stops. The lock is now released, the key can no longer be removed from the battery.
- Pull the battery [1] upwards about 20.0 mm along the slide rail [6] and then remove it completely.
- Store the battery in a clean place.

#### ⚠ CAUTION!

Ensure that the socket [3] does not come into contact with any metal

**Charging the battery:**

Fully charge the battery before its first use. The battery capacity upon delivery is generally 30%. The battery can be charged at any capacity without adversely impacting its service life. The battery achieves its maximum service life when it is charged at an ambient temperature of between 10 °C and 30 °C

**Connecting the battery charger:**

The battery [1] does not need to be removed from the Pedelec for charging, it can stay where it is. Only the motor cable plug [9] needs to be removed (see section 4.2). Then, proceed as before.

- Insert the plug [9] of the battery charger into the socket [3] on the battery.
- The two parts are correctly aligned and interlocked automatically by means of a Magnetic latch.
- Carry out the charging process according to the specifications in the charger operating manual. In addition, follow the instructions on the charging process given in section 2.3.

Before inserting the plug [6] in the socket [3], ensure that both parts are clean and there are no metallic particles on them. If you detect such particles, use a dry and clean cloth to remove them.

**Charging process:**


When charging the battery follow the instructions of the charger operating manual. Also observe the safety information and precautions provided in section 6.


**LED displays during the charging process:**

The following table shows the LED display [5] of the battery during charging.

**Fig.7.8**

**Fig.7.9**

LED lights Steadily	LED Flashes	Battery Charge
-	a	about 0-19 %
a	b	about 20-39 %
a, b	c	about 40-59 %
a, b, c	d	about 60-79 %
a, b, c, d	e	about 80-99 %
a, b, c, d, e	-	<b>Fully Charged 100%</b>

**⚠ CAUTION!**

If an error occurs during the charging process, all LEDs will light up. Check if all criteria (e.g. ambient temperature, properly connected charging plug, etc.) for the charging process have been met in accordance with this operating manual and the operating manual that is included with the charger. Before using the device, always check the charge status of the battery. The battery should be fully charged prior to use in order to provide motorised support at all times.

**⚠ WARNING!**

Do not leave the charger connected to the mains socket any longer than necessary. When the charging process is complete, disconnect the charger from the mains socket first before disconnecting it from the battery. The battery may only be charged in a dry room at temperatures between 0° and 40° degrees centigrade.

 Please observe the instructions in the operating manual enclosed with the charger. Observe the safety information and precautions regarding the battery provided in Part 2, Neodrives Battery, of this operating manual.

**Keys:**

The battery is supplied with two keys for locking in the battery slide rail. Your specialist dealer should note the identifier engraved on the key in the supplied add-on bike documentation so that these can be reordered if required. Please therefore check whether the key identifier is entered in the documentation. If this is not the case, please enter it. Keys can only be reordered by the specialist dealer from AXA (as of November 2013).

**Cleaning the battery:**

When cleaning the device, never use cleaning benzine, thinner, acetone or similar agents. Instead, use only conventional household cleaning agents and disinfectants (isopropyl alcohol).

- The plug of the charging cable on the battery charger [9], the plug of the motor cable [7] and the charger socket [3] on the battery may only be cleaned using a dry cloth.
- The battery must never be sprayed with a steam cleaner or similar device.

**Battery Storage:**

We recommend storing the unused battery in a dry place at temperatures between 10°C and 30°C.

- Do not keep your battery in places where the temperature exceeds 45°C or is less than -20°C.
- Avoid direct sunlight.
- Avoid storage in humid areas to prevent corrosion on the plug contacts.
- When the battery is not in use, recharge it at least every 12 weeks to prevent damage to the battery.”

**8.0 Disposal**

**Disposal:**


Electrical and electronic equipment need to be disposed of separately to general household waste at specific state-provided locations. The correct disposal and separate collection of used appliances serves to prevent potential damage to health and the environment. It is a requirement for the re-utilisation and recycling of used electrical and electronic equipment. Detailed information on the disposal of your used equipment can be obtained from your local authority, your waste disposal service, the specialist dealer from which you purchased the product, or your sales contact. This information only applies to equipment which is installed and sold in the countries of the European Union and which is subject to the European directive 2002/96/EC. In countries outside of the European Union, deviating conditions apply to the disposal of electrical and electronic waste.

**9.0 Liability**

**Liability:**

The manufacturer is not liable in any instance if the battery has been / is incorrectly handled.  
 the battery has been / is commissioned in contravention to the instructions in this operating manual.  
 the battery has been / is operated with insufficient battery charge.  
 the battery has been / is repaired or otherwise modified by a person not authorised to carry out such work.  
 the battery has been / is used contrary to the intended use.

Sunrise Medical GmbH  
 Kahlbachring 2-4  
 69254 Malsch/Heidelberg  
 Deutschland  
 Tel.: +49 (0) 7253/980-0  
 Fax: +49 (0) 7253/980-222  
 kundenservice@sunrisemedical.de  
 www.SunriseMedical.de

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