JAY Anterior Trunk Supports

EN

Important Consumer Information.

NOTE: This manual contains important instructions that must be passed on to the user of this product. Please do not remove this manual before delivery to the end user.

SUPPLIER: This manual must be given to the user of this product.

USER: Before using this product, read this entire manual and save for future reference.

WARRANTY REGISTRATION: To validate the warranty on this product, please complete and return the postcard attached to the enclosed instruction and warranty booklet.

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
Fig. 7

Fig. 8

The intended lifetime of this product is 5 years.

SHIPPING:

Check all parts are free from shipping damage. In case of damage DO NOT use. Contact Sunrise Medical or your supplier for further instructions.

⚠ WARNING!

- Jay Positioning harnesses and shoulder straps were designed for positioning, not to be used as transportation safety belt!
- JAY POSITIONING HARNESSSES, AND SHOULDER STRAPS may be used while in transit but are not crash tested or approved to provide protection during a crash.
- Any harness or support strap must not interfere with breathing.
- An appropriate clinical assessment by a clinician or qualified wheelchair professional is highly recommended to guarantee safety and correct clinical usage.

The products described in this manual were designed to be fitted and installed exclusively by:

- An authorized Sunrise Medical supplier.
- A qualified and competent healthcare professional who has undergone appropriate certified training by either Sunrise Medical or their approved supplier.

ANTERIOR TRUNK SUPPORTS INSTALLATION:

INSTALLATION REQUIREMENTS

To be able to mount any anterior trunk support / shoulder harness you have to ensure that the following are in place:

- A wheelchair with a solid backrest set approximately at shoulder height.
- An installed pelvic belt.

Installation steps (Fig. 2-6)

1. Fitting and installation should be performed with the user sitting in the wheelchair.
2. Place the desired anterior trunk support /shoulder harness on chest with upper straps over the backrest.
3. Position straps between neck and shoulder joint.

Top Straps

- i. Attach to the solid backrest using the fasteners/ strap positioners (Fig. 1, sold separately from harness).

Lower Straps

- ii. Identify optimal position for the lower straps on the backrest without interfering with other positioning devices (lateral supports, etc.).
- iii. Install end mounting between the top of the pelvis and the seat.
- iv. Attach using the fasteners supplied.

Installation steps:Trunk Strap Style (Fig. 7)

- Mounts to backrest with strap anchors.
- Anterior Velcro and "D" ring closure provide quick user access.

Installation steps: Y-Style Harness (Fig.8)

1. Place the upper part of the Y-Strap above and below the rotated shoulder and install belt anchors to the backrest's rear face.
2. Stretch the lower strap diagonally across the trunk and attach the lower part of the Y-Strap to the backrest's rear face.
3. Using the tension straps on the top portion of the Y-Strap, pull the rotated shoulder back until a proper and comfortable position for the user is achieved.

SAFETY CHECK

Once the installation is finished, carry out the following checks to ensure that the user is appropriately and safely positioned.

If possible, ask the user to lean forward and to the sides to test the fit.

Check that:

1. The buckles and adjustment straps are operating properly.
2. The user sits comfortably in the seating and the harness.
3. The harness is not too low or too high. The support must not rub or cut into the users neck.
4. The mounting hardware is not interfering with any other positioning accessories/devices on the wheelchair.

HYGIENE MEASURES WHEN BEING RE-USED:

Prior to the wheelchair being re-used, it must be carefully prepared. All surfaces which come into contact with the user must be treated with a disinfection spray.

To do this, you must use a disinfectant as authorised/ recommended in your country, for rapid alcohol-based disinfection for medical products and medical devices, which must be disinfected quickly.

Please be aware of the manufacturer's instructions for the disinfectant you are using.

In general, a complete disinfection cannot be guaranteed on seams. We therefore recommend that you dispose of seat and back slings to avoid micro-bacterial contamination with active agents according to your local infection protection law.

CONDITIONS FOR SAFE STORAGE

- Products should be stored in a dry place: at 20 – 75% Relative Humidity.
- Protect from direct sunlight and dust for long periods.
- Away from direct heat and at a temperature between 5°C (40°F) to 40°C (104°F)
- Position and stack the packaged products according to the "this side up" arrows on the packaging.
- If a package is high and narrow and therefore potentially unstable, ensure it is well secured on the pallet or in the racking, to prevent it from falling down.
- Do not stack anything on top of the packaged products when the "do not stack" mark on the packaging is shown.
- Product should never be stored outside and subject to the elements.

24 MONTH LIMITED WARRANTY

Each Product is carefully inspected and tested to provide peak performance. Every Product is guaranteed to be free from defects in materials and workmanship for a period of 24 months from the date of purchase, provided normal use. Should a defect in materials or workmanship occur within 24 months from the original date of purchase, Sunrise Medical will, at its option, repair or replace it without charge. This warranty does not apply to punctures, tears or burns, nor to any removable cover.

Removable covers are guaranteed to be free from defects in materials and workmanship for a period of six (6) months from the date of purchase, provided normal use. Should a defect in materials or workmanship occur within six (6) months from the original date of purchase, Sunrise Medical will, at its option, repair or replace it without charge.

Claims and repairs should be processed through the nearest authorised supplier. Except for express warranties made herein, all other warranties, including implied warranties of merchantability and warranties of fitness for a particular purpose are excluded.

Type:	Product Name/ SKU Number	MD	This symbol means Medical Device
LOT	Lot		Manufacturer's address
xxxx-xx-xx	Date of Manufacture	EC REP	European Authorised Representative
CE	CE Mark	UK RP	UK Responsible Person
UK CA	UKCA Mark	CH REP	Swiss Representative's address
	Consult instructions for use		Importer's address

ISO 7010-M002
Instruction manual/
booklet must be read!
(Blue Icon)

The management system of SUNRISE MEDICAL is certified to EN ISO 13485 and ISO 14001.

We at SUNRISE MEDICAL have been awarded the ISO-13485 certificate, which affirms the quality of our products at every stage, from R & D to production. This product complies with the standards set forth in EU and UK regulations.

Indications

The varieties of fitting variants, as well as the modular design, mean that it can be used by those who cannot walk or have limited mobility because of:

- Joint contractures
- Muscle Dystrophy
- Neurological disease pattern
- Spasticity
- Assymetric Postures
- Hypo- Hyper tonic
- Cerebral Palsy
- MS
- SCI

Contraindication

No contraindications are known or available.

NOTE: General user advice. Not following these instructions may result in physical injury, damage to the product or damage to the environment!

A notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

As the manufacturer, SUNRISE MEDICAL, declares that this product conforms to the Medical Device Regulation (2017/745).

As the manufacturer, SUNRISE MEDICAL, declares that the product conforms to the UK Medical Devices Regulation 2002 No. 618.

B4Me special adaptations

Sunrise Medical strongly recommends that in order to ensure that your B4Me product operates, and performs as intended by the manufacturer; all the user information supplied with your B4Me product is read and understood, before the product is first used. Sunrise Medical also recommends that the user information is not discarded after reading it, but it is kept safely stored for future reference.

Medical Device Combinations

It may be possible to combine this Medical device with one or more other Medical Device or other product. Information on which combinations are possible can be found at www.Sunrisemedical.co.uk. All combinations listed have been validated to meet the General Safety and Performance Requirements, section 14.1 of the Medical Device Regulation 2017/745.

Guidance on the combination, such as mounting, can be found at www.SunriseMedical.co.uk.

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Sunrise Medical B.V.
Sunrise Medical HCM B.V.
Sunrise Medical S.A.S