JAY[®]

JAY Pelvic Positioning Belts

JAY Pelvic Positioning Belts

EN

Important Consumer Information.

NOTE: This manual contains important instructions that must be passed on to the user of this product. Please do not remove this manual before delivery to the end user.

SUPPLIER: This manual must be given to the user of this product.

USER: Before using this product, read this entire manual and save for future reference.

WARRANTY REGISTRATION: To validate the warranty on this product, please complete and return the postcard attached to the enclosed instruction and warranty booklet.

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

14

15

The intended lifetime of this product is 5 years.

SHIPPING:

Check all parts are free from shipping damage. In case of damage DO NOT use. Contact Sunrise Medical or your supplier for further instructions.

WARNING!

- JAY Positioning belts were designed for positioning, not to be used as transportation safety belt!
- JAY POSITIONING BELTS may be used while in transit but are not crash tested or approved to provide protection during a crash.
- An appropriate clinical assessment by a clinician or qualified wheelchair professional is highly recommended to guarantee safety and correct clinical usage.
- Any belt must not interfere with breathing.
- Never place a pelvic belt, or positioning support across the waist or abdomen.

The products described in this manual were designed to be fitted and installed exclusively by:

- An authorized Sunrise Medical supplier.
- A qualified and competent healthcare professional who has undergone appropriate certified training by either Sunrise Medical or their approved supplier.

PELVIC BELT INSTALLATION:

Initial Considerations:

The angle that the pelvic belt is mounted to the wheelchair has a direct impact on the pull angle on the pelvis (Fig. 1).

A 60° attachment angle to the seat rail positions the pelvic belt anterior and inferior to the axis, which reduces the tendency of the user with a posterior pelvic tilt to slide under the belt.

A 45° attachment angle on the back post encourages a posterior pelvic tilt, which assists in positioning users with an anterior pelvic tilt.

Pelvic belts can be installed in different ways according to user's seating system, wheelchair type and hardware used (Fig. 2-3).

All JAY POSITIONING BELTS include belt anchors for attachment (Fig. 4).

SAFETY:

For both types of belt, when the installation is finished, carry out the following procedure to ensure that the user is appropriately and safely positioned:

1. Check if the user is sitting comfortably.
2. Check if there are any areas of irritation. This usually manifests as a reddening of the skin especially in soft tissue.
3. Check that the user does not slide down through the belt, if they do, adjust the belt as required. If the user still slides under the belt and the belt is a two point model, assess if there is a need for the four point belt, this will provide additional security and adjustment.
4. Check if the pelvic belt mounting hardware is interfering with any other positioning accessories/ devices on the wheelchair and reposition it if necessary.

HYGIENE MEASURES WHEN BEING RE-USED:

Prior to the wheelchair being re-used, it must be carefully prepared. All surfaces which come into contact with the user must be treated with a disinfection spray.

To do this, you must use a disinfectant as authorised/ recommended in your country, for rapid alcohol-based disinfection for medical products and medical devices, which must be disinfected quickly.

Please be aware of the manufacturer's instructions for the disinfectant you are using.

In general, a complete disinfection cannot be guaranteed on seams. We therefore recommend that you dispose of seat and back slings to avoid micro-bacterial contamination with active agents according to your local infection protection law.

CONDITIONS FOR SAFE STORAGE

- Products should be stored in a dry place: at 20 – 75% Relative Humidity.
- Protect from direct sunlight and dust for long periods.
- Away from direct heat and at a temperature between 5°C (40°F) to 40°C (104°F)
- Position and stack the packaged products according to the "this side up" arrows on the packaging.
- If a package is high and narrow and therefore potentially unstable, ensure it is well secured on the pallet or in the racking, to prevent it from falling down.
- Do not stack anything on top of the packaged products when the "do not stack" mark on the packaging is shown.
- Product should never be stored outside and subject to the elements.

24 MONTH LIMITED WARRANTY

Each Product is carefully inspected and tested to provide peak performance. Every Product is guaranteed to be free from defects in materials and workmanship for a period of 24 months from the date of purchase, provided normal use. Should a defect in materials or workmanship occur within 24 months from the original date of purchase, Sunrise Medical will, at its option, repair or replace it without charge. This warranty does not apply to punctures, tears or burns, nor to any removable cover.

Removable covers are guaranteed to be free from defects in materials and workmanship for a period of six (6) months from the date of purchase, provided normal use. Should a defect in materials or workmanship occur within six (6) months from the original date of purchase, Sunrise Medical will, at its option, repair or replace it without charge.

Claims and repairs should be processed through the nearest authorised supplier. Except for express warranties made herein, all other warranties, including implied warranties of merchantability and warranties of fitness for a particular purpose are excluded.

Type:	Product Name/ SKU Number	MD	This symbol means Medical Device
LOT	Lot		Manufacturer's address
xxxx.xx.xx	Date of Manufacture	EC REP	European Authorised Representative
CE	CE Mark	UK RP	UK Responsible Person
UK CA	UKCA Mark	CH REP	Swiss Representative's address
	Consult instructions for use		Importer's address

ISO 7010-M002
Instruction manual/
booklet must be read!
(Blue Icon)

The management system of SUNRISE MEDICAL is certified to EN ISO 13485 and ISO 14001.

We at SUNRISE MEDICAL have been awarded the ISO-13485 certificate, which affirms the quality of our products at every stage, from R & D to production. This product complies with the standards set forth in EU and UK regulations.

Indications

The varieties of fitting variants, as well as the modular design, mean that it can be used by those who cannot walk or have limited mobility because of:

- Joint contractures
- Muscle Dystrophy
- Neurological disease pattern
- Spasticity
- Assymetric Postures
- Hypo- Hyper tonic
- Cerebral Palsy
- MS
- SCI

Contraindication

No contraindications are known or available.

NOTE: General user advice. Not following these instructions may result in physical injury, damage to the product or damage to the environment!

A notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

As the manufacturer, SUNRISE MEDICAL, declares that this product conforms to the Medical Device Regulation (2017/745).

As the manufacturer, SUNRISE MEDICAL, declares that the product conforms to the UK Medical Devices Regulation 2002 No. 618.

B4Me special adaptations

Sunrise Medical strongly recommends that in order to ensure that your B4Me product operates, and performs as intended by the manufacturer; all the user information supplied with your B4Me product is read and understood, before the product is first used. Sunrise Medical also recommends that the user information is not discarded after reading it, but it is kept safely stored for future reference.

Medical Device Combinations

It may be possible to combine this Medical device with one or more other Medical Device or other product. Information on which combinations are possible can be found at www.Sunrisemedical.co.uk. All combinations listed have been validated to meet the General Safety and Performance Requirements, section 14.1 of the Medical Device Regulation 2017/745.

Guidance on the combination, such as mounting, can be found at www.SunriseMedical.co.uk.

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Sunrise Medical B.V.
Sunrise Medical HCM B.V.
Sunrise Medical S.A.S