 **Wheelchair**
 **Rollstuhl**
 **Fauteuil roulant**

Nitrum Nitrum Pro

Adjusting castor fork angle, directional stability and castor fork height

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66009215

The castor fork adjustment mechanism has three settings:



- Castor fork angle
- Directional stability
- Castor fork height

The settings for directional stability and the castor fork height are combined; by twisting the slightly angled castor axis (A) (Fig. 1), the fork moves inwards (B) or outwards (C) (Fig. 2). At the same time, the distance of the fork from the frame also changes by 1 mm for every full rotation of the castor axis.

The detailed instructions for making the adjustments are described in the following section.

NOTE: The following adjustments may only be made by Sunrise Medical dealers

NOTE: You can also find a service video about this on the Sunrise Medical YouTube channel

WARNING:

- Faulty adjustments, or adjustments that have not been carried out correctly, may cause damage to the wheelchair or result in a fall and severe injuries to the user and other people.
- Wherever torque settings are specified, it is strongly recommended that a torque meter (not included), is used to verify correct torque specification is achieved.
- Some of the screws that are used during manufacture are fitted with threadlock, (blue dot on the thread) and can be undone and re-tightened up to three times before they must be replaced with new screws and threadlock. Alternatively, you can put Loctite™ 243 threadlock on the screws and re-fit them.

Replacement of the barrel nut

- In case there is a need to replace the barrel nut shown in Fig. 1-D, it is strongly recommended to add carbon assembly paste (part number 71208300), to ensure an optimal adjustment. It is highly recommended this adjustment to be performed by an authorized specialist.

Replacement of the barrel nut, please consider to remove the safety screw as shown in Chapter 5. Correct orientation after replacing the barrel nut, the countersink of the barrel nut must be rotated downwards.

1. Setting the directional stability:

- Insert special Sunrise Medical tool (C spanner, item no. 66009203-ET) into a ¼" ratchet or torque wrench
- Undo the round nut between the fork and castor connection with the C spanner (Fig. 3)
- Rotate the castor inwards, place a 90° setting gauge on the front of the fork and grasp the castor axis between the round nut and the fork with a flat 19 mm open-end spanner (item no. 66009203-ET) and rotate it until the front of the fork is level with the setting gauge (Fig. 4 and Fig. 5).

To do this, rotate it by a maximum of 90° in one direction, to avoid accidentally changing the height.

- After the adjustment has been made, keep the castor axis in position with the open-end spanner while tightening the round nut to a torque of 28 Nm.

Note: The fork angle must be set after adjustment has been made (see 3. Setting the castor fork angle)


Fig. 3


Fig. 4

Fig. 5

2. Adjusting the height:

If the wheelchair does not stand properly on 4 wheels, i.e. one of the castors is not completely touching the floor, the height of this castor should be adjusted as follows:

- Place the wheelchair on a flat surface, such as a worktop or similar.
- Load the footrest with a weight of approx. 2.5 kg, make sure that the air pressure is identical in both drive wheels, and rotate both castors into the straight line position (castors pointing to the back)
- Insert special Sunrise Medical tool (C spanner, item no. 66009203-ET) into a ¼" ratchet or torque wrench
- Undo the round nut between the fork and castor connection with the C spanner
- Grasp the castor axis between the round nut and fork with a flat 19 mm open-end spanner (item no. 66009203-ET) and rotate it anticlockwise. The castor will be lowered by 1 mm with each complete turn (360°)
- If the castor axis is rotated by one complete turn (360°), this will not influence the directional stability. In all other cases, the directional stability must be re-adjusted. See "1. Setting the directional stability" for more on this.
- Adjust the castor axis until all 4 wheelchair wheels are standing firmly on the ground
- After the adjustment has been made, keep the castor axis in position with the open-end spanner while tightening the round nut to a torque of 28 Nm.

3. Adjusting the castor fork angle

This adjustment is necessary if the seat height/centre of gravity has been changed, or the directional stability has been adjusted.

WARNING: an incorrectly adjusted castor fork angle may cause the castors to wobble, and lead to the user falling or severe injuries.

NOTE: setting the fork to toe-out (fork is inclined towards the back in the direction of travel) will increase the probability of the castors wobbling. Setting it to toe-in (fork is inclined towards the front in the direction of travel) can reduce the probability of the castors wobbling. The angle of the castor fork is set by the factory with a slight toe-in.

Loosen the Allen screws (A) (Fig. 6). The castor fork angle can now be adjusted to the front (toe-in) or to the back (toe-out).

To ensure that the castor fork is in an exactly horizontal position (90° to the floor, neither toe-in or toe-out), you will need a setting gauge (e.g. 90° flat angle, set square). Rotate the castors in the travelling position and place the setting gauge on the flat front side of the fork (Fig. 7). The castor fork will be at an angle of 90° to the floor when the setting gauge covers the entire surface.

If you want to set the fork slightly to toe-in, e.g. to prevent the castor fork from wobbling, set the angle so that a slight light gap can be seen between the setting gauge and contact area.

When you have set the correct position, tighten the allen screws (A) using the following procedure:

Tighten the outer allen screw first to a torque of 7Nm, than tighten the inner allen screw to a torque of 7Nm.

Repeat this procedure up to 6 times until both allen screws are evenly tightened to a torque of 7Nm

4. Replacing the castor fork

- Remove the locking clip (A) at the lower end of the castor axis (Fig. 8)
- Remove the nut (B) with a 13 mm open-end or ring spanner
- Now replace the fork
- Reassemble by following the instructions in reverse order. Observe the nut's tightening torque of 10 Nm (B)


Fig. 8


Fig. 6

Fig. 7

5. Replacing a Frog Leg Fork

To remove the Frog Leg Fork, please follow the below instructions:

- Remove the M6 screws shown in position 1.4 (Fig. 9) with an Allen key.
- Remove the cap shown in position 1.6 (Fig. 9)
- Remove the safety screws in position 1.5 (Fig. 9) with an Allen key.
- Open lock nut as shown in position 1.7 (Fig. 9) with tool kit 66009203-ET.
- Remove caster bolt in position 1.3 (Fig. 9) as shown in Position 1.1.

Once you are ready to assemble the new frog leg fork, please follow the below instructions always using the recommended torques in each step:

- Assemble the caster bolt in Pos 1.3 (Fig. 9)
- Lock the safety screws in Pos 1.5 (Fig. 9)
- Lock the cap shown in position 1.6 (Fig. 9)
- Adjust for caster fork angle, camber and height as described in the service manual.
- Lock the nut in Pos 1.7 (Fig. 9) using a torque of 28 Nm.
- Adjust the M6 screws shown in position 1.4 (Fig.9) using a torque of 7 Nm.


Fig. 9

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Sunrise Medical HCM B.V.
Sunrise Medical SAS