

# Q700-UPM

## HOW TO SET BLUETOOTH CONTROLS



# BLUETOOTH CONTROLS INTRODUCTION



Bluetooth controls are an option on the Sedeo Ergo order form that allow the client to pair their joystick of input device to a Smart phone, PC or Tablet. Once paired the client can control that device from their chair.

**This is not a standard option – It must be selected on the Sedeo Order Form**

You can select two Bluetooth Modules to control Android devices and one to control an apple device.

R-NET OPTIONAL CONTROLS / ELECTRONICS			£
JMV110031	● R-NET ISM-L (lighting module)	(required if L&I is selected)	
JMV110033	● R-NET CxSM (8-channel seating module)	(includes memory seat programming)	
JMV090068	□ R-NET BLUETOOTH MODULE (Controls 1 x Android and Windows device)		195
JMV090101	□ 2nd R-NET BLUETOOTH MODULE (Controls second Android and Windows device)		195
JMV090102	□ R-NET BLUETOOTH MODULE IOS (Controls Apple Device)		195
JMV090103	□ R-NET IO MODULE	(7 Switched controllable output commands, only available via B4Me)	252

If you wish to control more devices you can order a **CJSM2** with integrated Bluetooth via Built-4-Me. This can control 2 Android devices and two Apple devices as standard.

This guide will show you how to pair and customise Bluetooth for a demo and handover.



# BLUETOOTH CONTROLS

## PAIRING YOUR DEVICE



**Step 1:** Enter Bluetooth Menu on the joystick and push forward for 15 seconds until you hear a beep



**Step 2:** Pull back for 15 seconds until you hear a beep



**Step 3:** You will now be able to find R-Net Controls in your device Bluetooth menu




# BLUETOOTH CONTROLS

## PAIRING YOUR DEVICE



**Step 4:** Pair the devices together by selecting Rnet form the list of devices in the Bluetooth menu.



**Step 5:** To check you have connected successfully turn power on and off and select Bluetooth from the joystick Menu

A curser should appear On the device screen if you are connected.

**To use Bluetooth effectively you need to customise the commands via the PC programmer**





# BLUETOOTH CONTROLS

## PAIRING YOUR DEVICE



**Step 1:** Ensure your device is in Bluetooth discoverable mode.

**Step 2:** Enter Bluetooth via the joystick menu and push forward for 15 seconds until you hear a beep, then pull back for 15 seconds until you hear a beep.

**Step 3:** You will now be able to find R-Net Controls in your device Bluetooth menu





# BLUETOOTH CONTROLS

## CUSTOMISING BLUETOOTH



Once you have paired your device you can control the mouse automatically for that device but now you need to customise how the client will select ICONS and make commands like typing and sending messages.

The next steps will allow you to customise commands like **home, select, open, close** and enable the user to effectively use their devices from their joystick.

There are two easy ways in which a client can do this; via

- **Nudge functions**
- **Buddy button**



# BLUETOOTH CONTROLS

## CUSTOMISING BLUETOOTH – NUDGE



“**Nudges**” are the most common way of making commands with a joystick. A nudge is a fast push of the joystick in a direction.

The user can push the joystick to an icon on the device and “nudge” to open it up.

Nudge commands need to be programmed.

Nudge speed can also be programmed to suit a clients ability

**Step 1:** Connect the R-Net Dongle, Open the PG programmer and read out the program.






# BLUETOOTH CONTROLS


## CUSTOMISING BLUETOOTH – NUDGE



**Step 2:**Select **Mouse 1** from the menu  
on the left hand side



**Step 3:** Open the **Mouse 1** menu  
and the settings appear









# BLUETOOTH CONTROLS

## CUSTOMISING BLUETOOTH – NUDGE



**Step 4:** double click on the Forward Nudge – show to left

**Step 5 –** It will reveal the list of commands you can program. Select the desired command.





# BLUETOOTH CONTROLS

## CUSTOMISING BLUETOOTH – NUDGE



### POPULAR TEMPLATE

Below is a popular template that will work for most people. Start with this and then customise it to your client.

**Android Phone / Tablet**

**Reverse Nudge = Android OS Home (Phone Home Screen)**

**Left Nudge = Right Click (Back Up)**

**Right Nudge = Left Click (Forward)**

Mouse 1	
Forward Nudge	None
Reverse Nudge	Android OS Home
Left Nudge	Right Click
Right Nudge	Left Click
Nudge Time	300 ms
Pointer Speed	5 %
Pointer Acceleration	2000 ms
Action Beeps	Yes
Deflection Beeps	No
Double Click Time	0.3 s
External Switch Input 1	None
External Switch Input 2	None
Device Name	Gebs HTC
Screen Graphic	Computer



# BLUETOOTH CONTROLS

## CUSTOMISING BLUETOOTH – NUDGE



### CUSTOMISING NUDGE SPEED

A client must be able to double tap the joystick to complete a nudge command. The speed of this can be programmed to suit the ability of a client.

**Default = 200ms**

If a client is not quick enough, it can be increased

Mouse 1		
Forward Nudge		None
Reverse Nudge		None
Left Nudge		None
Right Nudge		None
Nudge Time		300 ms
Pointer Speed		10 %

If a client accidentally keeps nudging (they are too quick with controls) this can be reduced

Mouse 1		
Forward Nudge		None
Reverse Nudge		None
Left Nudge		None
Right Nudge		None
Nudge Time		100 ms
Pointer Speed		10 %



# BLUETOOTH CONTROLS


## CUSTOMISING BLUETOOTH – BUDDY BUTTON



Some clients can find it difficult to make the nudge commands and using a buddy button mounted to the armrest can be easier.

A buddy button can be programmed to “select” or “back up”.

Clients can use the joystick to navigate the mouse over icons and use the button to select (short press) or back up (long press)





# BLUETOOTH CONTROLS

## CUSTOMISING BLUETOOTH – BUDDY BUTTON



**Step 1:** Plug a buddy button into the external jack socket port via the joystick



**Step 2:** Connect the R-Net Dongle. Open the PG programmer and read out the program.






# BLUETOOTH CONTROLS

## CUSTOMISING BLUETOOTH – BUDDY BUTTON



**Step 3:** Select Assignable buttons form the Rnet menu shown below



**Step 4:** double click on first function / external jack socket 1 as shown below


Assign Buttons		First Function	Second Function
	Second Function Time	1.0 s	
	Profile Button	Profile Up	Profile Up
	Mode Button	Mode Up	Mode Up
	External Profile Jack 1	Mode Up	Mode Up
	External Profile Jack 2	Mode Up	Mode Up
	Speed Down Button	Speed Down	Speed Down

# BLUETOOTH CONTROLS

## CUSTOMISING BLUETOOTH – BUDDY BUTTON



**Step 5:** In the assignable button menu select “Left Mouse Button” or another function



**Step 6:** Repeat for second function and set to “Android OS Back”

Assign Buttons		First Function	Second Function
	Second Function Time	1.0 s	
	Profile Button	Profile Up	Profile Up
	Mode Button	Mode Up	Mode Up
	External Profile Jack 1	Mode Up	Mode Up
	External Profile Jack 2	Mode Up	Mode Up
	Speed Down Button	Speed Down	Speed Down

- ☐ Double Rgt Click Mouse
- ☐ Android OS Home
- ☒ Android OS Back
- ☐ Android OS Volume Up



# BLUETOOTH CONTROLS

## CUSTOMISING BLUETOOTH – BUDDY BUTTON



The program should look like this.

Program the chair and whilst hovering over the ICON use a short press to select / type and a long press to back up.

Remember time can be programmed via “second function time”

**Short Press = Left Mouse button = Enter**

**Long press = Right Mouse button = Back Up**

Assign Buttons		First Function	Second Function
	Second Function Time	1.0 s	
	Profile Button	Profile Up	Profile Up
	Mode Button	Mode Up	Mode Up
	External Profile Jack 1	Left Mouse Button	Right Mouse Button
	External Profile Jack 2	Mode Up	Mode Up
	Speed Down Button	Speed Down	Speed Down
	Speed Up Button	Speed Up	Speed Up
	Hom Button	Hom	Hom
	Lights Button	Lights	Lights
	Left Indicator Button	Left Indicator	Left Indicator
	Right Indicator Button	Right Indicator	Right Indicator
	Hazards Button	Hazards	Hazards
	Softkey 1	Hazards	Settings
	Softkey 2	Left Indicator	Left Indicator
	Softkey 3	Lights	Lights
	Softkey 4	Right Indicator	Right Indicator
	AB Latched Seating Timeout	45 s	



# Q700-UPM

## HOW TO SET I-DEVICE CONTROLS



# BLUETOOTH CONTROLS

## PAIRING YOUR DEVICE



**Step 1:** Enter I-Device Menu on the joystick and push forward for 15 seconds until you hear a beep



**Step 2:** Pull back for 15 seconds until you hear a beep




**Step 3:** You will now be able to find R-Net Controls in your device Bluetooth menu



# BLUETOOTH CONTROLS

## PAIRING YOUR DEVICE

**Step 4:** From the Apple device select: **Settings-General-Accessibility – Switch Control - Switches**




Within Accessibility there are Auto Scanning and Auto Hide options. It is recommended that these options are turned off for the initial set up.

# BLUETOOTH CONTROLS

## PAIRING YOUR DEVICE

### Step 5: Select Add New Switch - External

- To add a switch hold the joystick forward, right, back or left for a short medium or long duration then bring it back to the neutral position.
- When you have done this for one movement e.g. forward short.
- The following will appear for you to name the command e.g. forward short
- Once named select save



# BLUETOOTH CONTROLS

## PAIRING YOUR DEVICE



**Step 6:** Once saved the screen to below will appear, select which command you wish to assign to the nudge.

Choose a switch action.

SCANNER

Select Item

Scanner Menu

Resume Auto Scanning

Move To Next Item

Move To Previous Item

Stop Scanning

SYSTEM

Tap

App Switcher


Home Button

Notification Centre

Siri

Decrease Volume

**Step 7:** Once you have selected the command, repeat step 5 and step 6 to set up more nudge commands. You will be able to set up to 12 nudge commands.



# BLUETOOTH CONTROLS


## PAIRING YOUR DEVICE



**Step 8:** Once you have assigned the commands to the nudges as required, come out of **switches** to the **switch control menu** (see right)

Finally **turn on** switch control.

The I-Device is now connected to the Apple device.





# BLUETOOTH CONTROLS CUSTOMISING I DEVICE– NUDGE



## CUSTOMISING NUDGE SPEED

A client must be able to tap the joystick to complete a nudge command. The speed of these nudges can be programmed to suit the ability of a client.

**Short Nudge = 50ms – 200ms**

**Medium Nudge = 200ms – 300ms**

**Long Nudge = 300ms - +**

If a client is not quick enough, it can be increased

Mouse 1		
Forward Nudge		None
Reverse Nudge		None
Left Nudge		None
Right Nudge		None
Nudge Time		300 ms
Pointer Speed		10 %

If a client accidentally keeps nudging (they are too quick with controls) this can be reduced

Mouse 1		
Forward Nudge		None
Reverse Nudge		None
Left Nudge		None
Right Nudge		None
Nudge Time		100 ms
Pointer Speed		10 %