









# Tango®

# Setting a new standard in powerchairs

The Quickie Tango is the hassle free powerchair with the most complete standard configuration and a wide range of manual and powered options.



# Great standard features

The Quickie Tango offers the most complete features in its standard configuration: Full adjustability, tilt and recline, comfort ride and secure transportation.





Mechanical seat tilt
The seat can be tilted and set
manualy within a range of 0° to 9°
(with tools).



Mechanical reclining back
The backrest recline can be adjusted
from -3° to 12° in 3° increments
(without tools)



Tension & height adjustable back Backrest uphostery is tension adjustable for improved user comfort, and can be adjusted to 45, 47,5 and 50 cm height



Adjustable seat width

Through the width adjustable armrest receiver the seat width can be adjusted from 40 to 50 cm

## Adjustable seat depth

The backrest can be moved backwards to offer seat dept adjustment from 42.5 to 50 cm in 2,5 cm



Fold down backrest Fold down backrest for easy transportation and storage

Height adjustable armrest: The single post armrests are also height adjustable as standard. Width adjustable hangers: Hangers offer 4cm of outward width adjustment.



### Drive wheel suspension

Drive wheel suspension as standard impoves comfort and minimises the effect of challenging environments

### Tie down points

Integrated and easy accessible tie down points at the rear and at the front of the frame for secure transportation in vehicles



# Comfort seating

Pro Seating comfort cushion and back is available as an option to improve user comfort.



Standard Headrest For standard upholstery back



Comfort Headrest For comfort back



# **Electronics**

Simple yet efficient electronics make the chair easy to handle and to service and can accommodate all the available powered options.

## VR2 programmable electronics

- VR2 controller controls up to two actuators via the joy stick
- 50 Amp controller as standard
- 60 Amp controller with two actuator outputs (optional)
- Dual control option
- Four actuator control box, optional



Four actuator control box (option)



Parallel Swing- Away contoller (standard)



Easy to service

Easy accesible batteries for service and mantenance

# Additional options:

- Three colour choices: Brilliant Silver, Brilliant Blue and Brilliant Red
- Two choices of wheel size (12" x 8" and 14" x 9") pneumatic or solid
- Crutch holder
- Kerb climber

# Wide range of Manual and Powered Options

In addition to the standard features the Quickie Tango offers a wide range of manual and powered options. Seat tilt, Back recline and Elevating legrest, to accomodate individual end user needs.



# Elevating legrest

- Manual elevating legrest
- Powered elevating legrests can be adjusted via the joystick

### Back recline

- Manual backrest recline angle with gas strut adjustable from -3° to -30°
- Powered backrest recline adjustable through the joystick from -3° to 30°

### Seat tilt

- Manual seat tilt with actuator lever, adjustable from 3° to 25°
- Powered seat can tilt from 3° to 25° by operating it through the control system



Lights & Indicators
Modern looking lights and
indicators available as
option for outdoor safety



Drive wheel mud guards
Optional mud guards improve the chair look and and protect the user

# Tango®



# Seating system adjutsments

Seat depth:	Adj. : 42.5, 45, 47.5, 50 cm	Backrest heights:	45 / 47.5 / 50 cm	
Seat width:	Adj. from 40 to 50 cm	Backrest angles:	Adj3° - 12°	
Seat height:	43 cm (with 3° tilt)	Armrest height:	Adj. from 27 to 35.5 cm	
Seat angles:	Adj. 3°, 6°, 9°			



# Technical data (Measurements taken with the standard configuration)

Max. speed:	6 kph	Total weight (with 50 Amph batteries)	100 kg
Max. range: 50 Amph batt.	25 km (ISO 7176-4)	Max user weight:	125 kg
Batteries:	40 / 50 Amph	Controllers:	VR2 std
According EN 12184 classified as Class B		Crash tested (ISO 7176-19)	
Minimum turning radious:	90 cm		
			CR40

For further information on the full specification and options and accessories please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text



Kahlbachring 2-4 69254 Malsch/HD

Tel.: +49 (0) 72 53 / 9 80-0 Fax: +49 (0) 72 53 / 9 80-222

E-Mail:international@sunrisemedical.de

www.SunriseMedical.eu