



SIMPLY · EFFECTIVE · SEATING







Contoured foam base



Applewood solid seat (option)

## High level comfort and stability

## **Clinical Application:**

- Client with mild symmetrical or mild asymmetrical posture and moderate postural support requirements
- Moderate risk of skin breakdown
- Client requires comfortable, no maintenance cushion

## Features:

- Exceptionally light and comfortable
- Sacral notch and Back cane cuts
- Incontinence cover
- Choice of flat or bevelled bases for solid or sling seats
- Optional solid seats
- Visco foam filled well

	Clinical Performance		
	Factors	Product Assessment	
	Skin Protection	2 cm ischial immersion into memory foam	_
	Shear Prevention	Ischial immersion and increased contouring to limit sliding	_
	Stability Provided	Thigh and trocanter support provides lateral and posterior pelvic stability	_
	Positioning	Stable foam base for moderate positioning	_
	Comfort	Exceptional softness, firmness and shape enabling prolonged periods seated	
	Clinical Benefit	High comfort and stability	

Order Part No: JEV Width, Depth, and Flat or Curved ie. JEV2525

Specifications				
Positioning	High level comfort and stability			
Purpose	Lightweight cushions with high level comfort and stability			
Skin Integrity (Max score is 5)	3			
Shear (Max score is 5)	2			
Stability (Max score is 5)	4			
Positioning (Max score is 5)	3			
Seating Tolerance (Max score is 5)	4			
User Weight	113kg			
Product Weight (MIN)	0.50 kg			
Product Weight (MAX)	2.70 kg			
Cushion Type	Open Cell Foam			
Cover Type	Incontinenece, Air Exchange / Microclimatic			
Width and Depth	25 - 60 cm			
Height Front and Rear	6 cm			
Guarantee and usage	2 years			



