JAY®

LITE[™]

The lightest cushion in its class, providing the perfect complement to an ultra lightweight wheelchair.



©Sunrise Medical. Jay Lite. 2013. Rev A.AU

JAY LITE

- The JAY Lite was designed through the careful consideration of the science behind seating. The result: An extremely lightweight cushion that features superior pressure reduction with microclimatic heat and moisture dissipation. The key features of the JAY Lite that provide these great benefits include 3DX[™] microclimatic spacer fabric, airflow base, and OptiWell[™] pelvic loading area. It's the perfect cushion for users on the go.

Extremely Lightweight:

- Product Weight: 771 g (16" x 18") Skin Protection cushion 790 g (16" x 18") - Skin Protection and Positioning cushion
- Perfect complement to ultra lightweight aluminum and titanium wheelchairs

Superior Pressure Distribution:

- Based on anthropometric data, the OptiWell[™] pelvic loading area is optimally shaped and sized to evenly distribute load under the pelvis
- Ischial cutouts and scooped reticulated foam insert minimize compression and reaction forces when the cushion is loaded

Microclimatic Heat and Moisture Dissipation:

- The 3DX[™] Microclimatic Spacer Fabric and airflow base provide ventilation and airflow beneath the user
- Protects the skin against heat and moisture build-up

Quickie Bundling:

• The JAY Lite can be ordered and shipped complete with a variety of Quickie manual and power chairs

• This helps to save time and money!





Sunrise Medical Australia Pty Ltd. 6 Healey Circuit, Huntingwood, NSW 2148

Email: enquiries@sunrisemedical.com.au

Phone: (02) 9678 6600

www.sunrisemedical.com.au

Fax: (02) 9678 6655

SPECIFICATIONS	
Jay Lite	
Weight	771 g (16" x 18")
Product Width	35.5 cm - 61 cm (14" – 24")
Product Length	35.5 cm - 61 cm (14" – 24")
Product Height	8.6 cm (3.4") - Skin Protection
	11.5 cm (4.5") - Skin Protection and Positioning
Weight Capacity	113.5 kg
Cover Options	Microclimatic 3DX [™] Spacer Fabric
Warranty	24 months
TGA #	145061

