WHITMYER®

= Head Support Solutions



Whitmyer is a product line of head supports and accessories that addresses the range of user needs from basic support to multiple-point control. A variety of pad shapes and sizes are available in single, two-pad or three-pad systems with each pad designed to address a specific support purpose.

TABLE OF CONTENTS

HEAD SUPPORTS







PLUSH

SPECIALTY PLUSH

CONTOUR CRADLE







S.O.F.T. DUAL SUB-OCCIPI-TAL Adult & Paediatric



HEADS UP



CUDDLES INFANT

HARDWARE SYSTEMS



AXYS



ONYX



LINX





COBRA XTRA Flip-Back

HEAD SUPPORTS

PLUSH Single Pad Headrest	4
SPECIALTY PLUSH Single Pad Headrest	
CONTOURED CRADLE Single Pad Headrest	
S.O.F.T. SINGLE SUB-OCCIPITAL	10
S.O.F.T. DUAL SUB-OCCIPITAL Adult & Paediatric	12
HEADS UP Anterior Head Support System	14
CUDDLES INFANT Contoured Cradle or Triple Pad Headrest	16

ACCESSORIES

LATERAL FACIAL SUPPORTS	17
ANTERIOR SUPPORT OPTIONS	17

HARDWARE SYSTEMS


AXYS	18
NYX	18
.INX	18
PRO	18
COBRA XTRA & COBRA XTRA Flip Back	18
CUSTOM MOUNTING OPTIONS	19

PLUSH Single Pad Headrest Systems utilize dual density foam layers on a contoured, curved panel. These single pad headrests are easily adjustable and accommodate changes or inconsistencies in body positioning. Unique "cupping" shape provides posterior support with a high degree of comfort and safety. This basic headrest is ideal for tilt/recline, and has the ability to add lateral support, anterior support, or switches as client's needs change. The larger size pads facilitate safety for clients with behavioral issues and active movements.



PLUSH ORDERING GUIDE

DETERMINE PAD SIZE



SELECT HARDWARE SYSTEM AND PAD COVER (BASED ON PAD SIZE)

COBRA XTRA AXYS ONYX LINX PRO **COBRA XTRA** Flip Back Pad Size (Cover selection: L=Lycra V=Dartex Reversed) Part Numbers ONYX-P4 L LINX-P4 L PRO-P4 L COBRA-P4 L FB-COBRA-P4 L AXYS-P4 L AXYS-P4 V ONYX-P4 V LINX-P4 V PRO-P4 V COBRA-P4 V FB-COBRA-P4 V ONYX-P2 L LINX-P2 L PRO-P2 L COBRA-P2 L AXYS-P2 L FB-COBRA-P2 L AXYS-P2 V ONYX-P2 V LINX-P2 V PRO-P2 V COBRA-P2 V FB-COBRA-P2 V ONYX-P1 L LINX-P1 L PRO-P1 L COBRA-P1 L FB-COBRA-P1 L AXYS-P1 L AXYS-P1 V ONYX-P1 V LINX-P1 V PRO-P1 V COBRA-P1 V FB-COBRA-P1 V AXYS-P3 L ONYX-P3 L LINX-P3 L PRO-P3 L COBRA-P3 L FB-COBRA-P3 L AXYS-P3 V ONYX-P3 V LINX-P3 V PRO-P3 V COBRA-P3 V FB-COBRA-P3 V AXYS-P8 L ONYX-P8 L LINX-P8 L PRO-P8 L COBRA-P8 L FB-COBRA-P8 L AXYS-P8 V ONYX-P8 V LINX-P8 V PRO-P8 V COBRA-P8 V FB-COBRA-P8 V

SELECT MOUNTING BRACKET FOR HARDWARE SYSTEM

3

			Moun	ting Brackets	
Detachable	0 0				C
ē	M2500	M2100	M2300	M2200	
Swing-Away					S
ау	M5500	M5100	M5300	M5200	

COBRA XTRA Hardware System includes detachable mounting bracket.

Select COBRA XTRA Flip Back Hardware System for swing-away capability.


Specialty PLUSH Single Pad Headrest Systems are formerly custom solutions that are now offered as part of the standard product line including Adjust-A-PLUSH (shown below), Dual Mount and 12" Narrow pads.

SPECIALTY PLUSH



ADJUST-A-PLUSH

Designed for use on clients that exhibit poor head control and require lateral support on at least one side to maintain an upright head posture. Assists in controlling passive and active lateral cervical flexion. Suitable for clients with poor targeting that require some degree of lateral confinement in their headrest. The Adjust-A-PLUSH has two tapered locking joints that can be narrowed or widened to provide maximum lateral support. Zippered cover hides adjustment points.



DUAL MOUNT PLUSH

An extra mount provides added strength making this head support ideal for high tone, high extension clients. The pad is offered in 19" and 14" widths.



NARROW PLUSH

Offers low profile support in a 12" width.





Check out the SPECIALTY PLUSH's versatility and adjustability in action. www.YouTube.com/user/JAYCushionsandBacks


Click on the playlist: "Whitmyer Headrest Adjustment Systems"

to view now!

Scan this QR code


SPECIALTY PLUSH ORDERING GUIDE

DETERMINE PAD SIZE



SELECT HARDWARE SYSTEM AND PAD COVER (BASED ON PAD SIZE)

2



SELECT MOUNTING BRACKET FOR HARDWARE SYSTEM

3

V I	CON					
				Moun	ting Brackets	
	Detachable	0 0				COBRA X
	ē	M2500	M2100	M2300	M2200	
	Swing-Away					Select Co
	ау	M5500	M5100	M5300	M5200	

OBRA XTRA Hardware System includes detachable mounting bracket.

Select COBRA XTRA Flip Back Hardware System for swing-away capability.

6


7

Contoured Cradle Head Support's positioning success is derived from its unique anatomic shape. Generically designed to cradle the occiput and provide lateral cervical support, this unit supplies a high degree of comfort and positioning in a simple-to-adjust head support. A single rotational collar adjustment allows positioning for maximum support surface contact. Four pad sizes fit infant to adult. Add lateral facial pads, switches, and/or anterior support for more versatility.

CONTOURED CRADLE



CONTOURED CRADLE ORDERING GUIDE



SELECT HARDWARE SYSTEM AND PAD COVER (BASED ON PAD SIZE)

COBRA XTRA **AXYS** ONYX LINX PRO **COBRA XTRA** Flip Back Pad Size Part Numbers (Cover selection: L=Lycra V=Dartex Reversed) AXYS-CCS4 L ONYX-CCS4 L LINX-CCS4 L PRO-CCS4 L COBRA-CCS4 L FB-COBRA-CCS4 L AXYS-CCS4 V ONYX-CCS4 V LINX-CCS4 V PRO-CCS4 V COBRA-CCS4 V FB-COBRA-CCS4 V ONYX-CCS3 L LINX-CCS3 L PRO-CCS3 L COBRA-CCS3 L AXYS-CCS3 L FB-COBRA-CCS3 L ONYX-CCS3 V AXYS-CCS3 V LINX-CCS3 V PRO-CCS3 V COBRA-CCS3 V FB-COBRA-CCS3 V AXYS-CCS2 L ONYX-CCS2 L LINX-CCS2 L PRO-CCS2 L COBRA-CCS2 L FB-COBRA-CCS2 L ONYX-CCS2 V COBRA-CCS2 V AXYS-CCS2 V LINX-CCS2 V PRO-CCS2 V FB-COBRA-CCS2 V AXYS-CCP4 L ONYX-CCP4 L LINX-CCP4 L PRO-CCP4 L COBRA-CCP4 L FB-COBRA-CCP4 L AXYS-CCP4 V ONYX-CCP4 V LINX-CCP4 V PRO-CCP4 V COBRA-CCP4 V FB-COBRA-CCP4 V AXYS-CCP3 L ONYX-CCP3 L LINX-CCP3 L COBRA-CCP3 L FB-COBRA-CCP3 L PRO-CCP3 L AXYS-CCP3 V ONYX-CCP3 V LINX-CCP3 V PRO-CCP3 V COBRA-CCP3 V FB-COBRA-CCP3 V AXYS-CCP2 L ONYX-CCP2 L LINX-CCP2 L PRO-CCP2 L COBRA-CCP2 L FB-COBRA-CCP2 L AXYS-CCP2 V ONYX-CCP2 V LINX-CCP2 V PRO-CCP2 V COBRA-CCP2 V FB-COBRA-CCP2 V

SELECT MOUNTING BRACKET FOR HARDWARE SYSTEM

3

			Moun	ting Brackets	
Detachable	0 0				СОВ
ē	M2500	M2100	M2300	M2200	
Swing-Away					Sele
ау	M5500	M5100	M5300	M5200	

OBRA XTRA Hardware System includes detachable mounting bracket.

Select COBRA XTRA Flip Back Hardware System for swing-away capability.

via a single curved sub-occipital pad. Unique to our single sub-occipital system is the independent

S.O.F.T. SINGLE SUB-OCCIPITAL



"Whitmyer Headrest Adjustment Systems"




S.O.F.T. SINGLE ORDERING GUIDE

DETERMINE PAD SIZE



SS10 10" PLUSH occipital w/Adult Single Sub-occipital



SS8 Adult occipital w/Adult Single Sub-occipital




Adult occipital w/Paediatric Paediatric occipital w/Paedi-Single Sub-occipital



atric Single Sub-occipital

SELECT HARDWARE SYSTEM AND PAD COVER (BASED ON PAD SIZE)



SELECT MOUNTING BRACKET FOR HARDWARE SYSTEM

			Moun	ting Brackets	
Detachable	0			M2200	
	M2500	M2100	M2300	M2200	
Swing-Away			Police		S
ау	M5500	M5100	M5300	M5200	

COBRA XTRA Hardware System includes detachable mounting bracket.

Select COBRA XTRA Flip Back Hardware System for swing-away capability.

S.O.F.T. Three Pad Head Supports solve the most challenging head positioning issues. This patented device combines two width and rotationally adjustable sub-occipital pads with an adjustable occipital pad for independent anterior/posterior movement. The added adjustability of the individual suboccipital pads allows adjustment to asymmetrical postures and deformities. This unique pad configuration supplies maximum occipital and proximal lateral cervical support for more aggressive positioning and comfort. Consistent, repeatable positioning of the client is required for best performance. Attaching lateral facial pads, switches, and/or anterior supports adds to this system's versatility.

S.O.F.T. DUAL SUB-OCCIPITAL

ADULT & PAEDIATRIC



Check out the S.O.F.T.'s versatility and adjustability in action.

www.YouTube.com/user/JAYCushionsandBacks Click on the playlist:

"Whitmyer Headrest Adjustment Systems"

Scan this QR code to view now!



S.O.F.T. DUAL ORDERING GUIDE

DETERMINE PAD SIZE



8– 8" PLUSH occipital w/ Size-1 Sub-occipitals



5– 10" PLUSH occipital w/ Size-1 Sub-occipitals



2— Adult occipital w/Size-2 Sub-occipitals



1- Adult occipital w/Size-1 Sub-occipitals



7– 8" PLUSH occipital w/ Size-3 Sub-occipitals



6- Paediatric occipital w/ Size-4 Sub-occipitals




4- Adult occipital w/Size-3 Sub-occipitals



3- Paediatric occipital w/ Size-3 Sub-occipitals

SELECT HARDWARE SYSTEM AND PAD COVER (BASED ON PAD SIZE)

2



SELECT MOUNTING BRACKET FOR HARDWARE SYSTEM

3

VI	MOUN ING BRACKET FOR HARDWARE SYSTEM							
				Moun	ting Brackets			
	Detachable	0				CC		
		M2500	M2100	M2300	M2200			
	Swing-Away			FBI		Se		
	⋖	M5500	M5100	M5300	M5200			

OBRA XTRA Hardware System includes detachable mounting bracket.

Select COBRA XTRA Flip Back Hardware System for swing-away capability.

12

HEADS UP

HEADS UP ORDERING GUIDE

The Heads Up features easy to create Smart Part numbers for configuring the headrest. Follow the 7 easy steps below to generate a part # that can be added into our online ordering system or ordered through customer service.

EXAMPLE SMART PART #: AXYS-HU-ITSSL Heads Up with AXYS Hardware - Infant Size, Tall Height, Short Arms, Standard Foam, Lycra (Neoprene) Cover

ANTERIOR AND LATERAL STABILIZATION OF SHOULDERS AND TRUNK FOR OPTIMAL **HEAD CONTROL**

- Innovative anterior stabilization arms
- Swing-out for easy transfer that maintains head support
- Lower lateral adjustments to accommodate varying body shapes


TWO DIFFERENT INDUSTRY-PROVEN

- Contoured Cradle
- Single Sub-Occipital

CHIN PAD INCLUDED


POSTERIOR SUPPORT OPTIONS:





www.YouTube.com/user/JAYCushionsandBacks Click on the playlist:

"Whitmyer Headrest Adjustment Systems"





SELECT HARDWARE SYSTEM Begin part # with hardware selection followed by -HU for Heads

AXYS	ONYX	LINX	PRO	COBRA XTRA	COBRA XTRA Flip Back
AXYS-HU-	ONYX-HU-	LINX-HU-	PRO-HU-	COBRA-HU-	FBCOB-HU-

SELECT NECK WIDTH (MAXIMUM) Add Selection to Smart Part #

Infant (2.8")	Small (3.5")	Medium (4.1")	Large (5.3")
1	S	M	L



SELECT PAD HEIGHT Add Selection to Smart Part #

Tall	Short
Т	S



shown with Occipital Pad (Step 7)

SELECT SWING-OUT SUPPORT ARM LENGTH Add Selection to Smart Part #

Short	Long
S	L



Infant & Small	Medium & Large
Short: 5.6"	Short: 7"
Long: 7.3"	Long: 9"

SELECT FOAM TYPE Add Selection to Smart Part #

Standard	Comfort
S	С

	Infant	Small	Medium	Large
Standard	0.4"	0.5"	0.6"	0.7"
Comfort	0.6"	0.6"	0.8"	0.9"

SELECT COVER Add Selection to Smart Part #

*Available on short height

Lycra (Neoprene)	Reversed Dartex		
L	V		

IF YOU DO NOT WANT TO ADD AN OCCIPITAL PAD, YOUR PART NUMBER IS COMPLETE. IF YOU DO, PLEASE PROCEED TO STEP 7.


SELECT OPTIONAL OCCIPITAL PADS* Add Selection to Smart Part #

	Paediatric 5" Pad (for Short Infant size)	Adult 6" Pad (for Short Small size)	Plush 8" Pad (for Short Medium size)	Plush 10" Pad (for Short Large size)
ts only.	-05	-06	-08	-10

SELECT MOUNTING BRACKET FOR HARDWARE SYSTEM

	Mounting Brackets for:						
	AXYS ONYX LINX PRO COBRA XTRA COBRA XTRA Flip Back						
Detach- able	M2500	M2100	M2300	M2200	COBRA XTRA Hardware System includes detachable mounting bracket.		
Swing- Away	M5500	M5100	M5300	M5200	Select COBRA XTRA Flip Back Hardware System for swing-away capability.		

CUDDLES INFANT



LATERAL FACIAL SUPPORTS



ORDERING GUIDE

DETERMINE PAD SIZE









CCS Infant Contoured Cradle Pad

SELECT HARDWARE SYSTEM AND PAD COVER (BASED ON PAD SIZE)




Cover selection: L=Lycra V=Dartex Reversed

SELECT MOUNTING BRACKET FOR HARDWARE SYS-



Check out the CUDDLES INFANT's versatility and adjustability in action.

www.YouTube.com/user/JAYCushionsandBacks Click on the playlist: "Whitmyer Headrest Adjustment Systems"





INFT500




INFT700






Part numbers for the pads and hardware are listed below. Hardware is swing-away and must be ordered with the pad part number.



Unilateral

Bilateral

NOTE: All lateral supports come with swing-away hardware



	SPOT	STANDARD	LONG	OVAL	LARGE-C	MED-C	SMALL-C	FSO-B
		°		(a)			8	20
				Part Nu	umbers			
	Plush & Spec.Plush							
	LFC10-PL	LFC15-PL	LFC20-PL	LFC25-PL	LFC30-PL	LFC35-PL	LFC40-PL	LFC80-PL
\Box	& LM115							
niLatera	Cont.Cradle/Heads	Cont.Cradle/Heads			Cont.Cradle/Heads	Cont.Cradle/Heads		Cont.Cradle/Heads
	Up	Up			Up	Up		Up
era	LFC10-CS	LFC15-CS	LFC20-CS	LFC25-CS	LFC30-CS	LFC35-CS	LFC40-CS	LFC80-CS
= =	& LM115							
	S.O.F.T. Single, Adult/							
_	Peds	Peds			Peds	Peds		Peds
	LFC10	LFC15	LFC20	LFC25	LFC30	LFC35	LFC40	LFC80
	& LM115	& LM115 Plush & Spec.Plush						
學	Plush & Spec.Plush	· ·			· ·			
atera	LFC10-PL2	LFC15-PL2	LFC20-PL2	LFC25-PL2	LFC30-PL2	LFC35-PL2	LFC40-PL2	LFC80-PL2
	& (2) LM115 Cont.Cradle/Heads							
	Up	Up			Up	Up		Up
	LFC10-CS2	LFC15-CS2	LFC20-CS2	LFC25-CS2	LFC30-CS2	LFC35-CS2	LFC40-CS2	LFC80-CS2
	& (2) LM115							


Anterior Supports are excellent for low tone individuals with poor head balance, individuals exhibiting forward flexion patterns, or individuals who require an upright head position for medical reasons, i.e. compromised respiration and swallowing. Successfully used for task oriented applications, such as classroom interaction or during meals. The addition of a Strap Cap brings a low profile solution to anterior support for individuals who may not be successful with the strap.

ANTERIOR SUPPORT OPTIONS

	DYNAMIC FOREHEAD STRAP	DFS ²	CTDAD CAD	STATIC FOREHEAD STRAP
	DINAMIC FORENEAD STRAP	DF3 ⁻	STRAP CAP	STATIC FUNEREAD STRAP
	Part Nui		mbers	
Adult	DFS10	DFSQ10	CAP10 (black)	HU10
Peds	DFS30	DFSQ30	CAP30-R (red) CAP30-B (blue)	HU10

MOUNTING HARDWARE OPTIONS

Each head support pad is available with most of the hardware options listed below. Part numbers for order can be found on the detailed head support product pages.



AXYS is a compact mounting system that features a black chrome finish combined with a ½" steel square rod. It was designed for rigorous applications, ease of adjustment, and simple modifications.

Height: 12" Horizontal Depth: 6" Pitch Adjustment: 30°

LINX




LINX is a contemporary styled mounting hardware that combines two simple. unique points of adjustment which allows full anterior/posterior fluid movement, creating an infinite range of support surface placement. It features satin black tubing with contrasting polished aluminum clamps and knurled mounting ball ensures stay-in-place adjustments.

Height adjustment: 7" 8.5" forward. Horizontal Depth: ¾" rearward

Pitch Adjustment: 30°

COBRA XTRA



The COBRA XTRA Mounting System incorporates a sleek and unobtrusive design that eliminates protrusions towards the rear of the chair. With superior strength and positioning range, the COBRA XTRA Mounting System is available in Small (Early Intervention and Paediatric), Medium (Paediatric), and Large (Adult). The components are modular and easily changed in the field to enable exact positioning of the head support pad – even if the user's head position changes from the time of the initial prescription.

Height adjustment: 6.5" Horizontal Depth Adj: 6.5" forward,

Pitch Adjustment: 30°

Paediatric

Height adjustment: 3.5" Horizontal Depth Adi: 5" forward.

1.5" rearward

2" rearward


Pitch Adjustment: 30°

Early Intervention

Height adjustment: 3.5" Horizontal Depth Adj: 4" forward,

1/2" rearward

Pitch Adjustment: 30° ONYX



PRO

8-1/4" Peds

10-1/4" Adult

ONYX is a compact mounting solution that features a black chrome finish combined with a small profile 7/8" steel rod. It was designed for rigorous applications, ease of adjustment, and simple modifications.

Paediatric

Height adjustment: 8" Horizontal Depth Adj: 1" to 6" Pitch Adjustment: 30°

Adult

Height adjustment: 11" Horizontal Depth Adj: 1" to 8" Pitch Adjustment: 30°

The PRO Mounting System is a proven hardware option for demanding applications. This mounting system features 1.25" aluminum extrusion and is ideal for use when client strength or strong tonal patterns are an issue.

Paediatric

Height adjustment: 6"

Horizontal Depth Adj: 1.5" forward, 6" rearward

Pitch Adjustment: 45°

Adult

Height adjustment: 6"

Horizontal Depth Adj: 1.5" forward, 8" rearward

Pitch Adjustment: 45°

COBRA XTRA Flip Back



The COBRA XTRA Flip Back Mounting System is operated by a single push-button on the side of the bracket and neatly flips the assembly back behind the wheelchair back. With superior strength and positioning range, the COBRA XTRA Flip Back Mounting System is available in Small (Early Intervention and Paediatric), Medium (Paediatric), and Large (Adult). The components are modular and easily changed in the field to enable exact positioning of the head support pad – even if the user's head position changes from the time of the initial prescription.

Adult

Height adjustment: 6.5" Horizontal Depth Adj:

6.5" forward, 2" rearward

Pitch Adjustment: 30°

Paediatric

Height adjustment: 3.5" Horizontal Depth Adj:


5" forward, 1.5" rear-

ward

30° Pitch Adjustment:

Early Intervention Height adjustment: 3.5"

Horizontal Depth Adj: 4" forward,



CUSTOM MOUNTING OPTIONS

FLEX INTERFACE BRACKET

The Flex Interface Bracket provides upward flex, travel, and return when the head support is extended against. This dynamic motion assists in maintaining head support contact throughout the user's exten-

pattern. The Flex Interface Bracket utilizes elastomeric bushings engineered to supply the correct resistance-to-force while allowing the head support to return to the adjusted position. Used between the head support mount and contact surface, the Flex Interface Bracket can be used with most Whitmyer head support accessories.


Part Number: M9450

1-7/8

UNIVERSAL SWING-AWAY MOUNT

Whitmyer detach brackets attach to our universal swing-away mount. When attached to the swing-away mount, the assembly becomes a combination detach/swing-away mounting bracket.


Part Number: M5000



SLIDE MOUNT

This slotted plastic plate allows the headrest/head support to be positioned in a more perpendicular fashion behind the head of individuals who present with a significant amount of asymmmetry. The plate bolts into the chair's existing hole pattern and allows 4" of movement left or right. Perfect for trying off-set head support locations, adjusting for misplaced T-nut locations and relocating the head support to be more in line with the direction of applied forces.


Part Number: M7010



PRO OFF-SET MOUNTING SYSTEM

PRO system off-set bracket provides 2-1/2" of rearward adjustment to the horizontal bar. Flipping the horizontal bar will also allow for 2-1/2" of forward adjustment.

Part Number: PRC7070



ONYX REARWARD OFF-SET VERTICAL ROD

ONYX off-set vertical rod provides 2-1/2" of rearward/forward adjustment to the horizontal rod.

Part Number: OX6030

19 18

Whitmyer Head Supports can be ordered with any Quickie*, Zippie* or Breezy* mobility product base and as a great addition to a JAY* standard or custom seating solution.

Order them together to ensure compatibility prior to shipment and save time and money with ONE ORDER and ONE INVOICE.

Sunrise MobilityLink



Adult Manual

























Paediatric







Kidz











Breezy Standard







BasiX 2

Rubix 2 XL

Breezy EC 2000HD





Custom Seating

