



Technical Bulletin

BULLETIN NO.:

TBUK030

DATE ISSUED

01/11/13

BULLETIN DETAILS

Product Group:	Coopers PVC Handle Elbow Crutches		
Model Affected:	8303C, 8249C, 8251C, 8253C & 9731C		
Effective Date:	November 2013	Serial No. Effective From:	N/A
Bulletin Subject:	Technical Update to Product		

1. Maximum user weight increase:

The Coopers PVC Handle Crutches have been retested and as a result the maximum user weight has been increased to 159Kg (previously 127kg).

No functional or dimensional specifications have changed on the product and performance is not affected in anyway. This update is purely an improvement to the weight carrying capacity of the Elbow Crutches through product retesting.

The product codes affected by this change are 8303C, 8249C, 8251C, 8253C & 9731C.



2. Product identification

In line with the introduction of this increased user weight all product labels have been updated with the new maximum user weight along with the user guide issued with each pair of crutches.

Products capable of carrying the new maximum user weight of 159Kg can be identified by the product label located under the crutch handle and on the outer packaging of the crutches. As in the example shown below.



If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88	+44 (0)1384 42 16 33	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk