



# Technical / News Bulletin

UK

BULLETIN NO.:

UK032

DATE ISSUED

14/03/2014

## BULLETIN DETAILS

**Product Group:** Power

**Model Affected:** Quickie Salsa M / R

**Effective Date:** Immediate effect      **Serial No. Effective From:** N/A

**Bulletin Subject:** Discontinuation of emergency freewheel option

### Subject of Bulletin:

#### Discontinuation of Low Push Force Emergency freewheel

As a result of customer feedback from April 1<sup>st</sup> 2014, the low push force Emergency freewheel will no longer be available on the Quickie Salsa R/M.

Final orders for this option can be placed up to and including the 31<sup>st</sup> March 2014.

If required, the emergency freewheel will still be available as a spare part however chairs in service fitted with the Emergency freewheel can be retro fitted with the standard drive wheels.

#### Improvements for the customer:

1. Simplified order forms
2. Simplified servicing using standard drive wheel
3. Reduced spares complexity

### Literature Affected:

**Sunparts:** No updates required

**Order Forms:** Updated to remove emergency freewheel option

**Website:** No change

**Owners Manual:** No change

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
<b>Telephone:</b>	<b>+44 (0)845 605 66 88</b> Option 2	<b>+44 (0)845 605 66 88</b> Option 3	enquiries@sunmed.co.uk
<b>Fax:</b>	<b>+44 (0)845 605 66 89</b>	<b>+44 (0)845 605 66 89</b>	help.technical@sunmed.co.uk