



# Technical / News Bulletin

UK

BULLETIN NO.:

0434

DATE ISSUED

14/3/2014

## BULLETIN DETAILS

<b>Product Group:</b>	Power		
<b>Model Affected:</b>	Quickie Salsa M / R2 and Jive R / R2 / F / M		
<b>Effective Date:</b>	Immediate effect	<b>Serial No. Effective From:</b>	N/A
<b>Bulletin Subject:</b>	Introduction of new Lift / Tilt module		

### Subject of Bulletin:

In an ongoing effort to improve the design, functionality and performance of our products the Lift & Tilt Module used on all Quickie and Zippie products has been improved to reduce the seat to floor height.

The reduction of the seat height when a Lift & Tilt option is selected aids with transfers, reduces head height requirements for vehicles and improves access under desks.

The new Lift & Tilt module provides a **20mm lower STFH** (seat to floor height) compared to the previous module. The seat height when a Lift/Tilt option is required is as follows:

- |                     |        |         |
|---------------------|--------|---------|
| • Quickie Salsa R2  | 47cm   | (18.5") |
| • Quickie R2 HD     | 43cm   | (16.5") |
| • Quickie Salsa M   | 48cm   | (19")   |
| • Quickie Salsa MHD | 44cm   | (17.3") |
| • Quickie Jive R    | 43cm   | (16.5") |
| • Quickie Jive F    | 44.5cm | (17.5") |
| • Quickie Jive M    | 44.5cm | (17.5") |

The new Lift/Tilt module is backwards compatible with Salsa M and all Jive variants.

The "old" Salsa (RWD) will change at a separate date to be announced.

### Literature Affected:

<b>Sunparts:</b>	Updated to include both old style and new style modules
<b>Order Forms:</b>	Seat to floor height tables updated to reflect changes in STFH
<b>Website:</b>	No change
<b>Owners Manual:</b>	No change

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
<b>Telephone:</b>	<b>+44 (0)845 605 66 88</b> Option 2	<b>+44 (0)845 605 66 88</b> Option 3	enquiries@sunmed.co.uk
<b>Fax:</b>	<b>+44 (0)845 605 66 89</b>	<b>+44 (0)845 605 66 89</b>	help.technical@sunmed.co.uk