

BULLETIN DETAILS

Product Group:	Power
Model Effected:	Quickie Salsa and Jive Families
Effective Date: 20.06.2014	Serial No. Effective From: N/A
Bulletin Subject:	Switch-It Mini Joystick and Powered swing-away chin control

Description of Bulletin:

As a result of the closure of the HMC factory in Belgium, Sunrise Medical has been required to source alternative solutions for all HMC components and as a result the HMC Mini Joystick and Powered Swing away chin control options have been upgraded.

1. Switch-It Mini JS (Micro guide):

The HMC Mini Joystick has been replaced by the Switch-It Mini Joystick (Micro Guide) on the Quickie Salsa and Jive powered wheelchairs. All mountings are compatible enabling existing chairs in the field to be easily upgraded to the new joystick without changing mounting systems.

All order Forms have been updated to reflect the new option and the Switch-It Mini JS (Micro guide) is available as a spare part under part code 21011083.

Benefits to customer

- Improved durability and easier mounting
- Increased compatibility with alternative mounting options e.g. bib mount, swivel joystick, bullet tray mount.



2. Powered Swing away chin control

The powered swing away control option on the Salsa and Jive powered wheelchairs. The option is backward compatible with all tubing, mounting dimensions and adjustment ranges remaining the same.



Powered swing-away mechanism



Stainless steel tubing

Benefit to customer

- Smaller actuator making it less obtrusive at the rear of the chair.

Literature affected:

- Parts Manual Updated

- Order Form Updated and available upon request from Customer Services or your Area Sales Manager

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 2	+44 (0)845 605 66 88 Option 3	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk