

## JIVE R<sup>2</sup>

HIGHENDPERFORMANCE  
WITH OUTSTANDING  
ADAPTABILITY





# PLEASE TAKE YOUR SEAT...



## **JIVE<sup>®</sup> R<sup>2</sup>**



The Jive R<sup>2</sup> is offered with the easily adjusted ASAP II seat frame. Featuring simple infield width, depth, armrest, legrest and back height adjustability. You can have confidence that the ASAP II seat frame will provide safe and comfortable support for years to come.

# EXPLORE YOUR ENVIRONMENT & GET IT ALL UNDER CONTROL



The Jive R<sup>2</sup> has impressive outdoor and indoor performance. The compact and durable base offers great stability and control. Powerful 4 pole motors can achieve up to 10 kph speed without compromising on manoeuvrability. The suspension allows comfortable climbing.

Jive is offered with **R-net**, highly adjustable electronics. Personalise your chair by assigning commonly used commands and settings to a single button.

For conversations on the move, keep tweets and status updates flowing by using the Jive's **Bluetooth** option to control your smartphone or tablet.

Or control your environment with Jive's **Omni** controller. Getting ready to watch a film? Dim the lights, select play and turn up the volume – all without your hand leaving the joystick. And, if a joystick's not suitable you can even customise the way you control your Jive with a variety of speciality drive controls.



## Ready for the trail:

Stylish 3 spoke wheels and 4 pole powerfull motors (speed up to 10 kph).

base available in blue, red & white

Easy battery access through the rear cover  
Simple battery and controller access while user is in the chair. Batteries available in 70 Ah.



## JIVE R<sup>2</sup> IS APPROVED FOR TRANSPORT IN A VEHICLE

The Jive R<sup>2</sup> is transit tested according to ISO 7176-19 and therefore approved for transportation in a vehicle. The tests are done for the 4 point tie-down system.



# THE JIVE R<sup>2</sup> OFFERS A WIDE RANGE OF OPTIONS CONFIGURE THE WHEELCHAIR TO YOUR INDIVIDUAL NEEDS – GET THE BEST FOR YOU!

FOR FURTHER INFORMATION PLEASE CHECK THE JIVE R<sup>2</sup> ORDER FORM.

## ASAP II Seat Frame:

- In-field width, depth, armrest, legrest, and back height adjustability
- Wide range of seat sizes to accommodate from paediatric to bariatric users
- Manual reclining backrest for increased function and sitting tolerance -option
- Optional fold down backrest for easy transport
- Extensive seating and positioning options
  - Multiple armrest options
  - Multiple legrest options
  - Long thigh supports
  - Fixed and articulating vent trays

## Power Options:

- 50 ° COG tilt
  - Super low seat height
- Combined tilt and lift module
  - Combination module offered at same price as tilt alone
  - Super low seat heights
  - 12" Seat elevation
  - 50 ° COG tilt
- Shear reduction recline
  - Integrated Jay 3 Back
  - Shear reduction for recline with minimal sliding
  - Super low seat heights
  - Adjustable seat depth

## TECHNICAL DATA

### JIVE R<sup>2</sup> power wheelchair

Max. user weight:	160 kg
Battery size:	70Ah
Caster wheel options:	9" and 10"
Max. dynamic stability:	10° (16%)
Drive wheel options:	14"
Overall width:	63 cm
Seat width:	12 -22"
Turning radius:	max. 120 cm
Seat depth:	12 -24"
Max. kerb climbing:	10 cm (with kerb climber)
Max. speed:	6 kph, 10 kph
Colour choices:	3 colours – red, blue & white
Max. range (ISO 7176-4):	(40 km with 70 Ah)
Transport (ISO 7176-19):	4 Point tie-down system
TGA Number:	100375



Sunrise Medical Australia Pty Ltd.

6 Healey Circuit

Huntingwood NSW 2148

Phone: (02) 9678 6600

Fax: (02) 9678 6656

Email: [orders@sunrisemedical.com.au](mailto:orders@sunrisemedical.com.au)

For further information on product options please check the order form on our website.

You can view this document also in a larger font on a PDF-file at [www.Sunrisemedical.com.au](http://www.Sunrisemedical.com.au)