



Technical / News Bulletin

UK

BULLETIN NO.:

TB 0442

DATE ISSUED

06/08/2014

BULLETIN DETAILS

Product Group:	Power		
Model Affected:	Quickie Tango, Rumba and Rumba Modular		
Effective Date:	23/06/2014	Serial No. Effective From:	N/A
Bulletin Subject:	New 8" Castor Wheel		

Subject of Bulletin:

As part of Sunrise Medicals commitment to improving our products and listening to quality and design feedback from our customer base, a new version of 8" castor wheel has been introduced.

The new wheel is used on the Quickie Rumba, Quickie Rumba Modular and Quickie Tango. It is retro-fittable and is available with immediate effect on production chairs and as a spare part. There are slight visual differences between the old and the new castor wheel, as shown below. The old version is no longer available, so where spares are required for existing products fitted with the older style wheel, it is recommended that both of the castor wheels should be replaced.



Old Castor wheel – C06141



New Castor wheel – C06278

Literature Affected:

Sunparts:	Updated with new part number
Order Forms:	No change
Website:	No change
Owners Manual:	No change

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 2	+44 (0)845 605 66 88 Option 3	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk