



Technical / News Bulletin

UK

BULLETIN NO.:

TB UK037

DATE ISSUED

04/02/2015

BULLETIN DETAILS

Product Group:	Power
Model Effected:	Quickie Salsa and Jive Families
Effective Date: 01.09.2014	Serial No. Effective From: N/A
Bulletin Subject:	Special Controls Update

Description of Bulletin:

As has been stated in a previous Technical Bulletin, the closure of the HMC factory in Belgium has led to Sunrise Medical sourcing alternative solutions for all HMC components.

Sunrise Medical has already introduced the Switch-It Micro Pilot and the Switch-It Micro Guide as replacement items for the HMC Mini Joystick. We now announce the introduction of the Switch-It Versaguide as a direct replacement for the HMC Compact Joystick. The Versaguide can be used as either hand or chin control, and offers the same benefits to the customer as its predecessor.



This introduction completes the update on standard special controls that are available directly via the order form on Salsa and Jive models.

For information regarding additional special control inputs such as foot controls, tablet drives etc, please contact your local Area Sales Manager or speak to the Technical Service Centre on the number detailed below.

Literature affected:

- Parts Manual	Updated
- Order Form	Updated and available upon request from Customer Services or your Area Sales Manager

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 2	+44 (0)845 605 66 88 Option 3	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk