



Technical / News Bulletin

UK

BULLETIN NO.:

TB 0457

DATE ISSUED

04/02/2015

BULLETIN DETAILS


Product Group:	Manual Wheelchairs		
Model Affected:	Quickie Xenon Fixed Front and Quickie Xenon Swing-away		
Effective Date:	Immediate	Serial No. Effective From:	Implemented since February 2014
Bulletin Subject:	Change to upholstery coating		

Description of issue

After a period in the field, and when subjected to direct sunlight or friction between the cushion and the upholstery, some upholsteries were becoming grey and fading

Solution – Refer to image below

A more durable coating has been applied to the upholstery to prevent fading



old coating

new coating

Benefits to customer

1. More durable upholstery
2. Upholstery does not fade
3. Chair keeps aesthetic appearance

Literature Affected:

Sunparts:	Running change. No updates required
Order Forms:	No changes required
Website:	No changes required
Owners Manual:	No changes required

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 2	+44 (0)845 605 66 88 Option 3	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk