



Technical / News Bulletin

UK

BULLETIN NO.:

TB 0460

DATE ISSUED

04-01-15

BULLETIN DETAILS

Product Group:	Manual Wheelchairs		
Model Affected:	Breezy Relax ²		
Effective Date:	Immediate	Serial No. Effective From:	Implemented since April 2014
Bulletin Subject:	Improvement to headrest bracket		

Description of issue

A small number of headrest brackets have been subject to breakages in the field.

Solution – Refer to image below

The material for the receiver of the standard and three wing headrests has been changed following extensive testing. Additionally, two extra welding points have been added to increase the strength by 25%.



Benefit to customer

1. Improved strength and protection against breakages.

Literature Affected:

Sunparts:	Running change. No updates required
Order Forms:	No changes required
Website:	No changes required
Owners Manual:	No changes required

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 2	+44 (0)845 605 66 88 Option 3	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk