



Technical / News Bulletin

UK

BULLETIN NO.:

TB 0468

DATE ISSUED

01/11/2015

BULLETIN DETAILS

Product Group:	Power
Model Effected:	Salsa M2
Effective Date: 30-04-2015	Serial No. Effective From: N/A
Bulletin Subject:	Salsa M2 with improved 7" castor wheel

Description of Bulletin:

As part of Sunrise Medicals after sales feedback mechanism, we have identified that the 7" castor wheel tyre used on the Salsa M2 is too soft and wears out too quickly.

We are now introducing, effective from 30.04.2015, a new updated 7" castor wheel with a harder type tyre which matches the existing 7" castor wheel as used on the Salsa M.

Spare wheels will also be exchanged with the updated wheel effective from the same date. Part number will not change.

The castor wheel is retro-fittable to existing Salsa M2 wheelchairs in the field.

Literature affected:

- Parts Manual	Updated
- Order Form	N/A

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 2	+44 (0)845 605 66 88 Option 3	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk