

## BULLETIN DETAILS

<b>Product Group:</b>	Manual
<b>Model Effected:</b>	Helium, Xenon <sup>2</sup> , Argon <sup>2</sup> , All Life and Life Teens
<b>Effective Date: 01.09.2014</b>	<b>Serial No. Effective From:</b> Various
<b>Bulletin Subject:</b>	Update to Exo backrest upholstery

### Description of Bulletin:

As part of the move to improve product options on Sunrise wheelchairs, the standard backrest options on the Active Wheelchair range has been updated to the new Exo backrest style upholstery.



The **Quickie Exo backrest** is available as standard on the Quickie Active wheelchair range and replaces the older traditional backrest upholstery. The optional **Quickie Exo Pro backrest** is made with 3DX<sup>TM</sup> vented fabric for extra ventilation.

The **Quickie Exo backrest** is an evolution in wheelchair upholstery design, incorporating technology from the Formula 1 and adrenaline sports industries onto your ultra-lightweight wheelchair and improving your everyday performance. Key features of the new style backrests are:

Slimline Design  
 Additional support and comfort  
 Lightweight  
 Easy Adjustment

Spare parts manuals have been updated to show the new backrest upholstery, which is completely retro-fittable to older chairs in the product range.

### **Literature affected:**

- Parts Manual	Updated
- Order Form	Updated and available upon request from Customer Services or your Area Sales Manager

**If you have any questions regarding this bulletin, please do not hesitate to contact:**

Customer Services	Customer Support	Technical Service Centre	E-mail
<b>Telephone:</b>	<b>+44 (0)845 605 66 88</b> <b>Option 2</b>	<b>+44 (0)845 605 66 88</b> <b>Option 3</b>	<a href="mailto:enquiries@sunmed.co.uk">enquiries@sunmed.co.uk</a>
<b>Fax:</b>	<b>+44 (0)845 605 66 89</b>	<b>+44 (0)845 605 66 89</b>	<a href="mailto:help.technical@sunmed.co.uk">help.technical@sunmed.co.uk</a>