



# Technical / News Bulletin

UK

BULLETIN NO.:

TB UK055

DATE ISSUED

01/11/2015


## BULLETIN DETAILS

<b>Product Group:</b>	Power
<b>Model Effected:</b>	Jive F/R2/M/FXL and Salsa HD family
<b>Effective Date: 01-07-2015</b>	<b>Serial No. Effective From:</b> N/A
<b>Bulletin Subject:</b>	Larger Lap belt option introduced

### Description of Bulletin:

As part of Sunrise Medicals response to customer feedback, we are delighted to offer the introduction of a longer lap belt option for Jive and Salsa HD wheelchairs. The longer lap belt has a length of 190cm and caters for the wider seat widths available on these products.

The belt is an aircraft buckle style as per the standard length belt that was originally the only option on the products



The shown drawing has no interactivity.

Pos. ▾	Item number ▾	Description ▾	get price
1	013905	STD AIRCRAFT BUCKLE, 130CM, INC FIXINGS	
1	20009064-190	STD AIRCRAFT BUCKLE, 190 CM, INC FIXINGS	
2	013907	PADDED AIRCRAFT BUCKLE, INC FIXINGS	

The new lapbelt has been added to the parts manuals for the affected products as shown above

### Literature affected:

- Parts Manual Updated
- Order Form N/A

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
-------------------	------------------	--------------------------	--------

<b>Telephone:</b>	<b>+44 (0)845 605 66 88 Option 2</b>	<b>+44 (0)845 605 66 88 Option 3</b>	<a href="mailto:enquiries@sunmed.co.uk">enquiries@sunmed.co.uk</a>
<b>Fax:</b>	<b>+44 (0)845 605 66 89</b>	<b>+44 (0)845 605 66 89</b>	<a href="mailto:help.technical@sunmed.co.uk">help.technical@sunmed.co.uk</a>