
S	U	N	I	2	S	E	
М	Е	D	I	С	А	L.	

Technical / News Bulletin



	DATE ISSUED					
	01/03/2016					
	BULLETIN DETAILS					
Product Group:	Manual					
Model Affected:	GPU, GPV & GP SA					
Effective Date: 1 st May 2016	Serial No. Effective From: N/A					
Bulletin Subject:	GPU, GPV & GP SA discontinuation					
Description of Bulletin:						
Technical Bulletin UK059 is confirmation GPV and Quickie GP SA commencing	on of the official discontinuation of the Quickie GPU, Quickie from 1 st May 2016.					
Since the introduction of the Quickie Li	d in functionality and features by the Quickie Life R in 2013. ife R, the demand for the GPU, GPV and GP SA has ched a point where the manufacturing for the GPU, GPV and					
The GPU, GPV and GP SA will continue to be available up to 1 st May 2016, after this date orders will no longer be processed.						
Spares or equivalents will be maintained and supplied for 5 years for the GPU, GPV and GP SA.						
• •	n the Quickie Life R, which is available with a fixed front and or more information on the benefits of switching to the					
Discontinued GPU, GPV & GP SA	Replacement Quickie Life R					

(GPV shown above)

(Fixed front shown above)

Benefits to customers moving to the more modern Quickie Life R:

The Quickie Life R brings a significant number of benefits to the customer:

1. The Life R is available in both a fixed front frame or with swing away frame.

- 2. The standard max. user weight for the Life R is 125kg (19.6 stones) which is greater than that of the GPU, GPV or GP SA.
- 3. It has a seat depth range of 36 50cm and seat width of 34 48cm. There is also the possibility of a smaller seat depth and width with the Quickie Life RT (teen version).
- 4. The Life R has a more modern style to the frame and the options and accessories that are available on the Life R have been rolled down from our flag ship products. For example the Exo backrest upholstery is a standard on the Life R, this upholstery allows greater breathability and comfort to customers, as well as better styling.
- 5. There is an extensive variety of colour options available on the Life R, from 32 frame colour choices to 5 upholstery trim colours. Customers are also able to customise their chair to their own unique style with the non-cost anodised fork and axle stem.
- 6. It is lighter than the GPV with the weight of the Life R from just 10.5kg.
- 7. There is a 5 year frame warranty and a 2 year guarantee on components on the Life R.
- 8. The chair is much easier to adjust with a new angle backrest plate which is also stronger and more durable. The camber can be adjusted from 0-4 degrees with no extra parts.
- 9. The Life R has been successfully transit tested.

Literature affected:

- Parts Manual N/A

- Order Form GPU, GPV & GP SA to be removed.

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail	
Telephone:	+44 (0)845 605 66 88 Option 2	+44 (0)845 605 66 88 Option 3	enquiries@sunmed.co.uk	
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed .co.uk	