



Technical / News Bulletin  
UK

**BULLETIN NO.:**

**TB 0507**

**DATE ISSUED**

25/07/16

**BULLETIN DETAILS**

<b>Product Group:</b>	Rehab Power	
<b>Model Affected:</b>	Quickie Jive Up	
<b>Effective Date: 11<sup>th</sup> July 2016</b>	<b>Serial No. Effective From:</b>	N/A
<b>Bulletin Subject:</b>	Jive Up 5 button switch box	

**Description of Issue / Solution:** In line with the acquisition of Switch-It, we have moved to their option when choosing the 5 button switch box on the Jive Up. This offers end users the additional benefits listed below.

**Old Style  
(Button)**



**New Button**



**New Toggle**



**USB Port**

**Benefits to the customer**

- No price difference
- USB port as standard (for charging devices)
- LED Light feedback
- LED light feedback
- Up to two switch boxes per chair
- Can be retro fitted to all R-net systems

**Literature affected:**

- Parts Manual - To be updated

**If you have any questions regarding this bulletin, please do not hesitate to contact:**

Customer Services	Customer Support	Technical Service Centre	E-mail
<b>Telephone:</b>	<b>+44 (0)845 605 66 88 Option 2</b>	<b>+44 (0)845 605 66 88 Option 3</b>	<a href="mailto:enquiries@sunmed.co.uk">enquiries@sunmed.co.uk</a>
<b>Fax:</b>	<b>+44 (0)845 605 66 89</b>	<b>+44 (0)845 605 66 89</b>	<a href="mailto:help.technical@sunmed.co.uk">help.technical@sunmed.co.uk</a>