



# Technical / News Bulletin

UK

BULLETIN NO.:

TB UK064

DATE ISSUED

15/08/2016

## BULLETIN DETAILS

**Product Group:** Seating

**Model Affected:** JAY GS Cushion

**Effective Date:** Immediate effect

**Serial No. Effective From:** N/A

**Bulletin Subject:** Improvements to GS Cushion

### Subject of Bulletin:

In an ongoing effort to improve the design, functionality and performance of our JAY GS (Growth System) cushion, a series of improvements have been introduced with immediate effect. These include a new fluid pad option, new style labels, Velcro attachment configuration of accessories onto base foam, new order form which utilises metric sizing and JAY SMART part numbers and a new GS accessory box.

#### Feature

##### New 'balance' Flow pad option



New 'balance' GS fluid pad



Old style 'Standard' fluid pad



Old Style 'Pressure Relief' pad

#### Benefit

Designed to meet the pressure care needs of clients who also require high skin protection and aggressive lower limb positioning;

- ✓ New "balance" pad is the same as the old Pressure Relief pad but **under-filled by 60%**
- ✓ It offers the **best compromise** between more aggressive positioning and high skin protection
- ✓ **3D construction** of the fluid pad allows immersion and envelopment around bony prominences for greater pressure relief
- ✓ Allows the user to '**sit in**' the cushion rather than '**on top**' when multiple positioning accessories are required

The old style 'Standard' and 'Pressure relief' pads are still available. Prescribers can now choose one of 3 pads to meet user need.

#### When would I use the standard fill fluid pad?

- 2D construction
- Excellent for aggressive positioning of the user for whom pressure is less of an issue e.g. paediatric users

#### When would I use the Pressure relief pad?

- 3D construction
- Large amount of JAY fluid in each segment
- Offers higher skin protection for those

users at higher risk of skin breakdown who do not require leg positioning accessories

**New style label**



A new style logo and rubber labelling gives the JAY GS cushion a distinct, clearly identifiable and modern look.

**New orientation of Velcro on base foam**



The foam base now has the soft loop Velcro attached to it instead of the old hook style.

**Benefit:**

- Allows easier attachment and repositioning of the accessories as the user grows. It is also more comfortable for the user to sit on when less accessories are used.

**New option of oversized covers**

The new JAY GS now has an **option of Oversized Incontinence or Microclimatic covers** which readily accept 2 or 3" positioning accessories:

**Benefit:**

- The Shower cap design of all JAY GS cushion covers allows pelvic and thigh positioning components to be added or removed as required
- The oversized cover option is designed to accommodate the components without over-stretching or 'hammocking' the cover

**New JAY SMART ordering code**

	Width	Depth	Fluid type	Outer Cover	fixed
JGSC	25	33	S	I	EU

Jay GS Cushion

The JAY GS is ordered using a **new Order Form and JAY Smart code** (similar to the JAY Balance and JAY Xtreme Active cushions).

**Benefits:**

- ✓ Metric measurements to more accurately match wheelchair width and depths
- ✓ More accurately meet user need by modularly constructing the cushion – both fluid pad and cover type
- ✓ More easily identify the construction of the cushion

	<p><i>The cushions are ordered in metric sizes but these will be the same as those currently used.</i></p> <p><i>E.g. JGSC2533SIEU = same as old JAY GS 10"x13"with standard fluid pad and incontinence cover</i></p> <p><b>*** Please see the accompanying <i>JAY GS Cushion Conversion Matrix</i> to aid selection of like-for-like replacements</b></p> <p><i>New order forms are available via customer Services or from our website at <a href="http://www.sunrisemedical.co.uk">www.sunrisemedical.co.uk</a></i></p>
<b>Increased User Weight Limit</b>	<p>In keeping with new cushions within the JAY cushion range, the maximum user weight has been <b>increased to 150kg</b> as standard.</p> <p>Benefit:</p> <ul style="list-style-type: none"> <li>✓ Provides an aggressive positioning cushion solution for heavier users over 18 stone (113kg).</li> </ul>
<b>New JAY GS Accessory Box</b>	<p>A new style JAY GS Accessory box featuring the full range of JAY GS positioning accessories is available.</p> <p>Benefit:</p> <ul style="list-style-type: none"> <li>✓ Ideal for use at assessment to establish the exact positioning accessories required for optimal pelvic and lower limb support</li> </ul>

Please refer to [www.sunrisemedical.co.uk](http://www.sunrisemedical.co.uk), or contact your area sales manager for more information

If you require any further technical assistance please contact:

Matthew Eveleigh, JAY Clinical Specialist on [matthew.eveleigh@sunmed.co.uk](mailto:matthew.eveleigh@sunmed.co.uk)

**Literature Affected:**

<b>Sunparts:</b>	Available from <a href="http://parts.sunrisemedical.eu">http://parts.sunrisemedical.eu</a>
<b>Order Forms:</b>	Updated and available via Customer Services
<b>Website:</b>	Available on <a href="http://www.sunrisemedical.co.uk">www.sunrisemedical.co.uk</a>
<b>Owners Manual:</b>	Available on <a href="http://www.sunrisemedical.co.uk">www.sunrisemedical.co.uk</a>

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
<b>Telephone:</b>	<b>+44 (0)845 605 66 88 Option 2</b>	<b>+44 (0)845 605 66 88 Option 3</b>	<a href="mailto:enquiries@sunmed.co.uk">enquiries@sunmed.co.uk</a>
<b>Fax:</b>	<b>+44 (0)845 605 66 89</b>	<b>+44 (0)845 605 66 89</b>	<a href="mailto:help.technical@sunmed.co.uk">help.technical@sunmed.co.uk</a>

# JAY GS Cushion Conversion Matrix

Please use matrix to convert a part code for old JAY GS cushion to the same specification of the new JAY GS cushion

Old style JAY GS Cushion Code	Description	New JAY GS Cushion Code	Description
JAY1111.1NF	JAY GS CUSHION 10x13" with pressure relief pad & incontinence cover	JGSC2533PRIEU	J GS CUSH 25X33 PR INC
JAY1111NF	JAY GS CUSHION 10x13" with standard pad & incontinence cover	JGSC2533SIEU	J GS CUSH 25X33 S INC
JAY1115.1NF	JAY GS CUSHION 12x15" with pressure relief pad & incontinence cover	JGSC3038PRIEU	J GS CUSH 30X38 PR INC
JAY1115NF	JAY GS CUSHION 12x15" with standard pad & incontinence cover	JGSC3038SIEU	J GS CUSH 30X38 S INC
JAY1110.1NF	JAY GS CUSHION 14x17" with pressure relief pad & incontinence cover	JGSC3643PRIEU	J GS CUSH 36X43 PR INC
JAY1110NF	JAY GS CUSHION 14x17" with standard pad & incontinence cover	JGSC3643SIEU	J GS CUSH 36X43 S INC
JAY1105.1NF	JAY GS CUSHION 16x19" with pressure relief pad & incontinence cover	JGSC4048PRIEU	J GS CUSH 40X48 PR INC
JAY1105NF	JAY GS CUSHION 16x19" with standard pad & incontinence cover	JGSC4048SIEU	J GS CUSH 40X48 S INC
JAY1100.1NF	JAY GS CUSHION 18x20" with pressure relief pad & incontinence cover	JGSC4550PRIEU	J GS CUSH 45X50 PR INC
JAY1100NF	JAY GS CUSHION 18x20" with standard pad & incontinence cover	JGSC4550SIEU	J GS CUSH 45X50 S INC
JAYB1111FR	FOAM BASE 10x13"	BJGS2533	BASE JAY GS 2533 VELCRO
JAYB1115FR	FOAM BASE 12x15"	BJGS3038	BASE JAY GS 3038 VELCRO

<b>JAYB1110FR</b>	FOAM BASE 14x17"	<b>BJGS3643</b>	BASE JAY GS 3643 VELCRO
<b>JAYB1105FR</b>	FOAM BASE 16x19"	<b>BJGS4048</b>	BASE JAY GS 4048 VELCRO
<b>JAYB1100FR</b>	FOAM BASE 18x20"	<b>BJGS4550</b>	BASE JAY GS 4550 VELCRO
<b>JAYC1111NF</b>	Incontinence resistant cover 10x13"	<b>CJGS2533I</b>	CVR JAY GS INC 25X33
<b>JAYC1115NF</b>	Incontinence resistant cover 12x15"	<b>CJGS3038I</b>	CVR JAY GS INC 33X38
<b>JAYC1110NF</b>	Incontinence resistant cover 14x17"	<b>CJGS3643I</b>	CVR JAY GS INC 36X43
<b>JAYC1105NF</b>	Incontinence resistant cover 16x19"	<b>CJGS4048I</b>	CVR JAY GS INC 40X48
<b>JAYC1100NF</b>	Incontinence resistant cover 18x20"	<b>CJGS4550I</b>	CVR JAY GS INC 45X50
<b>JAYF1111.1</b>	Pressure relief pad 10x13"	<b>FPRJGS25</b>	PR FLUID PAD JGS 25
<b>JAYF1115.1</b>	Pressure relief pad 12x15"	<b>FPRJGS30</b>	PR FLUID PAD JGS 30
<b>JAYF1110.1</b>	Pressure relief pad 14x17"	<b>FPRJGS36</b>	PR FLUID PAD JGS 36
<b>JAYF1105.1</b>	Pressure relief pad 16x19"	<b>FPRJGS40</b>	PR FLUID PAD JGS 40
<b>JAYF1100.1</b>	Pressure relief pad 18x20"	<b>FPRJGS45</b>	PR FLUID PAD JGS 45
<b>JAYF1111</b>	Standard fluid pad 10x13"	<b>FSJGS25</b>	STD FLUID PAD JGS 25
<b>JAYF1115</b>	Standard fluid pad 12x15"	<b>FSJGS30</b>	STD FLUID PAD JGS 30

<b>JAYF1110</b>	Standard fluid pad 14x17"	<b>FSJGS36</b>	STD FLUID PAD JGS 36
<b>JAYF1105</b>	Standard fluid pad 16x19"	<b>FSJGS40</b>	STD FLUID PAD JGS 40
<b>JAYF1100</b>	Standard fluid pad 18x20"	<b>FSJGS45</b>	STD FLUID PAD JGS 45
<b>JAYB118FR</b>	2" ABDUCTOR BUILD-UP	<b>JGSABDBUP2</b>	2" ABDUCTOR BUILD-UP
<b>JAYB1118FR</b>	3" ABDUCTOR BUILD-UP	<b>JGSABDBUP3</b>	3" ABDUCTOR BUILD-UP
<b>JAYB118WFR</b>	2" ABDUCTOR WEDGE	<b>JGSABW2</b>	2" ABDUCTOR WEDGE
<b>JAYB119JFR</b>	1" ADDUCTOR WEDGES (PR)	<b>JGSADDW2</b>	1" ADDUCTOR WEDGES (PR)
<b>JAYB119JFR</b>	2" ADDUCTOR WEDGES (PR)	<b>JGSADDW5</b>	2" ADDUCTOR WEDGES (PR)
<b>JAYB11113FR</b>	3" ADDUCTOR WEDGES (PR)	<b>JGSADDW7</b>	3" ADDUCTOR WEDGES (PR)
<b>JAY1155FR</b>	JAY GS ACCESSORY BOX	<b>JGSBOX</b>	JAY GS ACCESSORY BOX
<b>JAYB163FR</b>	2" CHILD ADDUCTOR WEDGES	<b>JGSCADDW</b>	2" CHILD ADDUCTOR WEDGES
<b>JAYB1167FR</b>	2" HIP GUIDES (PR)	<b>JGSHG2</b>	2" HIP GUIDES (PR)
<b>JAYB118FR</b>	2" JR/CHILD ABDUCTOR	<b>JGSJABDBUP2</b>	2" JR/CHILD ABDUCTOR
<b>JAYF1119</b>	MINI FLUID SUPPLEMENT PAD	<b>JGSMFSP</b>	MINI FLUID SUPPLEMENT PAD
<b>JAYB1171FR</b>	PRE-ISCHIAL SHELF EXTENDERS 10x13"	<b>JGSPISEXT2533</b>	PRE-ISCHIAL SHELF EXTENDERS 25x33

<b>JAYB1173FR</b>	PRE-ISCHIAL SHELF EXTENDERS 12x15"		<b>JGSPISEXT3038</b>	PRE-ISCHIAL SHELF EXTENDERS 30x38
<b>JAYB1174FR</b>	PRE-ISCHIAL SHELF EXTENDERS 14x17"		<b>JGSPISEXT3643</b>	PRE-ISCHIAL SHELF EXTENDERS 36x46
<b>JAYB1175FR</b>	PRE-ISCHIAL SHELF EXTENDERS 16x19"		<b>JGSPISEXT4048</b>	PRE-ISCHIAL SHELF EXTENDERS 40x48
<b>JAYB1176FR</b>	PRE-ISCHIAL SHELF EXTENDERS 18x20"		<b>JGSPISEXT4550</b>	PRE-ISCHIAL SHELF EXTENDERS 45x50
<b>JAYK1115</b>	1/2" WELL BUILD-UP 12X15"		<b>JGSWBUP3038</b>	1/2" WELL BUILD-UP 30X38
<b>JAYK1110</b>	1/2" WELL BUILD-UP 14X17"		<b>JGSWBUP3643</b>	1/2" WELL BUILD-UP 36X43
<b>JAYK1105</b>	1/2" WELL BUILD-UP 16X19"		<b>JGSWBUP4048</b>	1/2" WELL BUILD-UP 40X48
<b>JAYK1100</b>	1/2" WELL BUILD-UP 18X20"		<b>JGSWBUP4550</b>	1/2" WELL BUILD-UP 45X50