



Technical / News Bulletin

UK

BULLETIN NO.:

TB UK066

DATE ISSUED

20/12/2016

BULLETIN DETAILS


Product Group:	Manual Wheelchairs		
Models Affected:	Life FT, Life RT, Youngster 3, Simba		
Effective Date:	November 15 th 2016	Serial No. Effective From:	N/A
Bulletin Subject:	Safari Brake Lever Update		

Subject of Bulletin: Changes to the Safari Brake Lever design


Description:

- Based on customer feedback we have changed the design of the 'Safari Brake' lever
- Customers expressed concern that the levers were too long and interfered with the user's arms
- We have reduced the length of the lever by 10mm, this change will be applied to all chairs utilizing this product with immediate effect and with no change to part numbers or order forms.

OLD



NEW



Benefits now given to the Customer:

- Greater freedom of movement resulting in longer and more efficient propulsion strokes
- Less interference when completing transfers

Please refer to www.sunrisemedical.co.uk, or contact your area sales manager for more information

Literature Affected:

Website: Available on www.sunrisemedical.co.uk

Owners Manual: Available on www.sunrisemedical.co.uk

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 2	+44 (0)845 605 66 88 Option 3	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk