



ZIPPIE SALSA FAMILY

POWERED
WHEELCHAIRS
THAT FIT YOUR CHILD



POWERED WHEELCHAIRS THAT FIT YOUR CHILD



Whether it's navigating around the home or exploring outdoors, Zippie's range of power wheelchairs have just as much energy for adventure as your child! Available with adjustable seating systems, powered options and control systems to interact with the surrounding world – they're electric wheelchairs designed for your child's independence.



ZIPPIE *SALSA m²*

// The power wheelchair for children with indoor manoeuvrability and outdoor performance.



ZIPPIE *SALSA m²*

// The narrowest mid-wheel drive powered wheelchair on the market!



ZIPPIE *SALSA R²*

// The compact and powerful wheelchair with outstanding outdoor performance.



ZIPPIE ***SALSA M²***

The power wheelchair for kids with superb indoor manoeuvrability and outdoor performance.

The ZIPPIE Salsa M²'s mid wheel drive base delivers the best possible indoor performance thanks to its ultracompact turning radius (670 mm) and small footprint. A narrow width of just 610 mm and low seat height make access very easy. Performing as good outdoors as it does indoors, the unique suspension system delivers a smooth, safe and stable ride even when travelling over slopes and cambers. The ZIPPIE Salsa M² powered wheelchair can overcome obstacles of up to 75 mm and can travel up to 32 km on a full charge.

AGILE. COMPACT. VERSATILE.

// 60 AH BATTERIES

Range of 32 km according to the ISO standards.

// 2 OR 4 POLE MOTORS

Performance motors with 6 and 10 kph options.

// ANTI-PITCH & SUSPENSION

Enjoy a smooth ride over uneven terrain with minimal body movement.



// FOUR BODY COLOUR OPTIONS

blue, red, white & matt black



ZIPPIE *SALSA M² Mini*

Say no more to compromising on indoor manoeuvrability and outdoor performance thanks to the new Zippie Salsa M² Mini.

Introducing the world's narrowest mid-wheel drive performance powerchair with a chassis width of only 52 cm – That's 9 cm smaller than our market leading Zippie Salsa M².

Unlike other indoor powerchairs you don't have to compromise on outdoor performance, the Salsa M² Mini uses exactly the same drive base as the Salsa M², just narrower – so all of the great features you get on a bigger outdoor powerchair are packed into the world's narrowest mid-wheel drive powerchair.



THE NARROWEST MID-WHEEL DRIVE

// 36 AH BATTERIES

Range of 23 km according to the ISO standards.

// TIE DOWN

4 point tie-down system for transportation in a vehicle. Optional with Dahl car docking system available.

// ANTI-PITCH & SUSPENSION

Enjoy a smooth ride over uneven terrain with minimal body movement.



only
52 cm



// FOUR BODY COLOUR OPTIONS

blue, red, white & matt black



ZIPPIE ***SALSA R²***

The compact and powerful wheelchair with outstanding outdoor performance.

Thanks to the ZIPPIE Salsa R²'s powerful motors, 70 Ah batteries and chunky wheels, it's a powered wheelchair with outstanding outdoor performance. Capable of climbing over obstacles of up to 100 mm and tackling inclines of up to 10°, the Salsa-R² can travel 32 km in comfort. In addition, the Zippie Salsa R²'s ultra-compact base and super-low seat height make access into vehicles, around confined spaces and under tables easy.

GO OUT AND EXPLORE THE WORLD



// SUSPENSION

Adjustable drive wheel suspension



// TIE DOWN

4 point tie-down system for transportation in a vehicle. Optional with Dahl car docking system available.



// STYLISH AND COMPACT BASE

in 3 colours available –
blue, red & white

// ACTIVE SWING AWAY HANGERS

Swing away hangers for easy transfer. Also optional centre mount footrest available.

// LED LIGHTS

New base integrated LED lights and indicators (optional)





ZIPPIE SALSAS SEATING SYSTEM

The Zippie Salsa is equipped with a fully adjustable seating system which can be quickly adjusted throughout the full range of width and depth, back angle and height without additional parts. The seat is designed for perfect adaptation to the growing child. Alternative backs and cushions can be easily fitted. The Zippie Salsa seating system can be fitted to both Salsa M² and Salsa R².

Swing-away and removable tray table



Standard flip up armrests for easy accessibility.



Powered lift and tilt module with centre of gravity adjustment to suit correct users balance

Comfort Seating System (optional) for good positioning.



TECHNICAL DATA

ZIPPIE SALSA M² / ZIPPIE SALSA M² MINI / ZIPPIE SALSA R² // PEDIATRIC POWERED WHEELCHAIRS

Max. user weight	75 kg	Max. reversing width	Salsa M ² & Salsa M ² Mini: 670 mm / Salsa R ² : 910 mm
Seat width	300 – 400 mm (in 20 mm steps)	Max. safe slope	8°
Seat depth	300 – 400 mm (in 20 mm steps)	Max. kerb climbing	Salsa M ² & Salsa M ² Mini: 70 mm / Salsa R ² : 100 mm
Seat height	Salsa M ² & Salsa M ² Mini: 400 – 455 mm Salsa R ² : 415 – 470 mm	Castor wheel options	Salsa M ² & Salsa M ² Mini: 13" drive wheels; 7" castor Salsa R ² : 12,5" or 14" drive wheels; 9" castor
Seat angle	0° – 30°; -5° to 25° (optional)	Colour inlay choices	Salsa M ² & Salsa M ² Mini: Red, Blue, White, Matt Black Salsa R ² : Red, Blue, White
Backrest height	410 mm (± 25 mm)	Max. range (ISO 7176-4)	Salsa M ² & Salsa R ² : 32 km Salsa M ² Mini: 23 km (dependant on distance profile, user weight & electrical options)
Backrest angle	86° – 106° (mechanical); 90° – 120° (powered)	Wheelchair weight	from 97 kg (Salsa M ² Mini) 120 kg (Salsa M ² & Salsa R ²); including batteries
Total base width	Salsa M ² : 600 mm Salsa M ² Mini: 520 mm Salsa R ² : 620 mm (14" wheels) / 600mm (12,5" wheels)	Crash tested	Yes (approved for transport in a vehicle according ISO 7176-19)
Total length	Salsa M ² & Salsa M ² Mini: 1120 mm / Salsa R ² : 1070 mm	Intended use and environment	For people with limited mobility; self-propelling or being pushed by an assistant, indoor and outdoor use. Wheelchair type class B
Total height	Salsa M ² & Salsa M ² Mini: 925 mm / Salsa R ² : 945 mm		
Max. speed	6 kph & 10 kph (optional)		
Battery size	Salsa M ² 60 Ah / Salsa M ² Mini 41 Ah Salsa R ² 60 Ah & 70 Ah (optional)		
Electronics	R-net & VR2 controller; lights & indicators		



Sunrise Medical GmbH & Co. KG

Kahlbachring 2-4
69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: international@sunrisemedical.de
www.SunriseMedical.eu



For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.
If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu.

