

## Little Gem<sup>®</sup> 2



# Little Gem<sup>®</sup> 2

12. 23

Let the Little Gem 2 make you smile!



The Little Gem 2 has been designed with your comfort and safety in mind – with Sterling Lift the conveniently located handles make lifting and handling the parts easier than ever before.

#### So easy...

To get your Little Gem 2 ready for transporting, simply lift off the seat, lift out the battery and lift the rear handle releasing the Sterling Lock that quickly separates the front and rear of the scooter. It really is that simple! And by the way, it goes back together just as smoothly!





#### So comfortable...

The Sterling Comfort package features a soft padded seat and generous legroom; we know you'll never want to travel without your Little Gem 2 again!





#### So accessible...


With an infinitely adjustable tiller that can be pushed up and away from you, flip-up armrests and a swivelling seat, you can get on and off your Little Gem 2 with ease. Your Little Gem 2 is so compact and manoeuvrable you can take it almost anywhere.





### Little Gem<sup>®</sup> 2 will make you smile!





#### Technical data sheet

Version	4 Wheels	Drive unit weight	10.7kg
Speed	6.4kph (4mph)	Weight of the heaviest part	15kg
Batteries	12 Ah	Turning radius	1150mm (45.3")
Maximum range	16km (10 miles)	Maximum safe slope	6°
Colour	Blue	Ground clearance	60mm
Maximum user weight	113kg (17.8 stone)	Tyres	Solid
Total weight	48/106kg	Front wheel size	175mm (7")
Batteries weight	4.31kg unit	Rear wheel size	190mm (7.5")
Seat weight	10.4kg	Lights and indicator package	No

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.co.uk or alternatively is available on request in large text.



#### Sunrise Medical Limited

Thorns Road Brierley Hill West Midlands DY5 2LD Tel: +44 (0) 845 605 66 88 Fax: +44 (0) 845 605 66 89