



## Quality and reliability

Sunrise produce the highest quality components at an affordable price. The key to succeeding at this has been our parts commonality. Breezy wheelchairs use key components and superior designs from our elite line of Quickie chairs. Castor housings, armrests, backrest posts and footrests are just some of the components that are integrated into the Breezy SL/TL.



High quality T in T extension tubes ensure absolute secure foot support and quick accurate height adjustment.



The backrest is height adjustable for correct posture, comfort and freedom of movement (using a telescoping tube in the rear frame)



Armrests flip-up for easy transfer and close in the blink of an eye

#### Breezy transit model available





8" x 2" Pneumatic or solid castor (option)




# Breezy SL/TL


Because quality of life matters

## Breezy SL / TL



40 - 50cm (16 - 20")





40 - 47.5cm (15.7 - 18.6")







TL 94cm (37") SL 103cm (40.5")



### Technical data sheet

Frame colour:	Silver texture, blue, forest green, red	Turning circle:	76.5cm (30")
Upholstery colour:	Black	Folded width:	31.4cm (12.3")
Approved user weight:	125kg (19.6 stone)	Folded height:	91.2cm (35.7")
Max. weight:	15kg	Folded length:	TL 94cm (37") / SL 103cm (40.5")
Weight of heaviest part:	13.5kg	Transit tested to ISO 7176 part 19 an	d passed Yes

For further information on the full specification, options and accessories, please refer to the order form.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.co.uk or alternatively is available on request in large text.



Sunrise Medical Limited Thorns Road Brierley Hill West Midlands DY5 2LD

Tel: +44 (0)845 605 66 88 Fax: +44 (0)845 605 66 89