



Technical / News Bulletin  
UK

**BULLETIN  
NO.:**

**TB UK079**

**DATE ISSUED**

29/08/2017

**BULLETIN DETAILS**

<b>Product Group:</b>	Manual	
<b>Model Affected:</b>	Quickie RXS	
<b>Effective Date: 01/03/2018</b>	<b>Serial No. Effective From:</b>	N/A
<b>Bulletin Subject:</b>	RXS discontinuation	

**Description of Bulletin:**

Technical Bulletin UK079 is confirmation of the official discontinuation of the Quickie RXS commencing from 1<sup>st</sup> March 2018.

The RXS was replaced in functionality and features by the Quickie Life in 2013. Since the introduction of the Quickie Life F, the demand for the RXS has significantly reduced and has now reached a point where the manufacturing for the RXS is no longer viable.

The RXS will continue to be available up to 1st March 2018, after this date orders will no longer be processed.

Spares or equivalents will be maintained and supplied for 5 years for the RXS.

A replacement product can be found in the Quickie Life. Please see below for more information on the benefits of switching to the Quickie Life F.

Discontinued  
RXS



Replacement  
Quickie Life



## Benefits to customers moving to the more modern Quickie Life F:

The Quickie Life F brings a significant number of benefits to the customer:

1. The standard **max user weight for the Life F is 140kg**, which is greater than that of the RXS.
2. It has a **seat depth range of 36 – 50cm and seat width of 36 – 54cm**. There is also the possibility of a smaller seat depth and width with the Quickie Life FT (teen version).
3. It is **lighter than the RXS** with the weight of the Life F from just 12kg.
4. There is a **5 year frame warranty and a 2 year guarantee on components** on the Life F.
5. The Life offers the **option of a folding back** - this gives a great value product utilising the high end technology that is currently being used by some of the really high end active manual chairs we offer. This gives a smaller folding package making the Life far easier for transportation and stowage and reducing upper body strain on the user.
6. It gives you the **option of 70, 80 and 90 degree hanger angles available**; this will cater for users wanting either a more compact set up offering a tighter turning circle or if they clinically require a tighter hanger angle due to tight hamstrings or restricted range of knee movement.
7. The Life includes a number of **Rehab and Lifestyle options** including headrest, tray table, one arm drive, **JAY Shape back**, hemi armrest and amputee pad to widen the portfolio and increase comfort for the user.
8. It has the **option of a larger 25 inch wheel for taller users** to enable a more optimal rear wheel set up. It also allows a greater reach around the wheel for improved energy efficiency. Users will travel further per push – an important benefit when propelling longer distances and to decrease the strain placed on arms and shoulders.
9. The Life offers **over 30 different colour options** for the frame and a new axle plate design in a choice of black, orange and silver, the new design not only has a more modern look, it's also lighter and increases durability. There is also a choice of coloured forks, wheels and castor rims improving the overall look of the chair and allowing the user to tailor it more to their liking.

### Literature affected:

- Parts Manual      N/A

- Order Form      RXS to be removed.

**If you have any questions regarding this bulletin, please do not hesitate to contact:**

Customer Services	Customer Support	Technical Service Centre	E-mail
<b>Telephone:</b>	<b>+44 (0)845 605 66 88 Option 2</b>	<b>+44 (0)845 605 66 88 Option 3</b>	<a href="mailto:enquiries@sunmed.co.uk">enquiries@sunmed.co.uk</a>
<b>Fax:</b>	<b>+44 (0)845 605 66 89</b>	<b>+44 (0)845 605 66 89</b>	<a href="mailto:help.technical@sunmed.co.uk">help.technical@sunmed.co.uk</a>