



Technical / News Bulletin

UK

BULLETIN NO.:

TB UK089

DATE ISSUED

15.03.18

BULLETIN DETAILS

Product Group:	Quickie Attitude		
Models Affected:	Quickie Attitude Power		
Effective Date:	March 2018	Serial No. Effective From:	
Bulletin Subject:	Quickie Attitude Power Max Speed		

Description:

The new 9-pin motor on the Attitude Power (introduced with the Attitude VIP on the 05/2017), has many advantages such as:

- Improved throttle response due to optical position sensor of the motor
- Better regulation of the throttle from the start
- Improved uphill-capacity due to new planetary gear
- New planetary gear is more reliable than the old friction-roller gear. This will result in less motor faults on the Attitude Power

The max speed of the Attitude Power with the 9-pin motor is slightly lower (21kph instead of 23kph). Although it's only 2kph, we want to be transparent with this and based on the benefits the new 9-pin motor offers in general, we see it as a big improvement overall.

Please refer to www.sunrisemedical.co.uk, or contact your Area Sales Manager for more information

Literature Affected:

Website:	Yes
Order forms:	Yes
Owner's Manual:	Yes
Parts Manual:	Yes

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 2	+44 (0)845 605 66 88 Option 3	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk