


BULLETIN DETAILS

| | | | |
|--------------------------|-----------------------|-----------------------------------|--|
| Product Group: | JAY | | |
| Models Affected: | JAY Belts & Harnesses | | |
| Effective Date: | November 2018 | Serial No. Effective From: | |
| Bulletin Subject: | JAY Belts & Harnesses | | |


In an attempt to improve, evolve and driven by customer feedback, we will be introducing two new different pelvic belt clamps, to attach JAY belts and other positioning items to the frame of our wheelchair's, (which have a tube diameter of 1" or an oval tube.)

The clamps will be readily available through the EDC and will be categorised as spares.

Fitting will be made easier by simply using a 4mm & 5mm Allen key.



Item 120750 = pair, Ø 25mm/1"



Item 116268 = each, oval shape

The above two options are due to be released on our latest order form.

Please refer to www.sunrisemedical.co.uk, or contact your Area Sales Manager for more information

Literature Affected:

| | |
|------------------------|-----|
| Website: | N.A |
| Order forms: | Yes |
| Owner's Manual: | N.A |
| Parts Manual: | N.A |

If you have any questions regarding this bulletin, please do not hesitate to contact:

| Customer Services | Customer Support | Technical Service Centre | E-mail |
|-------------------|-----------------------------|-----------------------------|-----------------------------|
| Telephone: | +44 (0)845 605 66 88 | +44 (0)845 605 66 88 | enquiries@sunmed.co.uk |
| Fax: | +44 (0)845 605 66 89 | +44 (0)845 605 66 89 | help.technical@sunmed.co.uk |