



# Technical / News Bulletin

UK

BULLETIN NO.:

TB 0578

DATE ISSUED

04/05/2019

## BULLETIN DETAILS

<b>Product Group:</b>	JAY		
<b>Models Affected:</b>	JAY Balance Air, JAY3 Cushion Air		
<b>Effective Date:</b>	04/05/19	<b>Serial No. Effective From:</b>	
<b>Bulletin Subject:</b>	ROHO Insert Valve		

At Sunrise Medical we always endeavour to offer the highest quality products to our customers, and to do this we are constantly reviewing our current range based off both internal and external feedback. To this end we have made a change to the ROHO inflation valve.

From the 4<sup>th</sup> May 2019 the ROHO inflation valve, used within our JAY Balance Air and JAY 3 Air cushions, will be made from Zytel instead of metal. Zytel is known for its high strength, abrasion and impact resistance and this improvement is designed to provide increased longevity to our products. To identify the change, new products will have a black valve, instead of silver.



Please refer to [www.sunrisemedical.co.uk](http://www.sunrisemedical.co.uk), or contact your Area Sales Manager for more information

### Literature Affected:

<b>Website:</b>	N.A
<b>Order forms:</b>	N.A
<b>Owner's Manual:</b>	N.A
<b>Parts Manual:</b>	N.A

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
<b>Telephone:</b>	<b>+44 (0)845 605 66 88</b>	<b>+44 (0)845 605 66 88</b>	enquiries@sunmed.co.uk
<b>Fax:</b>	<b>+44 (0)845 605 66 89</b>	<b>+44 (0)845 605 66 89</b>	help.technical@sunmed.co.uk