



WHERE YOU WANT TO BE.

Diamond[®]





Diamond[®]

Speed and power






Be inspired by the Diamond, a beautifully packaged 8mph scooter.

If you are looking for speed, style and comfort in a compact package, choose the Sterling Diamond.

Combining a classic look with a powerful 8mph maximum speed.



Fully adjustable captain seat



Full lights and indicators package



Compact design for easy manoeuvrability

The full suspension chassis offers an exciting ride with great stability. Large wheels and narrow profile tyres offer easy manoeuvrability on all surfaces including rough terrain.



Suspension for added comfort




Stylish alloy wheels



Dual drive for finger or thumb control

Diamond®

Speed and power



Technical data sheet

Version	4 Wheels	Weight of the heaviest part	89kg (196lbs) (base)
Speed	8mph (12.5kph)	Turning radius	145cm (57")
Batteries	50Ah	Maximum safe slope	12° - 100kg user weight (15.7 stone)
Maximum range	24 miles (38km)	Maximum safe slope	8° - 150kg user weight (23.6 stone)
Colour	Black	Ground clearance	15cm (6")
Maximum user weight	150kg (23 stone)	Tyres	Pneumatic
Total weight	108kg (238lbs)	Front wheel size	33cm (13")
Batteries weight	22kg (48.5lbs) unit	Rear wheel size	33cm (13")
Seat weight	16.9kg (37lbs)	Lights and indicators package	Yes

Find out more today!

Visit www.sterlingscooters.co.uk for more information, to find your local dealer or to request a demonstration.

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.co.uk or alternatively is available on request in large text.



Sunrise Medical Limited

Thorns Road
Brierley Hill
West Midlands
DY5 2LD

Tel: +44 (0) 845 605 66 88
Fax: +44 (0) 845 605 66 89

www.SunriseMedical.co.uk