



Technical / News Bulletin
UK

BULLETIN NO.:

TB0631

DATE ISSUED

01/04/2020

BULLETIN DETAILS

Product Group:	QUICKIE		
Model Affected:	Q-Series Family		
Effective Date: See Below	Serial No. Effective From:	N/A	
Bulletin Subject:	QUICKIE Q-Series Crash Test Dahl Update		

Description of Issue / Solution:

Since the introduction of the Q-Series powerchair range we have endeavoured to ensure that each of these chairs are crash tested with a dock and lock car system to support clients who require these.

It is for this reason we are pleased to announce that the following range of chairs have been crash tested and passed with the DAHL system:

- QUICKIE Q100 R
- QUICKIE Q200 R
- QUICKIE Q400 M Lite
- QUICKIE Q500 M Pro
- QUICKIE Q500 R Pro
- QUICKIE Q500 F Pro
- QUICKIE Q700 M Pro
- QUICKIE Q700 R Pro

For more information please use the following link to the Dahl engineering website.

<https://dahleengineering.dk/en/products/dahl-docking-systems/approved-wheelchairs-and-manuals/>

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 1-1	+44 (0)845 605 66 88 Option 1-2	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk