



Technical / News Bulletin
UK

BULLETIN NO.:

TB0636

DATE ISSUED

01/04/2020

BULLETIN DETAILS

Product Group:	QUICKIE	
Model Affected:	Attitude Manual & Hybrid	
Effective Date: See Below	Serial No. Effective From:	N/A
Bulletin Subject:	QUICKIE Attitude Grip Shift Update	

Description of Issue / Solution:

With immediate effect the 10 gear grip-shift shifter has been discontinued on the Attitude Manual and Hybrid version.

This is due to the manufacture of this part discontinuing this option. Upon review alternative suppliers unfortunately we have found that there is no short term solution as this requires fundamental changes to the Attitude bike design.

Therefore with immediate effect the Shifter grip-shift 10 speed (ATT070042) is discontinued, an alternate option can be found in the trigger shifter (ATT070045).

10 speed cassette with derailleur for Hybrid

- ATT070045 ● Shifter, trigger
- ATT070042 ☐ Shifter, grip shift, 10 speed
- ATT070048 ● Shifter mount right
- ATT070047 ○ Shifter mount left
- ATT070044 ☐ Shifter, universal, 1 lever

(base to add levers)

Therefore the QUICKIE Attitude order form will be updated and will be available on the Sunrise Medical website under the QUICKIE Attitude product page.

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 1-1	+44 (0)845 605 66 88 Option 1-2	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk