



BULLETIN DETAILS

Product Group:	Zippie		
Model Affected:	Simba		
Effective Date: See Below	Serial No. Effective From:	N/A	
Bulletin Subject:	Jay Zip Preparation on Simba		

In September 2019 we introduced a change that allowed us to interface the J3 compact mounting hardware to the JAY Zip backrest shell (TB 0604). This improved the durability of the JAY Zip bracket and allowed it to be retro-fitted with existing products.

Since that introduction, we received feedback that when fitted onto a Simba the compact hardware can sometimes clash with the push handles.

Therefore, to alleviate this problem we will be adding additional spacers onto the push handle and this will allow the backrest to be freely adjusted.



Benefits to the Customer

- Jay Compact Hardware is now more stable and durable
- Incorporating proven design
- Improved adjustability
- JAY Zip is cross compatible with J3 compact hardware

For any queries on any of these updates please contact your local Area Sales Manager.

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 1-1	+44 (0)845 605 66 88 Option 1-2	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	technical@sunmed.co.uk