



Technical / News Bulletin  
UK

**BULLETIN NO.:**

**TB 0640**

**DATE ISSUED**

15/05/2020

**BULLETIN DETAILS**

<b>Product Group:</b>	Sterling		
<b>Model Affected:</b>	Trophy		
<b>Effective Date: See Below</b>	<b>Serial No. Effective From:</b>	N/A	
<b>Bulletin Subject:</b>	Sterling Trophy Handbrake Update		

Following customer feedback we will be adding a handbrake function on the Trophy as a standard option to provide an extra safety feature for the user.

The handbrake is activated by the operation of the handle this results in braking of the scooter and at the same time an inhibitor is activated to brake.



The brake position will be mounted as follows:

- Operating handle right hand side -> brake lever right hand side
- Operating handle left and right hand side -> brake lever right hand side
- Operating handle left hand side -> brake lever left hand side

**For any queries on any of these updates please contact your local Area Sales Manager.**

Customer Services	Customer Support	Technical Service Centre	E-mail
<b>Telephone:</b>	<b>+44 (0)845 605 66 88</b> <b>Option 1-1</b>	<b>+44 (0)845 605 66 88</b> <b>Option 1-2</b>	<a href="mailto:enquiries@sunmed.co.uk">enquiries@sunmed.co.uk</a>
<b>Fax:</b>	<b>+44 (0)845 605 66 89</b>	<b>+44 (0)845 605 66 89</b>	<a href="mailto:help.technical@sunmed.co.uk">help.technical@sunmed.co.uk</a>