



Technical / News Bulletin

UK

BULLETIN NO.:

TB 0647

DATE ISSUED

31/07/20

BULLETIN DETAILS

Product Group:	Empulse		
Models Affected:	Wheeldrive		
Effective Date:	June 2020	Serial No. Effective From:	N/A
Bulletin Subject:	New Wheeldrive Battery		

The R20 Powerpack was introduced in September 2019 as the first product in the Sunrise Medical E-Mobility Portfolio, Empulse. Following on from this we have introduced our Wheeldrive Powered wheels into the portfolio also. Whilst incorporating the Wheeldrive system we have also taken the opportunity to introduce a new battery that not only contains the Empulse branding, but also introduces an improved battery across the range.

This new battery contains an improved Battery Management System (BMS) that offers the following benefits:

- Increased stability of the internal electronics for more reliable and accurate function
- More accurate counter of charging cycles
- Improved temperature control for a longer life
- Internal error log for service purposes

The improvements gained by the BMS offer improved reliability and a longer battery to the user.

The new battery EMPULSE E1801-9120 has the new spare part number 1020824 which replaces the old version.



Please refer to www.sunrisemedical.co.uk, or contact your Area Sales Manager for more information

Literature Affected:

Website:	N.A
Order forms:	Yes
Owner's Manual:	N.A
Parts Manual:	N.A

If you have any questions regarding this bulletin, please do not hesitate to contact:

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88	+44 (0)845 605 66 88	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk

