



SELF
INFLATING,
NO PUMP
REQUIRED

Equazone

Instant support. Wherever you sit.

The immediate issue, self-inflating static air cushion that can reduce pressure for all risk categories which takes away the worry of regular air cushions with its integrated foam back up.

- ✓ Outperforms more expensive competitors
- ✓ Active pressure management thanks to air-filled element, with a foam back up
- ✓ Portable and easily transferred between settings
- ✓ Split into 4 quadrants for better pressure redistribution and air flow

For more information please contact:
sales@helpinghand.co.uk

 **Helping
Hand**

Equazone Seat

Designed to enhance the performance of existing chairs, the Equazone Seat delivers full-body pressure relief with no pump required for setup.

- ✓ **EXTREMELY STABLE SUPPORT**
- ✓ **QUICK AND EASY SETUP – NO PUMP NEEDED!**
- ✓ **FOAM BACKUP FOR ADDED PEACE OF MIND**
- ✓ **REUSABLE AND EASY TO CLEAN**



FREE-MOVING AIR BACKREST

Enhances comfort while reducing pressure for optimal support.



PADDED LEGREST

Provides enhanced protection and support for comfortable leg elevation.



CROSS-SECTION CUSHION

Enhances airflow to minimise moisture buildup, while individual air quadrants provide superior stability and comfort.



EASY SET UP

The Equazone Seat arrives wrapped with clear, printed instructions for easy installation.

TECHNICAL DETAILS

- Weight limit **254kg (40 stone)**
- **Sizes** Back - 43 x 71cm, Seat - 43 x 43cm, Leg - 43 x 43cm
- **Seat Height** under 2"
- **Apaptations** flat as standard, no covers, available in 2 or 3* sections

*3 sections as standard. 2 section removes the leg cushion

The ready to issue air cushion with superb stability and no pump needed.

Equazone Cushion

The immediate issue, self-inflating static air cushion that can reduce pressure for all risk categories which takes away the worry of regular air cushions with its integrated foam back up.

- ✓ **NO PUMP NEEDED**
- ✓ **4 EXTREMELY STABLE QUADRANTS**
- ✓ **STRAPS SECURELY TO ALL CHAIRS**
- ✓ **IT PRESSURE REDUCTION & TOTAL COCCYX RELIEF**



- **TARGETED CUTAWAYS FOR OPTIMAL PRESSURE RELIEF**

Designed to relieve pressure on bony areas, including the coccyx, for enhanced comfort and support.

- **CROSS-SECTION CUSHION**

Enhances airflow to minimise moisture buildup, while individual air quadrants provide superior stability and comfort.

- **REDISTRIBUTES WEIGHT FOR COMFORT**

Pressure shifts to the back of the legs, relieving stress on sensitive areas.

- **HASSLE-FREE SETUP**

The Equazone Cushion arrives wrapped with clear, printed instructions for easy installation.

- **VERSATILE REMOVABLE STRAP**

Securely attaches the cushion to a seat or riser recliner and rolls up for easy storage.



TECHNICAL DETAILS


- Weight limit **254kg (40 stone)**
- **Medium size** 43 x 43 x 5cm (17 x 17 x 2")
- **Large size** 53 x 48 x 5cm (21 x 19 x 2")
- **Apaptations** flat as standard, no cover

Smart Pressure Distribution


Equazone distributes pressure more evenly over the full area, even when compared with far more expensive alternatives, ensuring that there are fewer areas of intense pressure. It is the reduction of these pressure points that makes Equazone so effective at reducing the risk of pressure injuries and sores.

CLIENT 1

- 80 year old female
- 12 stone
- Polio sufferer




Medium risk foam cushion




Equazone cushion

CLIENT 2

- 93 year old male
- 9 stone
- Age related




Basic issue pump up air cushion




Equazone cushion

CLIENT 3

- 66 year old female
- 24 stone
- Stroke patient




High pressure reducing gel cushion




Equazone cushion

CLIENT 4

- 57 year old female
- 17 stone



High pressure relieving air cushion (over £200)




Equazone cushion

Case Study


- Female Stroke Patient
- North Midlands
- 66 years old
- 24 stone

This client has used over ten cushions in the last few years since her stroke, between the Wheelchair Therapist and District Nursing Team they have been unable to find a solution that she finds comfortable as well as reducing the increased pressure on her right side due to her lean.


Several trials were carried out over a weekly basis on a gel cushion, static air cushion and the Equazone. All cushions pressure mapped ok for the lady but the Equazone managed to significantly reduce the pressure under the lady's right IT as well as provide her with the comfort she hadn't felt for a few years. She is now happily using the Equazone on her wheelchair and armchair at home every day, there are no longer any signs of redness on her skin but most of all she is no longer in pain when sitting.



Client's current gel cushion



Static air cushion trial



Equazone cushion



01531 635678
sales@helpinghand.co.uk
helpinghand.co.uk

