



# LIFET & LIFE RT

MATCHING A TEENAGER'S  
ACTIVE LIFE

# TWO WHEELCHAIRS – ONE FAMILY

---

Live life to the full! Meeting your friends at the mall? Having a fun time at the park? Whatever your plans are for the day, the LIFE RT and LIFE T are two wheelchairs that perfectly match the needs of a teenager's active life.

The LIFE RT is strong, reliable and offers the support you need to enjoy all kinds of cool activities. Need a folding frame? The LIFE T is not only easy to transport, but is also the perfect all-rounder for your daily routine. Regardless of your selection, both wheelchairs are designed for an adolescent's body size and needs. So whether you require a fixed-front frame (with a choice of straight, abduction or inset front ends) or you prefer a swing-away hanger version – you're sure to find a LIFE to suit you.

COOL,  
STYLISH,  
MODERN.

**LIFE**



## LIFE T THE FOLDABLE

Folding frame,  
fixed & swing away front



## LIFE RT THE SPORTY

Rigid frame,  
fixed & swing away front





→ **One-hand operating footplate lock**  
This footboard always remains in position even when you are driving over rough terrain.

# LIFE T & LIFE RT – AN EXPRESSION OF A TEENAGER’S INDIVIDUALITY

---

Customise your wheelchair to your personal taste and needs. Define your style with a choice of more than 30 frame colours, three different style packages for the fork & axle plate and the full range of stylish spoke guards. Select from a wide range of side guards, armrests, hangers and footrests. You can even fine-tune your wheelchair with a sporty camber, stylish wheels and handrims or stunning Carbon options.




→ **Wide camber range for Life T & Life RT**  
Camber improves the lateral stability of the wheelchair, protects the fingers when passing through narrow doorways and gives a sporty look. Both wheelchairs offer a camber range of 6° in 1° steps.





**// INDIVIDUAL  
CUSTOMIZATION**

Define your style with a choice of three different anodized style packages for the fork & axle plate.



**// ANODIZED 24"  
REAR WHEELS**

Accentuate your colour scheme further by matching the rim and hub of your 24" rear wheel with 3 anodized colour selections.

# LEONARDO'S LIFE



→ LIFE RT – The Sporty



→ LIFE T – The Foldable

Options and features of the LIFE T and LIFE RT – Configure your wheelchair to your individual needs. For further information on product options please check the LIFE T and LIFE RT order form.

→ For LIFE T & LIFE RT



**// ABDUCTION FRAME**

This frame version offers 2,5 cm abduction on each side. The abduction already starts below the seat upholstery.



**// CARBON PLATFORM FOOTPLATE**

Lightweight and stylish footplate with angle and depth adjustment.



**// THE CONVINCING SAFARI WHEEL LOCK**

With its intuitive handling, the safari sideguard wheel lock helps you stop quickly and safely.



**// SEAT HEIGHT ADJUSTABLE IN SEVEN POSITIONS**

With a camber from 0° to 6° and 4 centre of gravity positions (from the very active 8,3 cm to the more stable 2,3 cm).

→ For LIFE RT



**// ANGLE ADJUSTABLE BACK**

This angle adjustable backrest provides an adjustment from 78° to 98° in 4° increments.



**// IT'S ALL ABOUT ERGONOMICS**

The push bar can be easily adjusted in angle and height for an enhanced ergonomic position.



**// DEPTH ADJUSTABLE SEAT UPHOLSTERY**

Depth adjustable seat upholstery offers growth in combination with a long frame.

→ For LIFE T



**// ANGLE ADJUSTABLE BACK**

Lightweight and simple to adjust. Provides angle adjustment from 78° to 102° in 4° increments.

# TECHNICAL DATA

	Life T	Life RT
<b>Max. user weight</b>	125 kg*	125 kg*
<b>Frame version and front angle</b>	Fixed front version: 75° (straight & abduction) // 85° (straight & inset) Swing-away front version: 90°/80°/70°	Fixed front version: 75° (straight & abduction) // 85° (straight & inset) Swing-away front version: 90°/80°/70°
<b>Seat width</b>	280 - 500 mm	280 - 480 mm
<b>Seat depth</b>	320 - 460 mm	320 - 460 mm
<b>Seat height front</b>	400 - 530 mm	360 - 500 mm
<b>Seat height rear</b>	400 - 500 mm	360 - 450 mm
<b>Backrest height</b>	250 - 475 mm	250 - 450 mm
<b>Backrest angle</b>	78° - 102°	78° - 98°
<b>Wheel camber</b>	0°/1°/2°/3°/4°/5°/6°	0°/1°/2°/3°/4°/5°/6°
<b>Total width // length // height</b>	720 mm // 1075 mm // 970 mm	720 mm // 1075 mm // 950 mm
<b>Max. pivot width</b>	860 mm	700 mm
<b>Max. safe slope</b>	7° (for use of parking brakes)	7° (for use of parking brakes)
<b>Front wheel // Rear wheel</b>	4"/5"/6"/7" // 22"/24"	4"/5"/6"/7" // 22"/24"
<b>Wheelchair weight</b>	from 10,5 kg	from 10 kg
<b>Dismantled weight</b>	from 7 kg	from 8 kg
<b>Mass of heaviest part</b>	Rear wheel 1,1 kg	Rear wheel 1,1 kg
<b>Removable parts</b>	Rear wheels, hangers & sideguards (for transport only)	Rear wheels, hangers & sideguards (for transport only)
<b>Colours</b>	32 colours	32 colours
<b>Crash tested</b>	Yes (approved for transport in a vehicle according ISO 7176-19)	
<b>Transport without user</b>	Fix the wheelchair safe in the vehicle (e.g. car or plane)	
<b>Intended use and environment</b>	For people with limited mobility (children & adults); self-propelling or being pushed by an assistant, indoor and outdoor use.	

\* When selecting a proton wheel or the abduction frame max. weight is 100 kg.



## Sunrise Medical GmbH

Kahlbachring 2-4  
69254 Malsch/HD  
Tel.: +49 (0) 72 53 / 9 80-0  
Fax: +49 (0) 72 53 / 9 80-222  
E-Mail: [export@sunrisemedical.de](mailto:export@sunrisemedical.de)  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.eu](http://www.SunriseMedical.eu).

