

BULLETIN DETAILS

Product Group:	Power
Model Effected:	Zippie Salsa R, M and M2
Effective Date: 01.09.2014	Serial No. Effective From: N/A
Bulletin Subject:	Zippie Salsa Laterals Update

Description of Bulletin:

As part of the ongoing product improvement process at Sunrise Medical, the Lateral options provided with the Zippie Salsa family have been updated. The resulting improvements provide a smaller pad, removal of any rotation under pressure, and increased positioning options.

There are two options now available with the new style laterals. These are the standard in-line option and the new 2" inset option, which provides 1" inset for the laterals on both sides.

1" Inset or straight version



The new laterals can be retro-fitted to older products, and are available in the spares catalogue.

The new part numbers for the laterals are as follows:

20003005-125 Swing-away lateral support rigid - small – Pad length 125mm, incl mounting accessories - Pair

20003006-125 Swing-away lateral support inset - small – Pad length 125mm, incl mounting accessories - Pair

Literature affected:	
- Parts Manual	Updated
- Order Form	Updated and available upon request from Customer Services or your Area Sales Manager
If you have any questions regarding this bulletin, please do not hesitate to contact:	

Customer Services	Customer Support	Technical Service Centre	E-mail
Telephone:	+44 (0)845 605 66 88 Option 2	+44 (0)845 605 66 88 Option 3	enquiries@sunmed.co.uk
Fax:	+44 (0)845 605 66 89	+44 (0)845 605 66 89	help.technical@sunmed.co.uk