

WheelDrive

kompatibilitetsguide

WheelDrive, neste generasjons motoriserte hjul for manuelle rullestoler, er enkel å montere på nesten alle de vanligste manuelle rullestoler - uten ekstra deler eller tilpasninger av rullestolen, "snap-on prinsipp".

Vennligst kontakt Sunrise Medical hvis din rullestolmodell ikke er inkludert i denne kompatibilitetsguiden.



Rullestolmerke / modell	Axle type	Drivhjulsbolt diameter	Rullestolmerke / modell	Axle type	Drivhjulsbolt diameter	Rullestolmerke / modell	Axle type	Drivhjulsbolt diameter
SUNRISE MEDICAL			Action 4	stand.	12 mm	Start4 M2		
Breezy UniX ²	stand.	12,7 mm	Action 4 NG			Ventus		
Breezy PariX ²	stand.	12,7 mm	Action 3			LEVO		
Breezy BasiX ²	stand.	12,7 mm	Action XT			LEVO active-easy LAE		
Breezy RubiX ²	stand.	12,7 mm	XLT			LEVO compact-easy LCEV		
Breezy HeliX ²	stand.	12,7 mm	SpinX			LEVO Summit EL		
Breezy HeliX ² Comfort	stand.	12,7 mm	Rea Focus	stand.	12 mm	LEVO Summit M		
Breezy Cirrus G5	stand.	12 mm	Top End Terminator	stand.	12,7 mm	DIETZ		
Breezy Cirrus 4	stand.	12 mm	KÜSCHALL			AS 01		
Breezy Emineo	stand.	12 mm	Advance			O⁴		
Breezy Exigo 10	stand.	12 mm	K-Series			Multihopper	stand.	12,7 mm
Breezy Exigo 20	stand.	12 mm	K-Series attract			OlympicHopper	stand.	12,7 mm
Breezy Exigo 30	stand.	12 mm	Champion			SwayHopper	stand.	12,7 mm
Breezy Ibis	stand.	12 mm	Ultra-Light			ETAC		
Quickie Impera Line	stand.	12 mm	Compact			Cross 5		
Quickie Argon	stand.	12,7 mm	Compact attract			KI MOBILITY		
Quickie Argon ²	stand.	12,7 mm	MEYRA			Rogue		
Quickie Argon Ti	stand.	12,7 mm	Hurricane			Little Wave Click		
Quickie Easy 200	stand.	12,7 mm	ZX1 / ZX3			SORG		
Quickie Easy 300	stand.	12,7 mm	X-Star			Vector (fra SB 280 mm)		
Quickie Easy ^{max}	stand.	12,7 mm	Avanti / Avanti Pro			Vector BSA		
Quickie Helium	stand.	12,7 mm	Eurochair ² 2.750	stand.	12 mm	Jump Alpha (from SW 280 mm)		
Quickie Easy 160 i	stand.	12,7 mm	Eurochair ² 2.850	stand.	12 mm	Jump Beta		
Quickie Life i	stand.	12,7 mm	Eurochair 1.750	stand.	12 mm	VELDINK4KIDS		
Quickie Life	stand.	12,7 mm	Eurochair 1.850	stand.	12 mm	Kiddo Classic		
Quickie Life R	stand.	12,7 mm	Service 3.600	stand.	12 mm	Kiddo Tilt		
Quickie Life T	stand.	12,7 mm	Smart ^r /Smart ^s			Adapt		
Quickie Life RT	stand.	12,7 mm	X1			Kiddo Xperience		
Quickie Neon ² SA	stand.	12,7 mm	VASSILLI			PRO ACTIV		
Quickie Xenon ² /Xenon ² SA	stand.	12,7 mm	Hi-Lo M			Traveler		
Quickie Xenon / Xenon SA	stand.	12,7 mm	Hi-Lo M Alu Line			Speedy		
PANTHERA			BEROLLKA					
S2	stand.	12,7 mm	Phoenix					
S2 Swing	stand.	12,7 mm	Cobra Sport					
U2	stand.	12,7 mm	Slade SLT/AR					
U2 Light	stand.	12,7 mm	Basic					
S3			Finess					
S3 Swing			Sprint AR/SLT					
U3			OTTO BOCK					
U3 Light			Motus					
Bambino			Motus CV/CS					
INVACARE			Avantgarde CV/CS/4/VR/T					
Myon HC			Start M2 Effect					
Action 5			Start M5 Comfort					

Merk:

Listen viser hvilke rullestolmodeller Wheeldrive kan monteres på. Likevel kan variasjoner forekomme i forhold til når stolen er produsert, hvilken tyngdepunktsetting du har og hjulstørrelsen på din stol. Forutsetningen er også at akslingen for drivhjul har en størrelse på 12 mm eller 12.7 mm i diameter. Derfor er det helt nødvendig å prøve og teste ut kombinasjonen for å være 100% sikker! Spør oss om råd!

Skriv til oss: post@sunrisemedical.no