



Deliver to:

Name:

Street:

Town:

E-Mail:

Tel:

Postal code:

Fax:



Deliver to:

Name:

Street:

Town:

E-Mail:

Tel:

Postal code:

Fax:

Order date:

Purchase order:

Number required:

Marked for:

Sunrise order no:

• = Standard

./ = Price reduction

□ = Cost option

○ = No cost option

Customer service fax number:

+497253 980-222

Max. user weight:

65 kg

Art. no. : 062001-000

USER WEIGHT

kg

Mandatory information to ensure correct configuration

Base Price

CLIENT MEASUREMENTS (all measurements without cushion)

SEAT WIDTH SELECTION - measured from inside sideguard to inside sideguard (mm)

| 220* | 240* | 260* | 280 | 300 | 320 | 340 | 360 |
|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |

* only abduction frame

Each seat width can be extended
by 20 mm without extra parts


SEAT DEPTH SELECTION - measured in mm

| 240 | 260 | 280 | 300 | 320 | 340 | 360 | 380 |
|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |

Each seat depth can be extended
by 40 mm without extra parts

FRONT SEAT HEIGHT SELECTION: in dependance of front castors, Carbotecture forks and frame versions

| | | Front seat height in mm | | | | | | | | | | | | | |
|--------------|------------|-------------------------|-------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----|
| Front castor | | Fork type | Frame | 370 | 380 | 390 | 400 | 410 | 420 | 430 | 440 | 450 | 460 | 470 | 480 |
| | 3" | 111 mm | Low | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | | | | | | | | | |
| | | 111 mm | High | | | | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | | | | | |
| | 4" | 111 mm | Low | <input type="radio"/> | | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | | | | | | | |
| | | 115 mm (*) | | | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | | | | | | | | |
| | | 111 mm | High | | | | <input type="radio"/> | <input type="radio"/> | | | <input type="radio"/> | <input type="radio"/> | | | |
| | | 115 mm (*) | | | | | | | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | | | |
| | 5" | 111 mm | Low | | | | <input type="radio"/> | <input type="radio"/> | | | | | | | |
| | | 115 mm (*) | | | | | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | | | | | | |
| | | 111 mm | High | | | | | | | | <input type="radio"/> | <input type="radio"/> | | | |
| 115 mm (*) | | | | | | | | | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | | | |
| 6" | 111 mm | Low | | | | | | <input type="radio"/> | <input type="radio"/> | | | | | | |
| | 115 mm (*) | | | | | | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | | | | | | |
| | 111 mm | High | | | | | | | | | | <input type="radio"/> | <input type="radio"/> | | |
| | 115 mm (*) | | | | | | | | | | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | |



REAR SEAT HEIGHT SELECTION: in dependance of rear wheel size

| | | Rear seat height in mm | | | | | | | | | | | | | |
|------------|--|------------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| Rear Wheel | | 290 | 300 | 310 | 320 | 330 | 340 | 350 | 360 | 370 | 380 | 390 | 400 | 410 | 420 |
| 20" (**) | | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 22" | | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 24" | | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |

The rear seat height can only be between 0 cm and 13 cm below the selected front seat height.

(*) fork compatible with additional drive
(**) only with Proton wheels/Right Run tyres

For more information
visit our webpage
www.SunriseMedical.eu
© Sunrise Medical 2025

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat width (SW):

- a) For the user: Allow space between sideguard and hip on each side of the chair (approx. 10 - 20 mm altogether).
- b) On the wheelchair: From inside sideguard to inside sideguard.

Seat depth (SD):

- a) For the user: When sitting upright: leave space for 2 - 3 fingers between front edge of the seat upholstery and back of knee. (approx. 30 - 60 mm). Care must be taken with adj. upholstery.
- b) On the wheelchair: From back tube to front of seat upholstery.

Front seat height:

- a) For the user: When sitting with 90° kneeangle:
From back of knee to bottom of shoe sole.
- b) On the wheelchair: From front of seat tube to the ground
(min. 30 mm longer than user measurement).
Note: Remember to include the cushion height.

Rear seat height:

- b) On the wheelchair: From back of seat tube to the ground.
Note: Max. seat height difference rear to front is 130 mm.

Back height:

- a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) On the wheelchair: From back of seat tube to upper edge of backrest upholstery.
Note: The user's ability to sit upright has a major influence on the back height and it can therefore differ considerably.

Back angle:

Important: The back height depends on the degree of disability.
0° is the angle 90° to the ground. Negative angles are <90°
and positive angles are >90°.

Technical Data Overview:

| | | | | |
|--------------------------|----------------|-----------------------|-----------------------|-----------------------|
| Max. overall height: | 1000 mm | | | |
| Max. overall length: | 850 mm | | | |
| Max. turning circle: | 650 mm | | | |
| Max. safe slope: | 10° | | | |
| Wheelchair weight: | from 9,0 kg | | | |
| Dismantled weight: | from 6,9 kg | | | |
| Weight of heaviest part: | 2,1 kg | | | |
| Max. overall width: | Camber | 20" rear wheel | 22" rear wheel | 24" rear wheel |
| | For 0° camber: | SW + 160 mm | SW + 160 mm | SW + 160 mm |
| | For 3° camber: | SW + 200 mm | SW + 210 mm | SW + 220 mm |
| | For 6° camber: | SW + 260 mm | SW + 270 mm | SW + 280 mm |
| | For 9° camber: | SW + 310 mm | SW + 320 mm | SW + 330 mm |

Product is **crash tested** according to **ISO 7176-19** and therefore approved for transportation in a vehicle.

For transport without a user the chair must be safe clamped in the vehicle.

All changes in the configuration might change the seat heights.

Manufacturing tolerances of +/- 10 mm

Zippie Simba

● = Standard ○ = No cost option □ = Cost option

Price in €

SEAT

Seat sling

| | | |
|--------------------------------------|--|--|
| SIB020002 | ● Seat sling, nylon, 1 utility bag | |
| SIB020004 | □ Seat sling, nylon, 1 utility bag, 1 catheter bag | |
| SIB020000 | ○ Without seat sling | |
| For Jay cushions see Jay order forms | | |

Seat cushion, cover with zipper, black

| | | | | | | |
|--------------------|------------|-----------|---------|-----------|---------|--|
| Foam (medium-soft) | Thickness: | SIB020101 | □ 30 mm | SIB020102 | □ 50 mm | |
|--------------------|------------|-----------|---------|-----------|---------|--|

Seat accessories

| | | |
|-----------|-------------------------------|--|
| SIB090018 | □ Lap belt with buckle, metal | |
|-----------|-------------------------------|--|

BACKREST

Back height range

Choose back height factory setting (mm).
Height adjustments of in total 5 cm possible.

| 200 - 250 | 250 - 300 | 300 - 350 | 350 - 400 |
|-----------|-----------|-----------|-----------|
| ○ 200 | ○ 250 | ○ 300 | ○ 350 |
| ○ 225 | ○ 275 | ○ 325 | ○ 375 |
| ○ 250 | ○ 300 | ○ 350 | ○ 400 |

| | | | |
|-----------|--|--|--|
| SIB030015 | ● Backrest, angle-adjustable, fold-down | (-16° to +20°) | |
| SIB030102 | □ Backrest, angle-adjustable, lock when folded | | |
| SIB030530 | □ Backrest angle 100° (Relax position) | (not in combination with double locking) | |

Push handles

| | | |
|-----------|--|---|
| SIB030200 | ○ Without push handles | |
| SIB030201 | ○ Push handles long | |
| SIB030203 | □ Push handles, fold-down | (only available from back height 300 mm) |
| SIB030204 | ● Push handles, height-adjustable | |
| SIB030210 | ○ Preparation for JAY Zip back | (not with Distance 20 mm for stabilising bar) |
| SIB030208 | □ Push bar, height-adjustable, width and angle-adjustable | (not with SW 220 mm) |
| SIB030210 | □ Preparation for JAY Zip back (available with height-adjustable push handles) | |

Backrest upholstery

| | | |
|-----------|--|--|
| SIB030316 | ● EXO Backrest sling, padded, adjustable, 1 accessory bag, black (EXO Evo version) Colour, binding tape: ● Black ○ Orange ○ Silver ○ Red ○ Blue | |
| SIB030317 | □ EXO PRO Backrest, padded, adj., breathable spacer fabric, 1 accessory bag, black (EXO Evo version) Colour, binding tape: ● Black ○ Orange ○ Silver ○ Red ○ Blue | |

Accessories for backrest

| | | | |
|-----------|--|---|----------------------------|
| | □ Headrest, depth- and angle-adjustable* | (only available from back height 300- 350 mm) | |
| SIB030411 | □ Small pad (200 x 110 mm) | SIB030410 | □ Large pad (230 x 110 mm) |

* Not available with JAY Zip Back, see back accessories on separate JAY Zip Back order form)

SIDEGUARDS AND ARMREST

Sideguards

| | | |
|-----------|--|--|
| SIB040103 | ● Aluminium sideguard, fender, wet weatherprotection (in frame colour) | |
|-----------|--|--|

Armrest (pair)

| | | | |
|-----------|------------------------------------|--|--|
| SIB040055 | □ Armpad, height-adjustable, short | (only with aluminum sideguard; not available with the Safari wheel lock) | |
|-----------|------------------------------------|--|--|

Accessories for sideguards and armrests

| | | | |
|-----------|---|---------------------------------|--------|
| SIB090111 | □ Tray, wood, elevated frame, push-fit, SW 220 - 360 mm | (not with aluminium sideguards) | |
| | □ Tray, plastic, swing-away, SW 260 - 380 mm | (not with aluminium sideguards) | |
| | Mounted: □ right | SIB090113 | □ left |
| | | SIB090112 | |

For more information
visit our webpage
www.SunriseMedical.eu
© Sunrise Medical 2025



Seat sling with 1 utility bag



Lap belt with buckle



Push handles
long



Push handles
fold-down



Push handles
height adj.



Push bar - height, angle
and width adjustable



Aluminium sideguard
with fender
(here also with Safari wheel lock)

Zippie Simba

● = Standard ○ = No cost option □ = Cost option

Price in €

FOOTBOARDS

Footboard, platform, lightweight

| | | | | |
|-----------|-------------------------------------|---|----------------------------|--|
| SIB050025 | <input type="checkbox"/> | Plastic, angle and depth-adjustable, calf strap, black | (LLL between 230 - 360 mm) | |
| SIB050056 | <input type="checkbox"/> | Aluminum, angle and depth-adjustable, calf strap, black | | |
| SIB050032 | <input type="checkbox"/> | Plastic, angle and depth-adjustable, flip-up (to the back), calf strap, black | (LLL between 230 - 360 mm) | |
| SIB050048 | <input checked="" type="checkbox"/> | Plastic, angle and depth-adjustable, high-mount, calf strap, black | (LLL 170 - 250 mm) | |

Accessories for footboards

| | | | | |
|-----------|--------------------------|--|---------------------------------|--|
| SIB050127 | <input type="checkbox"/> | Adjustable side protection, to create a side-stop on the footboard | (not with swing back footplate) | |
|-----------|--------------------------|--|---------------------------------|--|

LOWER LEG LENGTH (LLL): LLL has to be 25 mm less than front seat height; all measurements are in mm

- 170 ○ 180 ○ 190 ○ 200 ○ 220 ○ 230 ○ 240 ○ 250
○ 260 ○ 270 ○ 280 ○ 290 ○ 300 ○ 310 ○ 320 ○ 330
○ 340 ○ 350 ○ 360

FRAME VERSIONS

| | | | | |
|-----------|-------------------------------------|--|------------------------------------|--|
| SIB100118 | <input checked="" type="checkbox"/> | Frame FF 95° (integrated footrest), straight | (available from seat width 280 mm) | |
| SIB100119 | <input type="checkbox"/> | Frame FF 95° (integrated footrest), abducted (300 mm on each side) | | |

CENTRE OF GRAVITY (COG)

(measured from front of back post in mm)

| | | | | |
|--------|--------------------------|--------|--|--|
| COG20 | <input type="checkbox"/> | 20 mm | | |
| COG30 | <input type="checkbox"/> | 30 mm | | |
| COG40 | <input type="checkbox"/> | 40 mm | | |
| COG50 | <input type="checkbox"/> | 50 mm | | |
| COG60 | <input type="checkbox"/> | 60 mm | | |
| COG70 | <input type="checkbox"/> | 70 mm | | |
| COG80 | <input type="checkbox"/> | 80 mm | | |
| COG90 | <input type="checkbox"/> | 90 mm | | |
| COG100 | <input type="checkbox"/> | 100 mm | | |

FRAME ACCESSORIES

(choose either tip assist or anti-tip tubes, not both)

| | | | | |
|-----------------------------------|--------------------------|--------------------------------------|--|--|
| SIB090100 | <input type="checkbox"/> | Tie down brackets for transportation | (includes SIB030405) | |
| Tip assist | | | | |
| SIB090001 | <input type="checkbox"/> | Left | SIB090002 <input type="checkbox"/> Right | |
| SIB090009 | <input type="checkbox"/> | Crutch holder, loop | (only in combination with tip assist) | |
| Anti-tip tubes, swing-away | | | | |
| SIB090004 | <input type="checkbox"/> | Left | SIB090005 <input type="checkbox"/> Right | |
| SIB090006 | <input type="checkbox"/> | Pair | | |

CASTOR FORK

(compare with your choice on the first page)

| | | | | |
|-----------|--------------------------|-----------------------|---------------|--|
| SIB080013 | <input type="checkbox"/> | Fork Carbotechnology® | (only 111 mm) | |
| SIB080029 | <input type="checkbox"/> | Fork Aluminum | | |

STANDARD CASTORS

(compare with your choice on the first page)

| | | | | |
|-----------|--------------------------|--|--------------------|--|
| SIB080301 | <input type="checkbox"/> | Castor, solid (hard), rim: composite | | |
| SIB080302 | <input type="checkbox"/> | Castor, solid (soft), rim: composite | (only for 5" & 6") | |
| SIB080303 | <input type="checkbox"/> | Castor, solid (hard), light-up, rim: composite | (not with 6") | |
| SIB080304 | <input type="checkbox"/> | Castor, solid (hard), rim: aluminium | (only for 4" & 5") | |
| SIB080306 | <input type="checkbox"/> | Castor, solid (soft), rim: aluminium | (not with 3") | |

REAR WHEELS

(compare with your choice on the first page; choose only one wheel type)

| | | | | |
|-----------|--------------------------|---|------------------|--|
| SIB070001 | <input type="checkbox"/> | Rear wheel, Universal, 36 spokes, cross-wise spokes | (only 22" & 24") | |
| SIB070002 | <input type="checkbox"/> | Rear wheel, Design, 36 spokes, straight spokes | (only 24") | |
| SIB070003 | <input type="checkbox"/> | Rear wheel, Lightweight, 24 spokes, straight spokes | (only 22" & 24") | |
| SIB070004 | <input type="checkbox"/> | Rear wheel, Proton, 24 spokes, straight spokes | | |
| SIB070009 | <input type="checkbox"/> | Rear wheel, Mountainbike, 36 spokes, straight spokes | (only 24") | |
| SIB070010 | <input type="checkbox"/> | Mountainbike wheel, replacement kit incl. tyre (extra coarse profile) | (only 24") | |

TYRES

(choose only one tyre type)

| | | | | |
|-----------|--------------------------|--|-----------------|--|
| SIB070101 | <input type="checkbox"/> | Tyre, pneumatic, fine profile, Schwalbe Right Run | | |
| SIB070102 | <input type="checkbox"/> | Tyre, puncture proof, tread | (only 22", 24") | |
| SIB070107 | <input type="checkbox"/> | Tyre, pneumatic, slick, Schwalbe One | (only 24") | |
| SIB070104 | <input type="checkbox"/> | Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution | (only 24") | |

| | | | | |
|-----------|--------------------------|--|------------------------|--|
| SIB090163 | <input type="checkbox"/> | Spares-kit: 1 extra tyre (same as on the chair) and 1 extra inner tube | (not with solid tyres) | |
|-----------|--------------------------|--|------------------------|--|

For more information
visit our webpage
www.SunriseMedical.eu
© Sunrise Medical 2025



Plastic, angle and depth-adjustable



Adjustable side protection



Tip assist



Anti tip tubes,
swing away



Rear Wheel
Universal



Rear Wheel
Design



Rear Wheel
Lightweight



Rear Wheel
Proton



Rear Wheel
Mountainbike



Right Run
Fine Profile



Schwalbe
One



Marathon Plus



Puncture Proof

| | | | | | | |
|----------------------------|---|--|---|---|------------|--|
| Zippe Simba | | | ● = Standard ○ = No cost option □ = Cost option | | Price in € | |
| HANDRIMS | | | | | | |
| SIB070202 | ○ | Handrims, aluminium, anodised, silver | | | | |
| SIB070201 | ○ | Handrims, aluminium, powder-coated, clear | | | | |
| SIB070204 | ○ | Handrims, stainless steel, silver | | | | |
| SIB070203 | □ | Handrims, aluminium, powder-coated, coloured | (select colour on last page) | | | |
| SIB070205 | □ | Handrims, titanium, untreated | (only 24") | | | |
| SIB070206 | □ | Handrims, titanium, powder-coated, clear | (only 24") | | | |
| SIB070207 | □ | Handrims, Supergrip ® | (only 22" & 24") | | | |
| SIB070208 | □ | Handrims, Maxgrepp ® | (only 22" & 24") | | | |
| Hand rim mounting | | | | | | |
| SIB070307 | □ | Handrim mount, super narrow (7 mm) | (only with stainless steel; not in combination with Spinergy, Proton and Lightweight Wheels and stainless steel Handrim, not with Maxgrepp) | | | |
| SIB070317 | ● | Handrim mount, narrow (17 mm) | Not with Maxgrepp | | | |
| SIB070327 | ○ | Handrim mount, wide (27 mm) | | | | |
| SIB070250 | □ | Handrim cover, rubber | (only 22" and 24") | | | |
| WHEEL CAMBER | | | | | | |
| (choose only one option) | | | | | | |
| SIB070600 | ○ | 0° | | | | |
| SIB070603 | ○ | 3° | | | | |
| SIB070606 | ○ | 6° | | | | |
| SIB070609 | ○ | 9° | | | | |
| SPOKE GUARDS | | | | | | |
| (available in 22" and 24") | | | | | | |
| SIB070831 | □ | Transparent | (also available in 20") | | | |
| SIB070813 | □ | Glaktica | | | | |
| SIB070823 | □ | Choose back height factory setting (mm). | SIB070818 | □ Jungle Fever | | |
| SIB070816 | □ | Happy Lions | SIB070821 | □ Pirates | | |
| SIB070811 | □ | Feather | SIB070811 | □ Feather | | |
| SIB070819 | □ | Monkeys | SIB070814 | □ Graffiti | | |
| SIB070806 | □ | Chrome | SIB070806 | □ Chrome | | |
| SIB070820 | □ | Oriental | SIB070802 | □ Basketball | | |
| SIB070803 | □ | Basketball 2 | SIB070803 | □ Basketball 2 | | |
| SIB070829 | □ | Superhero | SIB070812 | □ Flying Water | | |
| SIB070828 | □ | Space | SIB070828 | □ Space | | |
| SIB070817 | □ | Hot Wheels | SIB070833 | □ Zebra | | |
| SIB070822 | □ | Pool | | | | |
| SIB070827 | □ | Silk | | | | |
| SIB070830 | □ | Tiger | | | | |
| WHEEL LOCKS | | | | | | |
| (choose only one option) | | | | | | |
| SIB060001 | ● | Standard brake | | | | |
| SIB060002 | □ | Knee-lever brake | | | | |
| SIB060010 | □ | Brake lever extension, knee-lever brake, fold-down | | | | |
| SIB060003 | □ | Compact brake | (only available on straight frame) | | | |
| SIB060013 | □ | Push to lock brake, integrated in sideguard, Safari | (not with solid tires) | | | |
| SIB060014 | □ | Additional brake actuation for attendant for Safari brake, integrated in sideguard | (only with safari wheel lock) | | | |
| ACCESSORIES | | | | | | |
| SIB090000 | ● | Tool kit | | | | |
| SIB090024 | □ | High pressure pump (0 - 14 bar) | | | | |
| SIB090025 | □ | Air pump | | | | |
| SIB090030 | □ | Back pack | | | | |
| SIB090016 | □ | Frame protection, black | (not with platform footplates highmounted) | | | |
| SIB090026 | □ | Mobile phone pocket | | | | |
| | | | | For more information visit our webpage www.SunriseMedical.eu © Sunrise Medical 2025 | | |



Handrims aluminium
Silver anodized



Handrims aluminium
powder coated, coloured



Handrim
Maxgrepp



Standard brake



Knee-lever brake



Compact brake



Push to lock brake
Safari



Tool kit



High pressure pump



Air pump



Back pack

Zippie Simba

● = Standard ○ = No cost option □ = Cost option


COLOURS

| | Frame | Handrim | Hub/rim | | Frame | Handrim | Hub/rim |
|------------------------|-------|---------|---------|----------------------------|-------|---------|---------|
| 02 Brilliant Silver | ○ | □ | | 27 Graphitemetallic | ○ | □ | |
| 03 Silver Matt | ○ | □ | | 29 Black | ○ | □ | □ |
| 06 Sun Yellow | ○ | □ | | 31 Matt Black | ○ | □ | □ |
| 14 Blue | ○ | □ | □ | 32 Signal White | ○ | □ | |
| 15 Brilliant Marine | ○ | □ | | 33 Pink | ○ | □ | |
| 16 Matt Marine | ○ | □ | | 34 Dark Chocolate Metallic | ○ | □ | |
| 19 Brilliant Purple | ○ | □ | | 35 Skyblue | ○ | □ | |
| 21 Red Ruby | ○ | □ | | 36 Cyan-Metallic | ○ | □ | |
| 22 Brilliant Red | ○ | □ | □ | 37 Marsred (Orange) | ○ | □ | |
| 23 Matt Red | ○ | □ | □ | 40 Kiwi | ○ | □ | |
| 24 Tutti Frutti Orange | ○ | □ | □ | 41 Flamingo-Glitter | ○ | □ | |
| | | | | 42 Matt Anthracite | ○ | □ | |
| | | | | 46 Mamba Green (new) | ○ | □ | |

STYLE PACKAGE - colours for forks

SIB100502 ● Black
SIB100503 □ Silver

Comments for special build - please apply for a quote; lead time conditions my change.



Simba - Base Price

| | |
|-------------------------|--|
| Option price | |
| Packaging and transport | |
| Exclusive VAT | |
| VAT. | |
| Total Price | |

Sunrise Medical Customer service

Address: Ul. Walczaka 114a,66 400, Gorsaw, Gorzów Wlkp. Polen

Tel. +48 95 737 25 14

E-mail: cs.pl@sunrise-medical.pl

Website: sunrisemedical.eu



www.SunriseMedical.eu

© Sunrise Medical 2025