

ASSEMBLY INSTRUCTION

External Joystick / External Joystick Lite



Description

Article number

External Joystick - Mounting on armrest, same side as Omni display

☐ 1016670

External Joystick Lite - Mounting on armrest, same side as Omni display

☐ 1016510



							
Included in kit	1015013	1015014	1016675	00355.0477	0035.0471	1016716	1016718
1016670	1		1	1	1	1	1
1016510		1	1	1	1	1	1



EN

1. Assemble the Omni screen bracket and the swing-away controller bracket.
2. Fasten the joystick to the swing-away bracket.
3. Fasten the black and red buddy buttons to the armrest.
4. Set up the cabling.

 In order to set up the system, see R-net Omni technical manual.

NO


1. Monter Omni skjermbrakett og utsvingbar styreboksbrakett.
2. Fest joysticken til den utsvingbare braketten.
3. Fest de svarte og røde Buddy-knappene til armlenet.
4. Sett opp kablene.


 Se den tekniske manualen for R-net Omni for oppsett av systemet.

SE

1. Montera fästet för Omni-skärmen och det undansvängbara controller fäste.
2. Fäst manöverboxen på det undansvängbara fästet.
3. Fäst de svarta och röda buddy-kontakterna till armstödet.
4. Installera kablaget.

 För att installera systemet, se R-net Omni teknisk manual.





Sunrise Medical HCM B.V.
Vossenbeemd 104
5705 CL Helmond
The Netherlands

T: +31 (0)492 593 888
customerservice@sunrisemedical.nl
www.SunriseMedical.eu